



## Upright Solid Door Laboratory Refrigerators

The **FUTURA Silver Series general purpose refrigerators** are intended for general purpose applications only, where precise temperature control and maintenance are not essential, and there are no requirements for temperature alarms. These economical units are available in several sizes that fit any storage requirements.



LABH-30-SD

### Standard Product Features:

- 18 & 30 Cu. Ft. Capacities
- Automatic defrost
- Adjustable mechanical thermostat
- 0°C to +10°C operating temperature; factory preset to +4°C
- Magnetic door gasket for positive seal
- White textured cabinet and door
- Right hand swing, solid door with heavy duty handle & key lock
- Leveling legs
- White cabinet and door
- Forced draft air circulation
- High density urethane foam cabinet and door insulation
- CFC-free refrigerant (R134a), foam insulation and packaging
- 115V, 60Hz operation
- Probe access port (3/8")
- UL/C-UL Listed
- Warranty: 2 years parts & labor, 5 years compressor parts only

### Optional Features:

- Casters
- Chart recorder
- Additional shelves
- Digital display

**LABRepCo**



**LABH-18-SD**

Model	LABH-18-SD	LABH-30-SD
Net Weight lbs.	157	312
Overall Height in.	67 3/4	72 1/2
Overall Width in.	34	35 3/4
Overall Depth in.	31 1/2	34 1/4
Interior Height in.	51 1/2	65
Interior Width in.	27	32
Interior Depth in.	17 3/4	22 1/2
Interior cu. ft.	18	30
Temperature	+2°C to +10°C	+2°C to +10°C
Number of Shelves	4 fixed	5 adjustable
Doors (Ext/Int)	1/0	1/0
Compressor Size HP	1/4	1/4
Receptacle	NEMA 5-15R	NEMA 5-15R
Total Amp Draw	5	6.8
Voltage	115V, 1PH, 60Hz	115V, 1PH, 60Hz
Listing	UL/C-UL	UL/C-UL

*Please Note: Interior dimensions are estimates.  
All specifications subject to change without notice.*



**LABH-18-SD (interior)**

Receptacle Required



115V/1PH/ 60Hz  
NEMA-5-15P



**LABRepCo**