

PLATE FUGE™

MICROPLATE MICROCENTRIFUGE



- Quickly spins down droplets in PCR plates
- Unique rotor design prevents sample spillage
- 50% smaller than traditional centrifuges
- Accepts all popular PCR plates (skirted, semi & non) & most microtiter plates



Before/After Centrifugation During Centrifugation

15% off

C2000 PlateFuge™ MicroCentrifuge with rotor and plate carrier ~~\$715.00~~ **\$607.75**

ORBI-SHAKER™



ORBI-SHAKER™ CO₂

All Environment Shaker



- Digital speed control up to 300 rpm
- TOUCH Screen control over time, speed and temperature
- Instantly exchange flask clamps with MAGIC Clamp™ accessories (PATENTED)
- Available in remote control version for use in CO₂ Incubators

BT3001
BT4001
BT3000-MR

Orbi-Shaker™ with rubber mat platform
Orbi-Shaker™ CO₂ w/remote controller and rubber mat platform
MAGIC Clamp™ universal platform (14x12") for flasks & tube racks

15% off

~~\$1599.00~~ **\$1359.15**
~~\$3199.00~~ **\$2719.15**
~~\$248.00~~ **\$210.80**

NEW! Roto-THERM™ Plus

INCUBATED TUBE ROTATOR

& Roto-Mini™ Plus

Specifications:

Temp Range: Ambient +5° to 60°C
Speed Range: 5 to 70rpm (Plus model)
Max. Capacity: 6 x 50ml, 16 x 15ml,
24 x 0.5ml, 44 x 1.5/2.0ml
96 x 0.2ml

- Temperature control up to 60°C (H2024 only)
- Includes tube holders for 0.5 to 50ml
- Available with variable speed & mix modes
Oscillating • Pausing • Rocking

Included Holders

- 11 x 1.5-2ml (2/pk.)
- 5 x 5-7ml (2/pk.)
- 4 x 10-15ml (2/pk.)
- 6 x 50ml* (1ea.)

*Also accepts 24 x 0.5ml

Includes Tube Holders for All Common Tube Sizes



15% off

H2024 Roto-Therm™ Plus w/included tube holders (variable speed) ~~\$699.00~~ **\$594.15**
R2024 Roto-Min™ Plus, w/included tube holders (variable speed) ~~\$425.00~~ **\$361.25**
R2020-PCR Optional tube holder for 96 x 0.2 (12 x 8 position PCR strips) ~~\$69.00~~ **\$58.65**



H2024

R2024

Additional tube holders available, please visit our website for a complete listing.

NEW! HOTPLATE STIRRER



- Large white ceramic top plate
- Digital speed and temp. control
- Safety LED indicates hot surface
- Stir capacity up to 15L

Max. Speed: 1500rpm
Max. Temp.: 400°C
Platform: 10x10"

Buy 2, Get 1 FREE!

H3710-HS HotPlate Stirrer, 10 x 10-inch with white ceramic top plate \$639.00

Use Promo Code: B2G1-HP

BENCHWAVER™

- Undulating 3D motion
- Adjustable speed and tilt
- Ideal for staining, washing and general mixing



Optional Stacking Platform

Speed: 5-105rpm
Platform: 13.5 x 13 in.
Tilt Angle: 0 to ±10°



B3D5000 BenchWaver™ Undulating Shaker \$1249.00
B3D5000-STK Optional stacking platform \$127.00
B3D5000-DIMP Optional dimpled mat \$84.00

Buy 2, Get 1 FREE!

Use Promo Code: B2G1-BW



NEW! BEADBLASTER™ 24R

REFRIGERATED MICROTUBE HOMOGENIZER

Use Promo Code: D2400-25
D2400R-R5

- Completely lyse, grind and homogenize a variety of sample types
- New! **POWERFUL REFRIGERATION** to protect delicate samples
 - True refrigeration - no dry ice or compressed air needed
 - Space saving design
- Process up to 24 samples simultaneously
- Wide variety of grinding media in prefilled 2ml and 5ml tubes

Specifications:

Capacity: 24 x 2ml tubes
 Speed range: 4.00 - 7.00 m/sec
 Cooling: Compressor*
 Temp. Range: -10°C to Ambient*
 Temp. Uniformity: ±2°C
 Warranty: 2 years

*D2400-R Only

D2400	BeadBlaster 24	\$8295.00	\$6221.25
D2400-R	BeadBlaster 24R with refrigeration	\$12990.00	\$9742.50
D2400-R5	Optional tube holder for 5ml tubes	\$359.00	\$322.15
D1032-05	Prefilled 2ml tubes, 0.5mm zirconium	\$155.00	\$131.75
D1032-15	Prefilled 2ml tubes, 1.5mm zirconium	\$155.00	\$131.75
D1032-30	Prefilled 2ml tubes, 3mm zirconium	\$155.00	\$131.75
D1033-28	Prefilled 2ml tubes, 2.8mm steel	\$155.00	\$131.75
D1034-30	Prefilled 5ml tubes, 3mm zirconium	\$219.00	\$186.15
D1034-28	Prefilled 5ml tubes, 2.8mm steel	\$207.00	\$175.95

Additional grinding media is available in prefilled tubes and loose bulk. Ask for details.

25% off

15% off



MC-24™ Touch

High Speed Microcentrifuge

- Full Color Touch Screen
- mySpin™ 1 button program recall
- COMBI-Rotor™ accepts tubes and strips
- Compact, only 9 x 12 x 7.75 in.
- Variable speed control, set in rpm or rcf

C2417 MC-24™ Touch Microcentrifuge with 24 place COMBI-Rotor ~~\$2115.00~~
\$1797.75 (15% Off)

Max. Speed: 13,500rpm / 16,800xg
 Max Capacity: 24 x 1.5/2.0ml
 2 x PCR Strips (16 x 0.2ml)



Touch Control

Unique Rotor for microtubes & PCR STRIPS



INCU-SHAKER™

SHAKING INCUBATORS

- Large flask capacity
- **MAGic Clamp** accessories for instant-magnetic clamp exchange
- Available with CO₂ control or cooling
- **NEW!** Adjustable clamp for 125, 250 and 500ml flasks



INCU-SHAKER Mini

Speed Range: 30 to 300rpm
 Temp Range: Amb. +5° to 70°C
 Platform Dims: 9.5 x 11.5 in.
 Universal Platform: H1000-MR
 Flask Capacity: 2 x 1L (H1000-MR-1000)
 5 x 500ml (H1000-MR-500)
 8 x 250ml (H1000-MR-250)
 12 x 125ml (H1000-MR-125)
~~\$2999~~
\$2549.15



Item: H1001-M

15% off

INCU-SHAKER 10L & 10LR

Speed Range: 30 to 300rpm
 Temp Range: (H1012, Cooling Model) Amb. -15° to 60°C
 Platform Dims: 17.8 x 17.8 in.
 Universal Platform: H1010-MR
 Flask Capacity: 5 x 2L (H1000-MR-2000)
 9 x 1L (H1000-MR-1000)
 16 x 500ml (H1000-MR-500)
 25 x 250ml (H1000-MR-250)
 36 x 125ml (H1000-MR-125)
~~\$6115 / \$8565~~
\$5197.75 / \$7280.25



Item: H2010 & H2012 (Cooling)

15% off

INCU-SHAKER Mini CO₂

Speed Range: 30 to 300rpm
 Temp Range: Amb. +5° to 60°C
 CO₂ Range: 0 to 20%
 Platform Dims: 9.5 x 11.5 in.
 Universal Platform: H1000-MR
 Flask Capacity: 5 x 500ml (H1000-MR-500)
 8 x 250ml (H1000-MR-250)
 12 x 125ml (H1000-MR-125)
~~\$8899~~
\$7564.15



Item: H3501

15% off



BENCHMIXER™ V2

- Reversible Top: Use "Cup" or "Flat/Dimpled" surface
- Powerful motor for INSTANT vortexing
- Variable speed (touch or continuous)

BV1003

\$279

Use Promo Code: B2G1-V2

Buy 2, Get 1
FREE!



PATENT PENDING

LABRepCo

One Company • Many Solutions

Ph: 1-800-521-0754

Email: info@labrepco.com