

Abeyance Cryo Solutions | Preserving Potential

SIMPLE & SECURE -190°C LN2 VAPOR STORAGE

DESIGNED FOR SAMPLES BY USERS

INCREASE CAPACITY | IMPROVE ERGONOMICS | STAY CONNECTED





Simple Controls, Secure Storage

Touchscreen controller with WiFi, LAN, text and email alerts automatically maintains the freezer conditions. Critical data is recorded and can be pushed to the cloud, emailed or saved to a USB flash drive. Supports independent temperature and level monitoring for redundant peace of mind.

Increase Capacity & Improve Ergonomics

Cryogenic LED and automatic fog clear enhance visibility and sample access. Optimized designs minimize liftover height and balance capacity and usability with sufficient workspace to maintain the cold chain.







Abeyance produces innovative cryo solutions for life science research and therapy.

Successful use of samples in the lab or clinic is the ultimate goal and cryopreservation is a critical link. Decades of experience across the cold chain, vertical integration and a commitment to quality make our unique freezers possible.

Increase capacity, improve ergonomics and stay connected with simple, secure cryo storage designed for samples by users.

Proudly made in the USA with domestic materials and skilled craftsmanship, we are dedicated to preserving your sample potential!

| | A220 | A440 | A700 | A1000 |
|---------------------------|-------------------|-------------------|--------------------|-----------------------|
| 2 ml Vial Capacity | | | | |
| 5-2 Slim Racks | 22,000 220 44 | 44,000 440 88 | 70,000 700 140 | 100,000 1,000 200 |
| 13-2 Square Racks | 15,600 156 12 | 41,600 416 32 | 67,600 676 52 | 98,800 988 76 |
| (Vials Boxes Racks) | | | | |
| | | | | |
| Outer Diameter - in (mm) | 34.0 (864) | 45.0 (1143) | 55.0 (1397) | 65.0 (1651) |
| Overall Height - in (mm) | 58.0 (1473) | 59.0 (1499) | 63.0 (1600) | 64.0 (1625) |
| Usable Height - in (mm) | 30.0 (762) | 30.0 (762) | 30.0 (762) | 30.0 (762) |
| Step Height - in (mm) | 9.5 (241) | 10.5 (267) | 10.0 (254) x 2 | 10.5 (267) x 2 |
| Liftover Height - in (mm) | 39.5 (1003) | 39.5 (1003) | 34.0 (864) | 34.0 (864) |
| | | | | |
| LN2 Capacity | 65 L | 150 L | 270 L | 390 L |
| Hold Time (< Tg) | > 10 Days | > 15 Days | > 25 Days | > 25 Days |
| | | | | |

Warranty

3 Year Parts and Craftsmanship | 7 Year Vacuum Insulation

Safety / Regulatory

Conforms to UL STD 61010-1 | CSA STD C22.2 # 61010-1 | LVD (2014/35/EU)



