


Smart storage solution designed for Pfizer's Covid-19 Vaccine Boxes An ideal tool to keep vaccines organized and ready for retrieval



UFD-252-T2
Stores (10) Pfizer
Covid-19 Vaccine Boxes

UFD-352-T6
Stores (15) Pfizer
Covid-19 Vaccine Boxes

All freezer racks come with our unconditional lifetime warranty

Freezer Racks are made from pharmaceutical grade stainless steel and feature double welds. All edges are smoothed for easy handling.

Part #	Dimensions (D x W x H) (in)	Dimensions (D x W x H) (in)	Capacity
UFD-252-T2	18.42 x 9.45 x 9.29 (in)	468 x 240 x 236 mm	10 Boxes
UFD-352-T6	27.48 x 9.45 x 9.29 (in)	698 x 240 x 236 mm	15 Boxes