




VALID FROM APRIL 1 TO JUNE 30, 2021

PURCHASE A KBF-S SOLID.LINE HUMIDITY TEST CHAMBER AND RECEIVE A FREE COMFORT WATER SUPPLY KIT

FEATURES:

- Temperature range: 0 °C to 70 °C
- Humidity range: 20% to 80% RH
- APT.line™ preheating chamber technology
- Humidity regulation with capacitive humidity sensor and steam humidification
- LCD to display temperature and humidity along with additional information and alarms
- Internal data logger, measured values can be read out in open format via USB
- Unit self-test for comprehensive status analysis

Series	115	240	720	1020
KBF-S	•	•	•	•



Model KBF-S 720



Comfort water supply
Flexible positioning of the water canister with magnetic holder, hose connection to the canister with practical quick coupling and water stop function.

