


LABRepCo Brand Laboratory Seating

Work hard,
sit smart.



- Industry leading 15-year warranty
- Laboratory Grade Vinyl and Polyurethane
- Entire range of Labrepco seating is available for QUICK SHIP
- Go to Labrepco.com/build-a-chair to build the perfect chair for your setting

LABRepCo
www.labrepco.com

Go from Boring Lab to Bright Lab


Trailblazer™

Trailblazer vinyl is inherently antibacterial due to its closed-cell nature.



Supernova™

Supernova textured vinyl has an additive for additional protection against bacteria.



A multitude of vinyl color options are available to help you compliment any laboratory setting. We can utilize your materials for further customizing.

For more information on color options, please contact your Labrepco Representative.

American Made - Laboratory Grade



Choice of caster or glides.

Choosing the Right Labrepco Brand Chair

Our seating is designed using ergonomic principles to create a full range of chairs and stools to be used in various environments. Here are a few chair selection tips to help you build the perfect laboratory chair.

1

Determine your work surface height:

Surface Height	Chair Height	
29" - 31"	16" - 21"	Desk (Low-Form)
32" - 36"	22" - 28"	Bench (Medium-Form)
37" - 42"	26" - 35"	High Bench (High-Form)

2

Determine the environment in which the chair will be used and pick the appropriate material:

- **Standard laboratory** (vinyl, polyurethane, wood)
- **Office** (vinyl, fabric, wood)
- **Clean Room and/or ESD** (vinyl, polyurethane)

3

Determine if a Foot Ring is needed. For work surfaces at or above 32" we recommend the use of an adjustable Foot Ring to provide for better circulation in the lower extremities.



Please go to Labrepco.com/build-a-chair

LabRepCo Laboratory Chairs and Stools



**Bench Height Vinyl Chair
Chrome Footring
LRC49993**

Fully adjustable and available with a variety of ergonomic options.




**High Bench Height
Polyurethane Chair
LRC44670**

This soil-resistant polyurethane chair features a large supportive back.



**Bench Height Tractor Stool
LRC45146**

New design may just be our most comfortable polyurethane stool yet.



Go to labrepco.com/build-a-chair to build your perfect chair in 30 seconds or less!



**High Bench Height Polyurethane
Chair
LRC44832**

Robust polyurethane chair with adjustable chrome foot ring. Soil, chemical, and puncture resistant.



**Vinyl Chair, Desk Height
LRC46006**

Supernova-Borscht

Our ergonomic chair is fully adjustable with a variety of ergonomic options.



**Bench Height Vinyl Stool with
Back
LRC51002**

Supernova - Marine Blue

These upholstered stools feature 3-inches of molded polyurethane foam seats, and an adjustable back.



WARNING: These products can expose you to chemicals, including Diisononyl Phthalate (DINP) or Indeno[1,2,3-CD]pyrene, which are known to the State of California to cause cancer. For more information go to www.p65warning.com.

Controlled Environment

Class 10 Clean Room



Class 10 Clean Room Chair, Medium Bench Height
LRC40388

Supernova - Sterling
Supernova Sterling vinyl with chrome footring and glides for stationary use.

Class 100 Clean Room Chairs



Desk Height Class 100 Polyurethane Clean Room Chair
LRC43442

ISO Class 5 (Class 100) Polyurethane Clean Room Chair features adjustable arms.



High Bench Height ISO Class 5 (Class 100) vinyl cleanroom chair
LRC43811

Supernova - Carbon
Our ergonomic chair features an adjustable chrome foot ring.

MADE ★ IN THE ★ USA

ESD Chairs



Fabric ESD Chair, Desk Height
LRC40455

Our black Fabric ESD Chair features proper grounding, conductive, and dissipating materials without the nuisance of a drag chain.



Polyurethane ESD Chair, Medium Bench Height
LRC40604

These Polyurethane ESD chairs combine the soil-, chemical-, and puncture-resistance of our polyurethane line with the grounding, conductive, and dissipating materials needed for sensitive electronic applications.

ESD/ Clean Room Chair



Bench Height Class 100 Clean Room/ESD Chair
LRC40436

These ISO 5 (Class 100) Clean Room/ESD Chairs provide both clean room and static dissipating protection.



WARNING: These products can expose you to chemicals, including Diisononyl Phthalate (DINP) or Indeno[1,2,3-CD] pyrene, which are known to the State of California to cause cancer. For more information go to www.p65warning.com.