

Entris[®] II Advanced Line

The new Entris[®] II balance is always the right choice

Our Entris[®] II Advanced Line gives additional value to the Entris[®] II portfolio with over 38 models offering benefits such as real-time level support, integrated protection systems, Cal Audit Trail, a graphic touch display and 13 built-in application programs.

Features & Benefits

- Simplifies balance leveling with built-in real-time level support
- Increases reliability of weighing results using integrated protection systems
- Controlled access to balance settings with user management
- Gapless documentation thanks to CalAuditTrail



Get a free Printer
and Cable

Buy a select Sartorius Entris® II Advanced Balance and receive a free YDP20- OCE printer, Power Supply (YEPS01-PS4) and cable (YCC-D09MF)



The printer meets the quality standards from the chemical and pharmaceutical industries, which are introduced in the research sector under the term GLP and in production under the term GMP. In addition to individual weighing data, the printer can also be used for statistical evaluation of weighing data. For this purpose, all weight values received from the balance are stored.

Special offer valid only for the NA until May 31, 2021. It cannot be used in combination with other special offers.

Entris® II Advanced Balance

Model No.	Capacity (g)	Readability (mg)
Entris® II Advanced Line Analytical Balance with External Calibration		
BCA64-1S	60	0.1
BCA124-1S	120	0.1
BCA224-1S	220	0.1
Entris® II Advanced Line Analytical Balance with Internal Calibration		
BCA64I-1S	60	0.1
BCA124I-1S	120	0.1
BCA224I-1S	220	0.1
BCA324I-1S	320	0.1
Entris® II Advanced Line Precision Balance with External Calibration		
BCA423-1S	420	1
BCA623-1S	620	1
BCA3202-1S	3200	10
BCA4202-1S	4200	10
BCA6202-1S	6200	10
BCA8201-1S	8200	100
BCA10201-1S	10200	100
BCA12201-1S	12200	100

Model No.	Capacity (g)	Readability (mg)
Entris® II Advanced Line Precision Balance with Internal Calibration		
BCA223I-1S	220	1
BCA423I-1S	420	1
BCA323I-1S	320	1
BCA623I-1S	620	1
BCA1203I-1S	1200	1
BCA1202I-1S	1200	10
BCA2202I-1S	2200	10
BCA3202I-1S	3200	10
BCA4202I-1S	4200	10
BCA6202I-1S	6200	10
BCA8201I-1S	8200	100
BCA10201I-1S	10200	100
BCA12201I-1S	12200	100

[Click here to view products](#)

