

**Vinyl Stool without Back
- Medium Bench Height
with Chrome Foot Ring
& Casters in Orange
Kist Trailblazer™ Vinyl**

SKU: LRC50875



Details

Vinyl Medium Bench Height Stool is adjustable from 19.75-27.25 inches from the top of the seat to the floor and features ergonomic adjustments, including 360-degree, all-around control, chrome adjustable foot ring, and black, nylon-reinforced fiberglass base with dual-wheel, self-braking, non-marring casters suitable for hard floors in antibacterial/antifungal Orange Kist Trailblazer™ Vinyl.

Specifications

Adjustment	Pneumatic Adjustment: 19.75-27.75 in
Seat	Diameter 17in, Thickness: 3.5in
Footring	18 in diameter round tube chrome foot ring with knob lock
Color	Orange Kist
Seat Material	Trailblazer™ Vinyl: Highly durable molded polyurethane foam cushioning covered with easily cleaned antibacterial/antifungal vinyl. Our upholstered vinyl chairs are protected by an industrial strength stain resistant finish.
Feet	Black dual-wheel, self-braking, non-marring casters (suitable for both hard and

carpeted floors)

Base	Nylon Reinforced Fiberglass Base
Base Width	24 in
Base Depth	24 in
Weight (lbs)	28
Weight Capacity	300 lbs
Warranty	Warranty: 15 year parts, lifetime on pneumatic lift, 5 year on upholstery

