



**Vinyl Stool with Back  
Chrome - Medium Bench  
Height with Chrome  
Foot Ring & Casters in  
Orange Kist  
Trailblazer™ Vinyl**

**SKU: LRC51001**

---

---

## **Details**

Vinyl Medium Bench Height Stool with Back in Chrome is adjustable from 22-29.75 inches from the top of the seat to the floor and features ergonomic adjustments including seat height and locking back tilt, chrome adjustable foot ring, and cast-aluminum base with chrome-hooded, dual-wheel, self-braking, non-marring casters suitable for hard floors in antibacterial/antifungal Orange Kist Trailblazer™ Vinyl.

---

## **Specifications**

<b>Adjustment</b>	Pneumatic Adjustment: 22-29.75 in
<b>Seat</b>	Diameter 17in, Thickness: 3.5in
<b>Back</b>	Width: 16.5in, Height: 11.5in, Thickness: 2.5in
<b>Footring</b>	20 in diameter round tube chrome foot ring with internal lock
<b>Color</b>	Orange Kist
<b>Seat Material</b>	Trailblazer™ Vinyl: Highly durable molded polyurethane foam cushioning covered with easily cleaned antibacterial/antifungal vinyl. Our upholstered vinyl chairs are

protected by an industrial strength stain resistant finish.

<b>Feet</b>	Chrome-hooded dual-wheel, self-braking, non-marring casters (suitable for both hard and carpeted floors)
<b>Base</b>	Polished Cast Aluminum Base
<b>Base Width</b>	24 in
<b>Base Depth</b>	24 in
<b>Weight (lbs)</b>	34
<b>Weight Capacity</b>	300 lbs
<b>Warranty</b>	Warranty: 15 year parts, lifetime on pneumatic lift, 5 year on upholstery

