



**Polyurethane Chair -
Medium Bench Height
with Chrome Foot Ring
& Casters in Black
Polyurethane**

SKU: LRC45016

Details

Black Polyurethane Medium Bench Height Chair is adjustable from 19.25-26.75 inches from the top of the seat to the floor and features ergonomic adjustments including seat height and locking back tilt, chrome adjustable foot ring, and black, nylon-reinforced fiberglass base with dual-wheel, self-braking, non-marring casters suitable for hard floors in chemical, soil, and puncture resistant polyurethane.

Specifications

Adjustment	Pneumatic Adjustment: 19.25-26.75 in
Seat	Width: 16in, Depth: 17in, Thickness: 1.5in
Back	Width: 15.5in, Height: 10.25in, Thickness: 1 in
Footring	20 in diameter round tube chrome foot ring with internal lock
Color	Black
Seat Material	Rugged polyurethane that is self-skinning, chemical- and puncture-resistant, easily cleaned, and soil resistant

Feet	Black dual-wheel, self-braking, non-marring casters (suitable for both hard and carpeted floors)
Base	Nylon Reinforced Fiberglass Base
Base Width	24 in
Base Depth	24 in
Weight (lbs)	34
Weight Capacity	300 lbs
Warranty	Warranty: 15 year parts, lifetime on pneumatic lift

