



MW200 Series



MW400 Series



MW600 Series



MW700 Series

MW Series Utility Carts — 12.01

- Consist of Super Erecta wire shelves with plastic split sleeves, two one-piece handles of matching finish, and designated casters with donut bumpers.
- Two- and three-tier models available.
- 375 lbs. (170kg) weight capacity per cart. MW carts use light-duty resilient casters.
- 18" (457mm) wide carts have 4" (102mm) casters. Cart is 38" (965mm) high.
- 21" and 24" (530 and 610mm) wide carts have 5" (127mm) casters. Cart is 39" (990mm) high.

Utility Carts with Stainless Steel Solid Shelves

Width/Length (in.) (mm)		Approx. Pkd. Wt. (lbs.) (kg)		Cat. No. 2-Shelf	Approx. Pkd. Wt. (lbs.) (kg)		Cat. No. 3-Shelf
18x24	457x610	40	18.1	LAB301WM	48	21.8	LAB302WM
18x30	457x760	45	20.4	LAB401WM	53	24.0	LAB402WM
18x36	457x914	48	21.8	LAB501WM	59	26.8	LAB502WM
21x36	530x914	54	24.5	LAB601WM	74	33.6	LAB602WM
24x36	610x914	60	27.2	LAB801WM	66	29.9	LAB802WM

Utility Carts with 1 Solid, 2 Wire Shelves

Width/Length (in.) (mm)		Approx. Pkd. Wt. (lbs.) (kg)		Shelf Material	Handles	Cat. No.
18x24	457x610	43	19.5	Stainless Steel, Chrome	Chrome	LAB104WM
18x30	457x760	47	21.3	Stainless Steel, Chrome	Chrome	LAB204WM
18x36	457x914	52	23.6	Stainless Steel, Chrome	Chrome	LAB304WM
21x36	530x914	59	26.8	Stainless Steel, Chrome	Chrome	LAB404WM
24x36	610x914	65	29.5	Stainless Steel, Chrome	Chrome	LAB604WM

Utility Carts with 2 Wire Shelves

Width/Length (in.) (mm)		Approx. Pkd. Wt. (lbs.) (kg)		Shelf Material	Handles	Cat. No.
18x24	457x610	34	15.4	Chrome	Chrome	LAB106WM
18x24	457x610	34	15.4	Stainless Steel	Stainless Steel	LAB206WM
18x30	457x760	37	16.8	Chrome	Chrome	LAB306WM
18x30	457x760	37	16.8	Stainless Steel	Stainless Steel	LAB406WM
18x36	457x914	40	18.1	Chrome	Chrome	LAB506WM
18x36	457x914	39	17.7	Stainless Steel	Stainless Steel	LAB606WM
21x36	530x914	44	20.0	Chrome	Chrome	LAB706WM
21x36	530x914	44	20.0	Stainless Steel	Stainless Steel	LAB806WM
24x36	610x914	47	21.3	Chrome	Chrome	LAB116WM
24x36	610x914	46	20.9	Stainless Steel	Stainless Steel	LAB216WM

Utility Carts with 3 Wire Shelves

Width/Length (in.) (mm)		Approx. Pkd. Wt. (lbs.) (kg)		Shelf Material	Handles	Cat. No.
18x24	457x610	39	17.6	Chrome	Chrome	LAB107WM
18x24	457x610	38	17.2	Stainless Steel	Stainless Steel	LAB207WM
18x30	457x760	43	19.5	Chrome	Chrome	LAB307WM
18x30	457x760	42	19.1	Stainless Steel	Stainless Steel	LAB407WM
18x36	457x914	47	21.3	Chrome	Chrome	LAB507WM
18x36	457x914	45	20.4	Stainless Steel	Stainless Steel	LAB607WM
21x36	530x914	50	22.7	Chrome	Chrome	LAB707WM
21x36	530x914	50	22.7	Stainless Steel	Stainless Steel	LAB807WM
24x36	610x914	60	27.7	Chrome	Chrome	LAB117WM
24x36	610x914	57	25.9	Stainless Steel	Stainless Steel	LAB217WM