


MVE 1536PD Dry Shipper




The MVE 1536PD Dry Shipper is specifically designed to provide you with peace of mind when shipping and/or relocating samples.

Features:


- Liquid Nitrogen absorbent can be charged for dry shipment
- Hold time up to 20 days when properly charged
- Rigid steel frame can be moved with pallet jack/forklift
- Sized to fit in a wide-body aircraft
- Inventory tray locks to prevent rotation during transit
- Locking lid
- Includes data logger to monitor and record sample temperature during shipment or relocation
- Inner vessel may be decontaminated after shipping if necessary.



Top View



Cabosil Under Tray



Locking Cabinet




Data Logger




Rack Layout

MVE 1536PD Dry Shipper



24, 13-2 Large Racks
16, 13-2 Mini Racks
36,400 Vials



MVE 1536PD

Maximum Storage Capacity

1.2 & 2 mL Vials (Internally Threaded)	36400
Number of Racks 100 cell boxes	24
Number of Racks 25 cell boxes	16
Number of Stages per Rack	13

Performance

Absorbed LN2 Capacity L	133
--------------------------------	-----

Unit Dimensions

Neck Opening in. (mm)	17.5 (444)
Usable Internal Height in. (mm)	30 (762)
Inner Diameter in. (mm)	38.5 (978)
Overall Height, including Frame in. (mm)	64 (1625)
Liftover Height in. (mm)	42.1 (1070)
Footprint in. (mm)	52x52 (1321x1321)
Step Height in. (mm)	20.5 (521)
Depth of Fold Down Step in. (mm)	10.8 (274)
Weight Empty lb. (kg)	1600 (726)
Weight Full lb. (kg)	2044 (927)

TWO Year Standard Warranty • THREE Year Vacuum Warranty

Chart Industries, Inc.
MVE, BioMedical Group
2200 Airport Industrial Drive, Ste 500
Ball Ground, GA 30107
Ph 770-721-7759 • Toll Free 1-800-482-2473
Fax 770-721-7758



Chart BioMedical, Ltd.
Unit 2, Maxdata Centre
Downmill Road, Bracknell
Berks RG12 1QS, United Kingdom
Ph +44(0) 1344 403100 • Fax +44(0) 1344 429224