


Aluminum Canisters

For Blood Bags, Bone Marrow, Stem Cell, Cord Blood & Plasma Bags

Canisters constructed of lightweight aluminum are made for most popular brands of bags including Baxter Fenwal, Chartered, Gambro and Stericon. Canisters have locking handle and thumb hole



For BAXTER/FENWAL BAGS

Canister Part Number	Blood Bag Number	Description	Dimension $L_c \times D_c \times H_c$	Frames	Platform Dividers
C-4R9951	4R-9951	BAXTER/FENWAL Blood Bag Canister	$6\frac{7}{16} \times 3\frac{9}{16} \times \frac{7}{16}$ inch 163 x 90 x 11 mm	F-C-4R9951-5 F-C-4R9951-7 F-C-4R9951-8	P20-4R9951 P24-4R9951 P31-4R9951 P39-4R9951
C-4R9953	4R-9953 4R-5461	BAXTER/FENWAL Blood Bag Canister	$7\frac{15}{16} \times 5\frac{1}{2} \times \frac{1}{2}$ inch 201 x 140 x 12 mm	F-C-4R9953-3 F-C-4R9953-4 F-C-4R9953-5	P20-4R9953 P24-4R9953 P31-4R9953 P39-4R9953
C-4R9955	4R-9955 4R-5462	BAXTER/FENWAL Blood Bag Canister	$9\frac{7}{8} \times 5\frac{1}{2} \times \frac{1}{2}$ inch 251 x 140 x 12 mm	F-C-4R9955-3 F-C-4R9955-4 F-C-4R9955-5	P20-4R9955 P24-4R9955 P31-4R9955 P39-4R9955
C-4R9957	4R-9957 4R-5463	BAXTER/FENWAL Blood Bag Canister	$10\frac{7}{8} \times 5\frac{1}{2} \times \frac{1}{2}$ inch 277 x 140 x 12 mm	F-C-4R9957-3 F-C-4R9957-4 F-C-4R9957-5	P20-4R9957 P24-4R9957 P31-4R9957 P39-4R9957
C-4R9959	4R-9959 4R-5464	BAXTER/FENWAL Blood Bag Canister	$13\frac{3}{4} \times 5\frac{1}{2} \times \frac{1}{2}$ inch 340 x 140 x 12 mm	F-C-4R9959-3 F-C-4R9959-4 F-C-4R9959-5	P20-4R9959 P24-4R9959 P31-4R9959 P39-4R9959
C-4R2111	4R-2111	BAXTER/FENWAL Blood Bag Canister	$8\frac{1}{2} \times 6\frac{1}{4} \times \frac{13}{16}$ inch 216 x 159 x 20 mm	F-C-4R2111-3 F-C-4R2111-4	P20-4R2111 P24-4R2111 P31-4R2111 P39-4R2111
C-4R2110	4R-2110	BAXTER/FENWAL Blood Bag Canister	$10 \times 8\frac{1}{16} \times \frac{13}{16}$ inch 254 x 204 x 20 mm	F-C-4R2110-3 F-C-4R2110-4	P20-4R2110 P24-4R2110 P31-4R2110 P39-4R2110

For STERICON BAGS

Part Number	Blood Bag Number	Description	Dimension $L_c \times D_c \times H_c$	Frames	Platform Dividers
C-91D	RCM-91D	STERICON Blood Bag Canister	$11 \times 5\frac{1}{2} \times \frac{1}{2}$ inch 280 x 140 x 12 mm	F-C-91D-3 F-C-91D-4 F-C-91D-5	P20-91D P24-91D P31-91D P39-91D

(Continued on next page)

Aluminum Canisters

For Bone Marrow, Stem Cell, Cord Blood & Plasma Bags

For CHARTERED BAGS

Part Number	Blood Bag Number	Description	Dimension $L_c \times D_c \times H_c$	Frames	Platform Dividers
C-7120	7120-2	CHARTERED Blood Bag Canister	$5\frac{1}{8} \times 4\frac{3}{16} \times 7\frac{7}{16}$ inch 130 x 107 x 11 mm	F-C-7120-4 F-C-7120-5 F-C-7120-6	P20-7120 P24-7120 P31-7120 P39-7120
C-7005	7005-2 7303-2 7039-2	CHARTERED Blood Bag Canister	$9\frac{1}{8} \times 7\frac{3}{8} \times 7\frac{7}{16}$ inch 232 x 188 x 11 mm	F-C-7005-2 F-C-7005-3	P20-7005 P24-7005 P31-7005 P39-7005
C-7400	7400-2 7403-2 7410-2 7025-2	CHARTERED Blood Bag Canister	$9 \times 9 \times 7\frac{7}{16}$ inch 229 x 229 x 11 mm	F-C-7400-2	P20-7400 P24-7400 P31-7400 P39-7400

For GAMBRO BAGS

Part Number	Blood Bag Number	Description	Dimension $L_c \times D_c \times H_c$	Frames	Platform Dividers
C-DF200	DF-200	GAMBRO Blood Bag Canister	$6\frac{7}{16} \times 4\frac{3}{16} \times 1\frac{13}{16}$ inch 163 x 106 x 21 mm	F-C-DF200-4 F-C-DF200-6	P20-DF200 P24-DF200 P31-DF200 P39-DF200
C-DF700	DF-700	GAMBRO Blood Bag Canister	$12 \times 6\frac{3}{16} \times 1\frac{13}{16}$ inch 304.5 x 157 x 21 mm	F-C-DF700-4 F-C-DF700-6	P20-DF700 P24-DF700 P31-DF700 P39-DF700

For MEDSEP BAGS


Part Number	Blood Bag Number	Description	Dimension $L_c \times D_c \times H_c$	Frames/Rack	Platform Dividers
CBBC-V1	791-05 791-05U 8-0346-0	Cord Blood Bag Canister	$3\frac{5}{16} \times 3\frac{3}{8} \times 7\frac{7}{16}$ inch 99.5 x 92 x 11 mm	F-C-CBBC-V1-6 F-C-CBBC-V1-7 F-C-CBBC-V1-8	P20-CBBC-V1 P24-CBBC-V1 P31-CBBC-V1 P39-CBBC-V1
			$3\frac{5}{16} \times 3\frac{3}{8} \times 7\frac{7}{16}$ inch 99.5 x 92 x 11 mm	CF-8-3-TN CF-9-3-TN CF-10-3-TN	

(Continued from last page)

Stainless Steel Canister Frames

Cannister frames are constructed with corrosion resistant stainless steel and include handles.

For BAXTER/FENWAL BAGS




Frame Part Number	Blood Bag Number	Description	Dimensions (D x H x L)	Canister	Platform Dividers
F-C-4R9951-5	4R-9951	Stainless Steel Frame for 5 BAXTER/ FENWALC-4R9951 Aluminum Canisters.	5 ¹¹ / ₁₆ x 18 ³ / ₄ x ⁹ / ₁₆ inch 145 x 476 x 15 mm	C-4R9951	P20-4R9951 P24-4R9951 P31-4R9951 P39-4R9951
F-C-4R9951-7		Stainless Steel Frame for 7 BAXTER/ FENWALC-4R9951 Aluminum Canisters.	5 ¹¹ / ₁₆ x 26 ¹ / ₁₆ x ⁹ / ₁₆ inch 145 x 662 x 15 mm		
F-C-4R9951-8		Stainless Steel Frame for 8 BAXTER/ FENWALC-4R9951 Aluminum Canisters.	5 ¹¹ / ₁₆ x 29 ³ / ₄ x ⁹ / ₁₆ inch 145 x 755 x 15 mm		
F-C-4R9953-3	4R-9953 4R-5461	Stainless Steel Frame for 3 BAXTER/ FENWALC-4R9953 Aluminum Canisters.	7 ⁵ / ₁₆ x 17 ⁵ / ₁₆ x ⁹ / ₁₆ inch 186 x 440 x 15 mm	C-4R9953	P20-4R9953 P24-4R9953 P31-4R9953 P39-4R9953
F-C-4R9953-4		Stainless Steel Frame for 4 BAXTER/ FENWALC-4R9953 Aluminum Canisters.	7 ⁵ / ₁₆ x 22 ¹⁵ / ₁₆ x ⁹ / ₁₆ inch 186 x 583 x 15 mm		
F-C-4R9953-5		Stainless Steel Frame for 5 BAXTER/ FENWALC-4R9953 Aluminum Canisters.	7 ⁵ / ₁₆ x 28 ⁹ / ₁₆ x ⁹ / ₁₆ inch 186 x 726 x 15 mm		
F-C-4R9955-3	4R-9955 4R-5462	Stainless Steel Frame for 3 BAXTER/ FENWALC-4R9955 Aluminum Canisters.	9 ³ / ₁₆ x 17 ⁵ / ₁₆ x ⁹ / ₁₆ inch 234 x 440 x 15 mm	C-4R9955	P20-4R9955 P24-4R9955 P31-4R9955 P39-4R9955
F-C-4R9955-4		Stainless Steel Frame for 4 BAXTER/ FENWALC-4R9955 Aluminum Canisters.	9 ³ / ₁₆ x 22 ¹⁵ / ₁₆ x ⁹ / ₁₆ inch 234 x 583 x 15 mm		
F-C-4R9955-5		Stainless Steel Frame for 5 BAXTER/ FENWALC-4R9955 Aluminum Canisters.	9 ³ / ₁₆ x 28 ⁹ / ₁₆ x ⁹ / ₁₆ inch 234 x 726 x 15 mm		
F-C-4R9957-3	4R-9957 4R-5463	Stainless Steel Frame for 3 BAXTER/ FENWALC-4R9957 Aluminum Canisters.	10 ⁵ / ₁₆ x 17 ⁵ / ₁₆ x ⁹ / ₁₆ inch 262 x 440 x 15 mm	C-4R9957	P20-4R9957 P24-4R9957 P31-4R9957 P39-4R9957
F-C-4R9957-4		Stainless Steel Frame for 4 BAXTER/ FENWALC-4R9957 Aluminum Canisters.	10 ⁵ / ₁₆ x 22 ¹⁵ / ₁₆ x ⁹ / ₁₆ inch 262 x 583 x 15 mm		
F-C-4R9957-5		Stainless Steel Frame for 5 BAXTER/ FENWALC-4R9957 Aluminum Canisters.	10 ⁵ / ₁₆ x 28 ⁹ / ₁₆ x ⁹ / ₁₆ inch 262 x 726 x 15 mm		
F-C-4R9959-3	4R-9959 4R-5464	Stainless Steel Frame for 3 BAXTER/ FENWALC-4R9959 Aluminum Canisters.	12 ¹¹ / ₁₆ x 17 ⁵ / ₁₆ x ⁹ / ₁₆ inch 322 x 440 x 15 mm	C-4R9959	P20-4R9959 P24-4R9959 P31-4R9959 P39-4R9959
F-C-4R9959-4		Stainless Steel Frame for 4 BAXTER/ FENWALC-4R9959 Aluminum Canisters.	12 ¹¹ / ₁₆ x 22 ¹⁵ / ₁₆ x ⁹ / ₁₆ inch 322 x 583 x 15 mm		
F-C-4R9959-5		Stainless Steel Frame for 5 BAXTER/ FENWALC-4R9959 Aluminum Canisters.	12 ¹¹ / ₁₆ x 28 ⁹ / ₁₆ x ⁹ / ₁₆ inch 322 x 726 x 15 mm		
F-C-4R2111-3	4R-2111	Stainless Steel Frame for 3 BAXTER/ FENWALC-4R2111 Aluminum Canisters.	7 ¹¹ / ₁₆ x 20 x ¹⁵ / ₁₆ inch 196 x 508 x 24 mm	C-4R2111	P20-4R2111 P24-4R2111 P31-4R2111 P39-4R2111
F-C-4R2111-4		Stainless Steel Frame for 4 BAXTER/ FENWALC-4R2111 Aluminum Canisters.	7 ¹¹ / ₁₆ x 26 ⁹ / ₁₆ x ¹⁵ / ₁₆ inch 196 x 674 x 24 mm		
F-C-4R2110-3	4R-2110	Stainless Steel Frame for 3 BAXTER/ FENWALC-4R2110 Aluminum Canisters.	9 ³ / ₁₆ x 25 ¹ / ₄ x ¹⁵ / ₁₆ inch 234 x 641 x 24 mm	C-4R2110	P20-4R2110 P24-4R2110 P31-4R2110 P39-4R2110
F-C-4R2110-4		Stainless Steel Frame for 4 BAXTER/ FENWALC-4R2110 Aluminum Canisters.	9 ³ / ₁₆ x 33 ¹ / ₂ x ¹⁵ / ₁₆ inch 234 x 851 x 24 mm		

(Continued on next page)

Stainless Steel Canister Frames

Canister frames are constructed with corrosion resistant stainless steel and include handles.



For CHARTERED BAGS


Part Number	Blood Bag Number	Description	Dimensions (D x H x L)	Canister	Platform Dividers
F-C-7120-4	7120-2	Stainless Steel Frame for 4 CHARTERED C-7120 Aluminum Canisters.	4 ⁵ / ₁₆ x 17 ³ / ₄ x ⁹ / ₁₆ inch 110 x 451 x 15 mm	C-7120	P20-7120 P24-7120 P31-7120 P39-7120
F-C-7120-5		Stainless Steel Frame for 5 CHARTERED C-7120 Aluminum Canisters.	4 ⁵ / ₁₆ x 22 ¹ / ₁₆ x ⁹ / ₁₆ inch 110 x 561 x 15 mm		
F-C-7120-6		Stainless Steel Frame for 6 CHARTERED C-7120 Aluminum Canisters.	4 ⁵ / ₁₆ x 26 ⁷ / ₁₆ x ⁹ / ₁₆ inch 110 x 671 x 15 mm		
F-C-7005-2	7005-2	Stainless Steel Frame for 2 CHARTERED C-7005 Aluminum Canisters.	8 ⁵ / ₈ x 15 ¹ / ₂ x ⁹ / ₁₆ inch 212 x 393 x 15 mm	C-7005	P20-7005 P24-7005 P31-7005 P39-7005
F-C-7005-3	7303-2 7039-2	Stainless Steel Frame for 3 CHARTERED C-7005 Aluminum Canisters.	8 ⁵ / ₈ x 22 ¹⁵ / ₁₆ x ⁹ / ₁₆ inch 212 x 582 x 15 mm		
F-C-7400-2	7400-2 7403-2	Stainless Steel Frame for 2 CHARTERED C-7400 Aluminum Canisters.	8 ¹ / ₄ x 18 ⁵ / ₁₆ x ⁹ / ₁₆ inch 209 x 475 x 15 mm	C-7400	P20-7400 P24-7400 P31-7400 P39-7400
F-C-7400-3	7410-2 7025-2	Stainless Steel Frame for 3 CHARTERED C-7400 Aluminum Canisters.	8 ¹ / ₄ x 27 ¹³ / ₁₆ x ⁹ / ₁₆ inch 209 x 707 x 15 mm		

For GAMBRO BAGS

Part Number	Blood Bag Number	Description	Dimensions (D x H x L)	Canister	Platform Dividers
F-C-DF200-4	DF200	Stainless Steel Frame for 4 GAMBRO C-DF200 Aluminum Canisters.	5 ⁷ / ₈ x 18 x 1 ¹ / ₈ inch 150 x 457 x 28 mm	C-DF200	P20-DF200 P24-DF200 P31-DF200 P39-DF200
F-C-DF200-6		Stainless Steel Frame for 6 GAMBRO C-DF200 Aluminum Canisters.	5 ⁷ / ₈ x 27 ⁷ / ₁₆ x 1 ¹ / ₈ inch 150 x 690 x 28 mm		
F-C-DF700-4	DF700	Stainless Steel Frame for 4 GAMBRO C-DF700 Aluminum Canisters.	11 ³ / ₁₆ x 25 ⁵ / ₈ x ¹⁵ / ₁₆ inch 284.5 x 651 x 24 mm	C-DF700	P20-DF700 P24-DF700 P31-DF700 P39-DF700
F-C-DF700-6		Stainless Steel Frame for 6 GAMBRO C-DF700 Aluminum Canisters.	11 ³ / ₁₆ x 38 ¹ / ₄ x ¹⁵ / ₁₆ inch 284.5 x 971 x 24 mm		

Stainless Steel Canister Frames

Cannister frames are constructed with corrosion resistant stainless steel and include handles.



For STERICON BAGS

Part Number	Blood Bag Number	Description	Dimensions (D x H x L)	Canister	Platform Dividers
F-C-91D-3	RCM-91D	Stainless Steel Frame for 3 STERICON C-91D Aluminum Canisters.	10 ³ / ₁₆ x 17 ⁵ / ₁₆ x ⁹ / ₁₆ inch 259 x 440 x 15 mm	C-91D	P20-91D P24-91D P31-91D P39-91D
F-C-91D-4		Stainless Steel Frame for 4 STERICON C-91D Aluminum Canisters.	10 ³ / ₁₆ x 22 ¹⁵ / ₁₆ x ⁹ / ₁₆ inch 259 x 583 x 15 mm		
F-C-91D-5		Stainless Steel Frame for 5 STERICON C-91D Aluminum Canisters.	10 ³ / ₁₆ x 28 ⁹ / ₁₆ x ⁹ / ₁₆ inch 259 x 726 x 15 mm		


For MEDSEP BAGS

Part Number	Blood Bag Number	Description	Dimensions (D x H x L)	Canister	Platform Dividers
F-C-CBBC-V1-6	791-OS 791-OSU 8-0346-O	Stainless Steel Frame for 6 Med Sep CBBC-V1 Aluminum Cord Blood Bag Canisters.	3 ¹ / ₄ x 23 ¹ / ₈ x ⁹ / ₁₆ inch 82 x 587 x 15 mm	CBBC-V1	P20-CBBC-V P24-CBBC-V P31-CBBC-V P39-CBBC-V
F-C-CBBC-V1-7		Stainless Steel Frame for 7 Med Sep CBBC-V1 Aluminum Cord Blood Bag Canisters.	3 ¹ / ₄ x 26 ⁷ / ₈ x ⁹ / ₁₆ inch 82 x 683 x 15 mm		
F-C-CBBC-V1-8		Stainless Steel Frame for 8 Med Sep CBBC-V1 Aluminum Cord Blood Bag Canisters.	3 ¹ / ₄ x 30 ¹¹ / ₁₆ x ⁹ / ₁₆ inch 82 x 780 x 15 mm		
CF-8-3-T		Stainless Steel Vertical Rack for 56 CBBC-V1 Canisters.	3 ⁹ / ₁₁ x 25 x 4 ³ / ₈ inch 97 x 635 x 111 mm		
CF-9-3-T		Stainless Steel Vertical Rack for 63 CBBC-V1 Canisters.	3 ⁹ / ₁₁ x 28 ¹ / ₁₆ x 4 ³ / ₈ inch 97 x 713 x 111 mm		
CF-10-3-T		Stainless Steel Vertical Rack for 70 CBBC-V1 Canisters.	3 ⁹ / ₁₁ x 31 ¹ / ₈ x 4 ³ / ₈ inch 97 x 791 x 111 mm		

(Continued on next page)

Aluminum Platform Dividers for Liquid Nitrogen Storage Tanks

For BAXTER/FENWAL BAGS




Part Number	Blood Bag Number(S)	Description	Dimensions (H x φ)	Canister	Frames
P20-4R9951	4R-9951	20" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F-C-4R9951.	20 x 20 inch 508 x 508 mm	C-4R9951	F-C-4R9951-5 F-C-4R9951-7 F-C-4R9951-8
P24-4R9951		24" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F-C-4R9951.	20 x 24 inch 508 x 610 mm		
P31-4R9951		31" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F-C-4R9951.	20 x 31 inch 508 x 787 mm		
P39-4R9951		39" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F-C-4R9951.	20 x 39 inch 508 x 991 mm		
P20-4R9953	4R-9953 4R-5461	20" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F-C-4R9953.	20 x 20 inch 508 x 508 mm	C-4R9953	F-C-4R9953-3 F-C-4R9953-4 F-C-4R9953-5
P24-4R9953		24" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F-C-4R9953.	20 x 24 inch 508 x 610 mm		
P31-4R9953		31" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F-C-4R9953.	20 x 31 inch 508 x 787 mm		
P39-4R9953		39" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F-C-4R9953.	20 x 39 inch 508 x 991 mm		
P20-4R9955	4R-9955 4R-5462	20" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F-C-4R9955.	20 x 20 inch 508 x 508 mm	C-4R9955	F-C-4R9955-3 F-C-4R9955-4 F-C-4R9955-5
P24-4R9955		24" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F-C-4R9955.	20 x 24 inch 508 x 610 mm		
P31-4R9955		31" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F-C-4R9955.	20 x 31 inch 508 x 787 mm		
P39-4R9955		39" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F-C-4R9955.	20 x 39 inch 508 x 991 mm		
P20-4R9957	4R-9957 4R-5463	20" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F-C-4R9957.	20 x 20 inch 508 x 508 mm	C-4R9957	F-C-4R9957-3 F-C-4R9957-4 F-C-4R9957-5
P24-4R9957		24" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F-C-4R9957.	20 x 24 inch 508 x 610 mm		
P31-4R9957		31" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F-C-4R9957.	20 x 31 inch 508 x 787 mm		
P39-4R9957		39" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F-C-4R9957.	20 x 39 inch 508 x 991 mm		
P20-4R9959	4R-9959 4R-5464	20" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F-C-4R9959.	20 x 20 inch 508 x 508 mm	C-4R9959	F-C-4R9959-3 F-C-4R9959-4 F-C-4R9959-5
P24-4R9959		24" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F-C-4R9959.	20 x 24 inch 508 x 610 mm		
P31-4R9959		31" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F-C-4R9959.	20 x 31 inch 508 x 787 mm		
P39-4R9959		39" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F-C-4R9959.	20 x 39 inch 508 x 991 mm		

(Continued on next page)

Aluminum Platform Dividers for Liquid Nitrogen Storage Tanks

For BAXTER/FENWAL BAGS



Platform Divider Part Number	Blood Bag Number(S)	Description	Dimension Height x Diameter	Canister	Frames
P20-4R2111	4R-2111	20" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F- C-4R2111.	20 x 20 inch 510 x 510 mm	C-4R2111	F-C-4R2111-3 F-C-4R2111-4
P24-4R2111		24" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F- C-4R2111.	20 x 24 inch 510 x 610 mm		
P31-4R2111		31" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F- C-4R2111.	20 x 31 inch 510 x 787 mm		
P39-4R2111		39" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F- C-4R2111.	20 x 39 inch 510 x 991 mm		
P20-4R2110	4R-2110	20" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F- C-4R2110.	20 x 20 inch 510 x 510 mm	C-4R2110	F-C-4R2110-3 F-C-4R2110-4
P24-4R2110		24" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F- C-4R2110.	20 x 24 inch 510 x 610 mm		
P31-4R2110		31" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F- C-4R2110.	20 x 31 inch 510 x 787 mm		
P39-4R2110		39" diameter Platform Divider for BAXTER/FENWAL bag Canister frames F- C-4R2110.	20 x 39 inch 510 x 991 mm		


For GAMBRO BAGS

Platform Divider Part Number	Blood Bag Number(S)	Description	Dimension Height x Diameter	Canister	Frames
P20-DF200	DF-200	20" diameter Platform Divider for GAMBRO bag Canister frames F-C-DF-200.	20 x 20 inch 510 x 510 mm	C-DF-200	F-C DF200-4 F-C DF200-6
P24-DF200		24" diameter Platform Divider for GAMBRO bag Canister frames F-C-DF-200.	20 x 24 inch 510 x 610 mm		
P31-DF200		31" diameter Platform Divider for GAMBRO bag Canister frames F-C-DF-200.	20 x 31 inch 510 x 787 mm		
P39-DF200		39" diameter Platform Divider for GAMBRO bag Canister frames F-C-DF-200.	20 x 39 inch 510 x 991 mm		
P20-DF700	DF-700	20" diameter Platform Divider for GAMBRO bag Canister frames F-C-DF-700.	20 x 20 inch 510 x 510 mm	C-DF-700	F-C DF700-4 F-C DF700-6
P24-DF700		24" diameter Platform Divider for GAMBRO bag Canister frames F-C-DF-700.	20 x 24 inch 510 x 610 mm		
P31-DF700		31" diameter Platform Divider for GAMBRO bag Canister frames F-C-DF-700.	20 x 31 inch 510 x 787 mm		
P39-DF700		39" diameter Platform Divider for GAMBRO bag Canister frames F-C-DF-700.	20 x 39 inch 510 x 991 mm		

(Continued on next page)

(Continued from last page)

Aluminum Platform Dividers for Liquid Nitrogen Storage Tanks



Inventory Control System
Canisters, Frames & Dividers


For CHARTERED BAGS

Platform Divider Part Number	Blood Bag Number(S)	Description	Dimension Height x Diameter	Canister	Frames
P20-7120	7120-2	20" diameter Platform Divider for CHARTERED bag Canister frames F- C-7120.	20 x 20 inch 510 x 510 mm	C-7120	F-C-7120-4 F-C-7120-5 F-C-7120-6
P24-7120		24" diameter Platform Divider for CHARTERED bag Canister frames F- C-7120.	20 x 24 inch 510 x 610 mm		
P31-7120		31" diameter Platform Divider for CHARTERED bag Canister frames F- C-7120.	20 x 31 inch 510 x 787 mm		
P39-7120		39" diameter Platform Divider for CHARTERED bag Canister frames F- C7120.	20 x 39 inch 510 x 991 mm		
P20-7005	7005-2 7303-2 7039-2	20" diameter Platform Divider for CHARTERED bag Canister frames F-7005.	20 x 20 inch 510 x 510 mm	C-7005	F-C-7005-2 F-C-7005-3
P24-7005		24" diameter Platform Divider for CHARTERED bag Canister frames F-7005.	20 x 24 inch 510 x 610 mm		
P31-7005		31" diameter Platform Divider for CHARTERED bag Canister frames F- C-7005.	20 x 31 inch 510 x 787 mm		
P39-7005		39" diameter Platform Divider for CHARTERED bag Canister frames F- C7005.	20 x 39 inch 510 x 991 mm		
P20-7400	7400-2 7403-2 7410-2 7025-2	20" diameter Platform Divider for CHARTERED bag Canister frames F- C-7400.	20 x 20 inch 510 x 510 mm	C-7400	F-C-7400-2
P24-7400		24" diameter Platform Divider for CHARTERED bag Canister frames F- C-7400.	20 x 24 inch 510 x 610 mm		
P31-7400		31" diameter Platform Divider for CHARTERED bag Canister frames F- C-7400.	20 x 31 inch 510 x 787 mm		
P39-7400		39" diameter Platform Divider for CHARTERED bag Canister frames F- C7400.	20 x 39 inch 510 x 991 mm		
P20-DF200	DF-200	20" diameter Platform Divider for CHARTERED bag Canister frames	20 x 20 inch 510 x 510 mm	C-DF-200	F-C DF200-4 F-C DF200-6
P24-DF200		24" diameter Platform Divider for CHARTERED bag Canister frames F-7005.	20 x 24 inch 510 x 610 mm		
P31-DF200		31" diameter Platform Divider for CHARTERED bag Canister frames	20 x 31 inch 510 x 787 mm		
P39-DF200		39" diameter Platform Divider for CHARTERED bag Canister frames	20 x 39 inch 510 x 991 mm		
P20-DF700	DF-700	20" diameter Platform Divider for CHARTERED bag Canister frames	20 x 20 inch 510 x 510 mm	C-DF-700	F-C DF700-4 F-C DF700-6
P24-DF700		24" diameter Platform Divider for CHARTERED bag Canister frames F-7005.	20 x 24 inch 510 x 610 mm		
P31-DF700		31" diameter Platform Divider for CHARTERED bag Canister frames	20 x 31 inch 510 x 787 mm		
P39-DF700		39" diameter Platform Divider for CHARTERED bag Canister frames	20 x 39 inch 510 x 991 mm		

(Continued from last page)

Inventory Control System

Aluminum Platform Dividers for Liquid Nitrogen Storage Tanks



For STERICON BAGS

Platform Divider Part Number	Blood Bag Number(S)	Description	Dimension Height x Diameter	Canister	Frames
P20-91D	RCM-91D	20" diameter Platform Divider for STERICON bag Canister frames F-C-91D.	20 x 20 inch 510 x 510 mm	C-91D	F-C 91D-3 F-C 91D-4 F-C 91D-5
P24-91D		24" diameter Platform Divider for STERICON bag Canister frames F-C-91D.	20 x 24 inch 510 x 610 mm		
P31-91D		31" diameter Platform Divider for STERICON bag Canister frames F-C-91D.	20 x 31 inch 510 x 787 mm		
P39-91D		39" diameter Platform Divider for STERICON bag Canister frames F-C-91D.	20 x 39 inch 510 x 991 mm		

For MED SEP BAGS

Platform Divider Part Number	Blood Bag Number(S)	Description	Dimension Height x Diameter	Canister	Fames
P20-CBBC-V1	791-OS 791-OSU 8-0346-O	20" diameter Platform Divider for MED SEP bag Canister frames F-CBBC-V1.	20 x 20 inch 510 x 510 mm	CBBC-V1	F-C-CBBC-V1-6 F-C-CBBC-V1-7 F-C-CBBC-V1-8
P24-CBBC-V1		24" diameter Platform Divider for MED SEP bag Canister frames F-CBBC-V1.	20 x 24 inch 510 x 610 mm		
P31-CBBC-V1		31" diameter Platform Divider for MED SEP bag Canister frames F-CBBC-V1..	20 x 31 inch 510 x 787 mm		
P39-CBBC-V1		39" diameter Platform Divider for MED SEP bag Canister frames F-CBBC-V1.	20 x 39 inch 510 x 991 mm		

(Continued on next page)

Cross Reference Tables for Blood Bag Storage Systems for Most Popular Liquid Nitrogen Storage Tanks

CHART-MVE ETERNE/HE SERIES

Frame Part Number (total frames total canisters)	BAXTER/FENWAL BAG NUMBER(S)						
	4R-9951	4R-9953 4R-5461	4R-9955 4R-5462	4R-9957 4R-5464	4R-9959 4R-5464	4R-2111	4R-2110
	C-4R9951	C-4R9953	C-4R9955	C-4R9957	C-4R9959	C-4R2111	C-4R2110
MVE 810HE/MVE-815P-150	F-C-4R9951-5(111/555) F-C-4R9951-7(111/777) F-C-4R9951-8(111/888)	F-C-4R9953-3(95/285) F-C-4R9953-4(95/380) F-C-4R9953-5(95/475)	F-C-4R9955-3(68/204) F-C-4R9955-4(68/272) F-C-4R9955-5(68/340)	F-C-4R9957-3(68/204) F-C-4R9957-4(68/272)		F-C-4R2111-3(44/132) F-C-4R2111-4(44/176)	F-C-4R2110-3(42/126)
MVE 1520HE/MVE-1536P-150	F-C-4R9951-5(248/1240) F-C-4R9951-7(248/1736) F-C-4R9951-8(248/1984)	F-C-4R9953-3(190/570) F-C-4R9953-4(190/760) F-C-4R9953-5(190/950)	F-C-4R9955-3(180/540) F-C-4R9955-4(180/720) F-C-4R9955-5(180/900)	F-C-4R9957-3(130/390) F-C-4R9957-4(130/520)	F-C-4R9959-3(96/288) F-C-4R9959-4(96/384)	F-C-4R2111-3(92/276) F-C-4R2111-4(92/368)	F-C-4R2110-3(80/240)
MVE1830HE/MVE-1892P-150	F-C-4R9951-5(480/2400) F-C-4R9951-7(480/3360) F-C-4R9951-8(480/3840)	F-C-4R9953-3(390/1170) F-C-4R9953-4(390/1560) F-C-4R9953-5(390/1950)	F-C-4R9955-3(276/828) F-C-4R9955-4(276/1104) F-C-4R9955-5(276/1380)	F-C-4R9957-3(240/720) F-C-4R9957-4(240/960)	F-C-4R9959-3(174/522) F-C-4R9959-4(174/696)	F-C-4R2111-3(180/540) F-C-4R2111-4(180/720)	F-C-4R2110-3(150/450)

Frame Part Number (total frames total canisters)	CHARTERMED		GAMBRO		STERICON	MED SEP	
	7120-2	7005-2 7303-2 7039-2	7400-2 7403-2 7410-2 7025-2	DF-200	DF-700	RCM-91D	791-05 791-05U 8-0346-0
	C-7120	C-7005	C-7400	C-DF-200	C-DF-700	C-91D	CBBC-V1
MVE 810HE/MVE-815P-150	F-C-7120-4(154/616) F-C-7120-5(154/770) F-C-7120-6(154/924)	F-C-7005-2(72/144) F-C-7005-3(72/216)	F-C-7400-2(82/164) F-C-7400-3(82/246)	F-C-DF200-4(52/108) F-C-DF200-6(52/312)	F-C-DF700-3(34/102) F-C-DF700-4(34/136)	F-C-91D-3(68/204) F-C-91D-4(68/272)	F-C-CBBC-V1-6(138/828) F-C-CBBC-V1-7(138/966) CF-8-3-T(26/1456)
MVE 1520HE/MVE-1536P-150	F-C-7120-4(350/1400) F-C-7120-5(350/1750) F-C-7120-6(350/2100)	F-C-7005-2(152/304) F-C-7005-3(152/456)	F-C-7400-2(150/300) F-C-7400-3(150/450)	F-C-DF200-4(104/416) F-C-DF200-6(104/624)	F-C-DF700-3(72/216) F-C-DF700-4(72/288)	F-C-91D-3(130/390) F-C-91D-4(130/520)	F-C-CBBC-V1-6(428/2568) F-C-CBBC-V1-7(428/2996) CF-9-3-T(32/3276)
MVE1830HE/MVE-1892P-150	F-C-7120-4(594/2376) F-C-7120-5(594/2970) F-C-7120-6(594/3560)	F-C-7005-2(306/612) F-C-7005-3(306/918)	F-C-7400-2(294/588) F-C-7400-3(294/882)	F-C-DF200-4(198/792) F-C-DF200-6(198/1188)	F-C-DF700-3(126/378) F-C-DF700-4(126/504)	F-C-91D-3(240/720) F-C-91D-4(240/960)	F-C-CBBC-V1-6(881/5286) F-C-CBBC-V1-7(881/6167) CF-9-3-T(108/6804)

CHART MVE 600/1400 SERIES

Frame Part Number (total frames total canisters)	BAXTER/FENWAL BAG NUMBER(S)						
	4R-9951	4R-9953 4R-5461	4R-9955 4R-5462	4R-9957 4R-5464	4R-9959 4R-5464	4R-2111	4R-2110
	C-4R9951	C-4R9953	C-4R9955	C-4R9957	C-4R9959	C-4R2111	C-4R2110
MVE 600/MVE-616C	F-C-4R9951-5(88/440) F-C-4R9951-7(88/616)	F-C-4R9953-3(66/198) F-C-4R9953-4(66/264) F-C-4R9953-5(66/330)	F-C-4R9955-3(56/168) F-C-4R9955-4(56/224) F-C-4R9955-5(56/280)	F-C-4R9957-3(46/138) F-C-4R9957-4(46/184) F-C-4R9957-5(46/230)	F-C-4R9959-3(36/108) F-C-4R9959-4(36/144) F-C-4R9959-5(36/180)	F-C-4R2111-3(28/84) F-C-4R2111-4(28/112)	F-C-4R2110-3(25/75)
Platform Divider Part NO. (For MVE 600/MVE-616C)	P24-4R9951	P24-4R9953	P24-4R9955	P24-4R9957	P24-4R9959	P24-4R2111	P24-4R2110
MVE 1400/MVE-1426CF	F-C-4R9951-5(120/600) F-C-4R9951-7(120/840)	F-C-4R9953-3(112/336) F-C-4R9953-4(112/448) F-C-4R9953-5(112/560)	F-C-4R9955-3(94/282) F-C-4R9955-4(94/376) F-C-4R9955-5(94/470)	F-C-4R9957-3(84/252) F-C-4R9957-4(84/336) F-C-4R9957-5(84/420)	F-C-4R9959-3(57/171) F-C-4R9959-4(57/228) F-C-4R9959-5(57/285)	F-C-4R2111-3(56/168) F-C-4R2111-4(56/224)	F-C-4R2110-3(25/75)
Platform Divider Part NO. (For MVE 1400/MVE-1426CF)	P24-4R9951	P24-4R9953	P24-4R9955	P24-4R9957	P24-4R9959	P24-4R2111	P24-4R2110

Frame Part Number (total frames total canisters)	CHARTERMED		GAMBRO		STERICON	MED SEP	
	7120-2	7005-2 7303-2 7039-2	7400-2 7403-2 7410-2 7025-2	DF-200	DF-700	RCM-91D	791-05 791-05U 8-0346-0
	C-7120	C-7005	C-7400	C-DF-200	C-DF-700	C-91D	CBBC-V1
MVE 600/MVE-616C	F-C-7120-4(88/352) F-C-7120-5(88/440) F-C-7120-6(88/528)	F-C-7005-2(48/96) F-C-7005-3(48/144)	F-C-7400-2(48/96) F-C-7400-3(48/144)	F-C-DF200-4(28/112) F-C-DF200-6(28/168)	F-C-DF700-4(24/96) F-C-DF700-6(24/144)	F-C-91D-3(46/138) F-C-91D-4(46/184) F-C-91D-5(46/230)	F-C-CBBC-V1-6(138/828) F-C-CBBC-V1-7(138/896) CF-9-3-T(22/1386)
Platform Divider Part NO. (For MVE 600/MVE-616C)	P24-7120	P24-7005	P24-7400	P24-DF200	P24-DF700	P24-91D	NONE
MVE 1400/MVE-1426CF	F-C-7120-4(148/592) F-C-7120-5(148/740) F-C-7120-6(148/888)	F-C-7005-2(96/192) F-C-7005-3(96/288)	F-C-7400-2(96/192) F-C-7400-3(96/288)	F-C-DF200-4(60/240) F-C-DF200-6(60/360)	F-C-DF700-4(40/160) F-C-DF700-6(40/240)	F-C-91D-3(84/252) F-C-91D-4(84/336) F-C-91D-5(84/420)	F-C-CBBC-V1-6(232/1392) F-C-CBBC-V1-7(232/1624) CF-9-3-T(35/2205)
Platform Divider Part NO. (For MVE 1400/MVE-1426CF)	P24-7120	P24-7005	P24-7400	P24-DF200	P24-DF700	P24-91D	NONE

(Continued on next page)

Cross Reference Tables for Blood Bag Storage Systems for Most Popular Liquid Nitrogen Storage Tanks

CHART MVE SERIES

Bag Manufacturer Bag Part Number Frame Part Number (total number of frames/total number of canisters) Canister Part Number Tank Model Number Old/New	BAXTER/FENWAL BAG NUMBER(S)						
	4R-9951	4R-9953 4R-5461	4R-9955 4R-5462	4R-9957 4R-5464	4R-9959 4R-5464	4R-2111	4R-2110
	C-4R9951	C-4R9953	C-4R9955	C-4R9957	C-4R9959	C-4R2111	C-4R2110
MVE 230/MVE-204	F-C-4R9951-5(36/180)	F-C-4R9953-3(29/87)	F-C-4R9955-3(21/63)	F-C-4R9957-3(18/54)	F-C-4R9959-3(13/39)		
MVE 230II/MVE-205	F-C-4R9951-5(36/180) F-C-4R9951-7(36/252)	F-C-4R9953-3(29/87) F-C-4R9953-4(29/116) F-C-4R9953-5(29/145)	F-C-4R9955-3(21/63) F-C-4R9955-4(21/84) F-C-4R9955-5(21/105)	F-C-4R9957-3(18/54) F-C-4R9957-4(18/72) F-C-4R9957-5(18/90)	F-C-4R9959-3(13/39) F-C-4R9959-4(13/52) F-C-4R9959-5(13/65)	F-C-4R2111-3(16/48) F-C-4R2111-4(16/64)	F-C-4R2110-3(13/39)
MVE 511/MVE-510	F-C-4R9951-5(56/280) F-C-4R9951-7(56/392) F-C-4R9951-8(56/448)	F-C-4R9953-3(43/129) F-C-4R9953-4(43/172) F-C-4R9953-5(43/215)	F-C-4R9955-3(32/96) F-C-4R9955-4(32/128) F-C-4R9955-5(32/160)	F-C-4R9957-3(30/90) F-C-4R9957-4(30/120) F-C-4R9957-5(30/150)	F-C-4R9959-3(21/63) F-C-4R9959-4(21/84) F-C-4R9959-5(21/105)	F-C-4R2111-3(20/60) F-C-4R2111-4(20/80)	F-C-4R2110-3(18/54)
Platform Divider Part No. (For MVE 511/MVE-510)	P20-4R9951	P20-4R9953	P20-4R9955	P20-4R9957	P20-4R9959	P20-4R2111	P20-4R2110
MVE 611/MVE-616	F-C-4R9951-5(88/440) F-C-4R9951-7(88/616)	F-C-4R9953-3(66/198) F-C-4R9953-4(66/264) F-C-4R9953-5(66/330)	F-C-4R9955-3(56/168) F-C-4R9955-4(56/224) F-C-4R9955-5(56/280)	F-C-4R9957-3(46/138) F-C-4R9957-4(46/184) F-C-4R9957-5(46/230)	F-C-4R9959-3(36/108) F-C-4R9959-4(36/144) F-C-4R9959-5(36/180)	F-C-4R2111-3(28/84) F-C-4R2111-4(28/112)	F-C-4R2110-3(25/75)
Platform Divider Part No. (For MVE 611/MVE-616)	P24-4R9951	P24-4R9953	P24-4R9955	P24-4R9957	P24-4R9959	P24-4R2111	P24-4R2110
MVE 1411/MVE-1426	F-C-4R9951-5(20/100) F-C-4R9951-7(20/140)	F-C-4R9953-3(112/336) F-C-4R9953-4(112/448) F-C-4R9953-5(112/560)	F-C-4R9955-3(94/282) F-C-4R9955-4(94/376)	F-C-4R9957-3(84/252) F-C-4R9957-4(84/336)	F-C-4R9959-3(57/171) F-C-4R9959-4(57/228)	F-C-4R2111-3(56/168) F-C-4R2111-4(56/224)	F-C-4R2110-3(52/156)
Platform Divider Part No. (For MVE 1411/MVE-1426)	P24-4R9951	P24-4R9953	P24-4R9955	P24-4R9957	P24-4R9959	P24-4R2111	P24-4R2110
MVE 1841/MVE-1839	F-C-4R9951-5(210/1050) F-C-4R9951-7(210/1470) F-C-4R9951-8(210/1680)	F-C-4R9953-3(156/468) F-C-4R9953-4(156/624) F-C-4R9953-5(156/780)	F-C-4R9955-3(144/432) F-C-4R9955-4(144/576) F-C-4R9955-5(144/720)	F-C-4R9957-3(136/408) F-C-4R9957-4(136/544) F-C-4R9957-5(136/680)	F-C-4R9959-3(106/318) F-C-4R9959-4(106/424) F-C-4R9959-5(106/530)	F-C-4R2111-3(92/276) F-C-4R2111-4(92/368)	F-C-4R2110-3(81/243)
Platform Divider Part No. (For MVE 1841/MVE-1839)	P39-4R9951	P39-4R9953	P39-4R9955	P39-4R9957	P39-4R9959	P39-4R2111	P39-4R2110

Bag Manufacturer Bag Part Number Frame Part Number (total number of frames/total number of canisters) Canister Part Number Tank Model Number	CHARTERMED			GAMBRO		STERICON	MED SEP
	7120-2	7005-2 7303-2 7039-2	7400-2 7403-2 7410-2 7025-2	DF-200	DF-700	RCM-91D	791-05 791-05U 8-0346-0
	C-7120	C-7005	C-7400	C-DF-200	C-DF-700	C-91D	CBBC-V1
MVE 230/MVE-204	F-C-7120-4(48/192)	F-C-7005-2(24/48)	F-C-7400-2(24/48)	F-C-DF200-4(16/64)		F-C-91D-3(18/54)	
MVE 230II/MVE-205	F-C-7120-4(48/192) F-C-7120-5(48/240) F-C-7120-6(48/288)	F-C-7005-2(24/48) F-C-7005-3(24/72)	F-C-7400-2(24/48) F-C-7400-3(24/72)	F-C-DF200-4(24/96) F-C-DF200-6(24/144)	F-C-DF700-3(10/30) F-C-DF700-4(10/40)	F-C-91D-3(18/54) F-C-91D-4(18/72) F-C-91D-5(18/90)	F-C-CBBC-V1-6(60/360) F-C-CBBC-V1-7(60/420) CF-9-3-T(8/448) CF-9-3-T(13/819)
MVE 511/MVE-510	F-C-7120-4(58/232) F-C-7120-5(58/290) F-C-7120-6(58/348)	F-C-7005-2(40/80) F-C-7005-3(40/120)	F-C-7400-2(40/80) F-C-7400-3(40/120)	F-C-DF200-4(20/80) F-C-DF200-6(20/120)	F-C-DF700-3(16/48) F-C-DF700-4(16/64)	F-C-91D-3(30/90) F-C-91D-4(30/120) F-C-91D-5(30/150)	F-C-CBBC-V1-6(90/540) F-C-CBBC-V1-7(90/630) CF-9-3-T(13/819)
Platform Divider Part No. (For MVE 511/MVE-510)	P20-7120	P20-7005	P20-7400	P20-DF200	P20-DF700	P20-91D	NONE
MVE 611/MVE-616	F-C-7120-4(88/352) F-C-7120-5(88/440) F-C-7120-6(88/528)	F-C-7005-2(48/96) F-C-7005-3(48/144)	F-C-7400-2(48/96) F-C-7400-3(48/144)	F-C-DF200-4(28/112) F-C-DF200-6(28/168)	F-C-DF700-3(24/72) F-C-DF700-4(24/96)	F-C-91D-3(46/138) F-C-91D-4(46/184) F-C-91D-5(46/230)	F-C-CBBC-V1-6(138/828) F-C-CBBC-V1-7(138/966) CF-9-3-T(22/1386)
Platform Divider Part No. (For MVE 611/MVE-616)	P24-7120	P24-7005	P24-7400	P24-DF200	P24-DF700	P24-91D	NONE
MVE 1411/MVE-1426	F-C-7120-4(148/592) F-C-7120-5(148/740) F-C-7120-6(148/888)	F-C-7005-2(96/192) F-C-7005-3(96/288)	F-C-7400-2(96/192) F-C-7400-3(96/288)	F-C-DF200-4(60/240) F-C-DF200-6(60/360)	F-C-DF700-3(40/120) F-C-DF700-4(40/160)	F-C-91D-3(84/252) F-C-91D-4(84/336)	F-C-CBBC-V1-6(232/1392) F-C-CBBC-V1-7(232/1624) CF-9-3-T(35/2205)
Platform Divider Part No. (For MVE 1411/MVE-1426)	P24-7120	P24-7005	P24-7400	P24-DF200	P24-DF700	P24-91D	NONE
MVE 1841/MVE-1839	F-C-7120-4(256/1024) F-C-7120-5(256/1280) F-C-7120-6(256/1536)	F-C-7005-2(148/296) F-C-7005-3(148/444)	F-C-7400-2(148/296) F-C-7400-3(148/444)	F-C-DF200-4(92/368) F-C-DF200-6(92/552)	F-C-DF700-3(68/204) F-C-DF700-4(68/272)	F-C-91D-3(136/408) F-C-91D-4(136/544) F-C-91D-5(136/680)	F-C-CBBC-V1-6(368/2208) F-C-CBBC-V1-7(368/2576) F-C-CBBC-V1-8(368/2944) CF-10-3-T(59/4130)
Platform Divider Part No. (For MVE 1841/MVE-1839)	P39-7120	P39-7005	P39-7400	P39-DF200	P39-DF700	P39-91D	NONE

(Continued on next page)

Cross Reference Tables for Blood Bag Storage Systems for Most Popular Liquid Nitrogen Storage Tanks

CHART MVE Bio Series

Frame Part Number (total frames total canisters)	Bag Manufacturer		BAXTER/FENWAL BAG NUMBER(S)						
	Bag Part Number	Canister Part Number	4R-9951	4R-9953 4R-5461	4R-9955 4R-5462	4R-9957 4R-5464	4R-9959 4R-5464	4R-2111	4R-2110
			Tank Model Number Old/New	C-4R9951	C-4R9953	C-4R9955	C-4R9957	C-4R9959	C-4R2111
BIO I	F-C-4R9951-5(36/180) F-C-4R9951-7(36/252)	F-C-4R9953-3(29/87) F-C-4R9953-4(26/116)	F-C-4R9955-3(21/63) F-C-4R9955-4(21/84)	F-C-4R9957-3(18/54) F-C-4R9957-4(18/72)	F-C-4R9959-3(13/39) F-C-4R9959-4(13/52)	F-C-4R2111-4(16/64)	F-C-4R2110-3(13/39)		
BIO II	F-C-4R9951-5(56/280) F-C-4R9951-7(56/392)	F-C-4R9953-3(43/129) F-C-4R9953-4(43/172)	F-C-4R9955-3(32/96) F-C-4R9955-4(32/128)	F-C-4R9957-3(30/90) F-C-4R9957-4(30/120)	F-C-4R9959-3(21/63) F-C-4R9959-4(21/84)	F-C-4R2111-3(20/60) F-C-4R2111-4(20/80)	F-C-4R2110-3(18/54)		
Platform Divider Part No. (For MVE BIO II)	P20-4R9951	P20-4R9953	P20-4R9955	P20-4R9957	P20-4R9959	P20-4R2111	P20-4R2110		
BIO III	F-C-4R9951-5(72/360) F-C-4R9951-7(72/504)	F-C-4R9953-3(66/198) F-C-4R9953-4(66/264) F-C-4R9953-5(66/330)	F-C-4R9955-3(48/144) F-C-4R9955-4(48/192) F-C-4R9955-5(48/240)	F-C-4R9957-3(46/138) F-C-4R9957-4(46/184) F-C-4R9957-5(46/230)	F-C-4R9959-3(36/108) F-C-4R9959-4(36/144) F-C-4R9959-5(36/180)	F-C-4R2111-3(28/84) F-C-4R2111-4(28/112)	F-C-4R2110-3(25/75) F-C-4R2110-4(25/100)		
Platform Divider Part No. (For MVE BIO III)	P24-4R9951	P24-4R9953	P24-4R9955	P24-4R9957	P24-4R9959	P24-4R2111	P24-4R2110		
BIO IV	F-C-4R9951-5(120/600) F-C-4R9951-7(120/840) F-C-4R9951-8(120/960)	F-C-4R9953-3(108/324) F-C-4R9953-4(108/432) F-C-4R9953-5(108/540)	F-C-4R9955-3(88/264) F-C-4R9955-4(88/352) F-C-4R9955-5(88/440)	F-C-4R9957-3(70/210) F-C-4R9957-4(70/280) F-C-4R9957-5(70/350)	F-C-4R9959-3(57/171) F-C-4R9959-4(57/228) F-C-4R9959-5(57/285)	F-C-4R2111-3(56/168) F-C-4R2111-4(56/224)	F-C-4R2110-3(52/156) F-C-4R2110-4(52/208)		
Platform Divider Part No. (For MVE BIO IV)	P31-4R9951	P31-4R9953	P31-4R9955	P31-4R9957	P31-4R9959	P31-4R2111	P31-4R2110		
BIO V	F-C-4R9951-5(210/1050) F-C-4R9951-7(210/1470) F-C-4R9951-8(210/1680)	F-C-4R9953-3(156/468) F-C-4R9953-4(156/624) F-C-4R9953-5(156/780)	F-C-4R9955-3(144/432) F-C-4R9955-4(144/578) F-C-4R9955-5(144/720)	F-C-4R9957-3(136/408) F-C-4R9957-4(136/544) F-C-4R9957-5(136/680)	F-C-4R9959-3(106/318) F-C-4R9959-4(106/424) F-C-4R9959-5(106/530)	F-C-4R2111-3(92/276) F-C-4R2111-4(92/368)	F-C-4R2110-3(81/243) F-C-4R2110-4(81/324)		
Platform Divider Part No. (For MVE BIO V)	P39-4R9951	P39-4R9953	P39-4R9955	P39-4R9957	P39-4R9959	P39-4R2111	P39-4R2110		

Frame Part Number (total frames total canisters)	Bag Manufacturer		CHARTERMED		GAMBRO		STERICON	MED SEP	
	Bag Part Number	Canister Part Number	7120-2	7005-2 7303-2 7039-2	7400-2 7403-2 7410-2 7025-2	DF-200	DF-700	RCM-91D	791-05 791-05U 8-0346-0
			Tank Model Number Old/New	C-7120	C-7005	C-7400	C-DF-200	C-DF-700	C-91D
BIO I	F-C-7120-4(48/192) F-C-7120-5(48/240) F-C-7120-6(48/288)	F-C-7005-2(24/48) F-C-7005-3(24/72)	F-C-7400-2(24/48)	F-C-DF200-4(24/96) F-C-DF200-6(24/144)	F-C-DF700-3(10/30) F-C-DF700-4(10/40)	F-C-91D-3(18/54) F-C-91D-4(18/72)	F-C-CBBC-V1-6(60/360) F-C-CBBC-V1-7(60/420) CF-8-3-T(8/48)		
BIO II	F-C-7120-4(58/232) F-C-7120-5(58/290) F-C-7120-6(58/348)	F-C-7005-2(40/80) F-C-7005-3(40/120)	F-C-7400-2(40/80) F-C-7400-3(40/120)	F-C-DF200-4(20/80) F-C-DF200-6(20/120)	F-C-DF700-3(16/48) F-C-DF700-4(16/64)	F-C-91D-3(30/90) F-C-91D-4(30/120)	F-C-CBBC-V1-6(90/540) F-C-CBBC-V1-7(90/630) CF-9-3-T(13/819)		
Platform Divider Part No. (For MVE BIO II)	P20-7120	P20-7005	P20-7400	P20-DF200	P20-DF700	P20-91D	NONE		
BIO III	F-C-7120-4(88/352) F-C-7120-5(88/440) F-C-7120-6(88/528)	F-C-7005-2(48/96) F-C-7005-3(48/144)	F-C-7400-2(48/96) F-C-7400-3(48/144)	F-C-DF200-4(28/112) F-C-DF200-6(28/168)	F-C-DF700-3(24/72) F-C-DF700-4(24/96)	F-C-91D-3(46/138) F-C-91D-4(46/184) F-C-91D-5(46/230)	F-C-CBBC-V1-6(138/828) F-C-CBBC-V1-7(138/966) CF-9-3-T(22/1386)		
Platform Divider Part No. (For MVE BIO III)	P24-7120	P24-7005	P24-7400	P24-DF200	P24-DF700	P24-91D	NONE		
BIO IV	F-C-7120-4(152/608) F-C-7120-5(152/760) F-C-7120-6(152/912)	F-C-7005-2(96/192) F-C-7005-3(96/288)	F-C-7400-2(96/192) F-C-7400-3(96/288)	F-C-DF200-4(56/224) F-C-DF200-6(56/336)	F-C-DF700-3(40/120) F-C-DF700-4(40/160)	F-C-91D-3(80/240) F-C-91D-4(80/320)	F-C-CBBC-V1-6(216/1296) F-C-CBBC-V1-7(216/1512) CF-9-3-T(33/2079)		
Platform Divider Part No. (For MVE BIO IV)	P31-7120	P31-7005	P31-7400	P31-DF200	P31-DF700	P31-91D	NONE		
BIO V	F-C-7120-4(256/1024) F-C-7120-5(256/1280) F-C-7120-6(256/1536)	F-C-7005-2(148/296) F-C-7005-3(148/444)	F-C-7400-2(148/296) F-C-7400-3(148/444)	F-C-DF200-4(92/368) F-C-DF200-6(92/552)	F-C-DF700-3(68/204) F-C-DF700-4(68/272)	F-C-91D-3(136/408) F-C-91D-4(136/544) F-C-91D-5(136/680)	F-C-CBBC-V1-6(368/2208) F-C-CBBC-V1-7(368/2576) CF-9-3-T(53/3339)		
Platform Divider Part No. (For MVE BIO V)	P39-7120	P39-7005	P39-7400	P39-DF200	P39-DF700	P39-91D	NONE		

TAYLOR-WHARTON LABS SERIES

Frame Part Number (total frames total canisters)	Bag Manufacturer		BAXTER/FENWAL BAG NUMBER(S)						
	Bag Part Number	Canister Part Number	4R-9951	4R-9953 4R-5461	4R-9955 4R-5462	4R-9957 4R-5464	4R-9959 4R-5464	4R-2111	4R-2110
			Tank Model Number	C-4R9951	C-4R9953	C-4R9955	C-4R9957	C-4R9959	C-4R2111
20K	F-C-4R9951-5(108/540) F-C-4R9951-7(108/756) F-C-4R9951-8(108/864)	F-C-4R9953-3(72/216) F-C-4R9953-4(72/288) F-C-4R9953-5(72/360)	F-C-4R9955-3(54/162) F-C-4R9955-4(54/216) F-C-4R9955-5(54/270)	F-C-4R9957-3(68/204) F-C-4R9957-4(68/272) F-C-4R9957-5(68/3340)		F-C-4R2111-3(50/150) F-C-4R2111-4(50/200)	F-C-4R2110-3(44/132)		
40K	F-C-4R9951-5(216/1080) F-C-4R9951-7(216/1512) F-C-4R9951-8(216/1728)	F-C-4R9953-3(174/525) F-C-4R9953-4(174/700) F-C-4R9953-5(174/870)	F-C-4R9955-3(118/354) F-C-4R9955-4(118/472) F-C-4R9955-5(118/590)	F-C-4R9957-3(96/288) F-C-4R9957-4(96/384) F-C-4R9957-5(96/480)	F-C-4R9959-3(72/216) F-C-4R9959-4(72/288) F-C-4R9959-5(72/360)	F-C-4R2111-3(90/270) F-C-4R2111-4(90/360)	F-C-4R2110-3(66/198)		
80K	F-C-4R9951-5(432/2160) F-C-4R9951-7(432/3024) F-C-4R9951-8(432/3456)	F-C-4R9953-3(324/972) F-C-4R9953-4(324/1296) F-C-4R9953-5(324/1620)	F-C-4R9955-3(240/720) F-C-4R9955-4(240/960) F-C-4R9955-5(240/1200)	F-C-4R9957-3(228/684) F-C-4R9957-4(228/912) F-C-4R9957-5(228/1140)	F-C-4R9959-3(162/486) F-C-4R9959-4(162/648) F-C-4R9959-5(162/810)	F-C-4R2111-3(174/522) F-C-4R2111-4(174/696)	F-C-4R2110-3(144/432)		

(Continued on next page)

Cross Reference Tables for Blood Bag Storage Systems for Most Popular Liquid Nitrogen Storage Tanks

Inventory Control Systems
Canisters, Frames & Dividers

Bag Manufacturer		CHARTERMED				GAMBRO		STERICON	MED SEP
		7120-2	7005-2 7303-2 7039-2	7400-2 7403-2 7410-2 7025-2		DF-200	DF-700	RCM-91D	791-05 791-05U 8-0346-0
Frame Part Number (total frames total canisters)		C-7120		C-7400	C-DF-200	C-DF-700	C-91D	CBBC-V1	
Tank Model Number		C-7120		C-7400	C-DF-200	C-DF-700	C-91D	CBBC-V1	
20K		F-C-4R7120-4(124/496) F-C-4R7120-5(124/620) F-C-4R7120-6(124/744)	F-C-4R7005-2(56/112) F-C-4R7005-3(56/168)	F-C-4R7400-2(56/112) F-C-4R7400-3(56/168)	F-C-DF200-4(48/384) F-C-DF200-6(48/576)	F-C-DF700-3(38/114) F-C-DF700-4(38/152)	F-C-91D-3(68/204) F-C-91D-4(68/272) F-C-91D-5(68/3340)	F-C-CBBC-V1-6(198/1188) F-C-CBBC-V1-7(198/1386) CF-9-3-T(26/1638)	
40K		F-C-4R7120-4(288/1152) F-C-4R7120-5(288/1440) F-C-4R7120-6(288/1728)	F-C-4R7005-2(138/276) F-C-4R7005-3(138/414)	F-C-4R7400-2(138/276) F-C-4R7400-3(138/414)	F-C-DF200-4(96/768) F-C-DF200-6(96/1152)	F-C-DF700-3(54/162) F-C-DF700-4(54/216)	F-C-91D-3(96/288) F-C-91D-4(96/384) F-C-91D-5(96/480)	F-C-CBBC-V1-6(402/2412) F-C-CBBC-V1-7(402/2814) CF-9-3-T(56/3528)	
80K		F-C-4R7120-4(570/2280) F-C-4R7120-5(570/2850) F-C-4R7120-6(570/3420)	F-C-4R7005-2(294/588) F-C-4R7005-3(294/882)	F-C-4R7400-2(294/588) F-C-4R7400-3(294/882)	F-C-DF200-4(192/1536) F-C-DF200-6(192/2304)	F-C-DF700-3(126/378) F-C-DF700-4(126/504)	F-C-91D-3(228/684) F-C-91D-4(228/912) F-C-91D-5(228/1140)	F-C-CBBC-V1-6(780/4680) F-C-CBBC-V1-7(780/5460) CF-9-3-T(112/7056)	

TAYLOR-WHARTON K SERIES

Bag Manufacturer		BAXTER/FENWAL BAG NUMBER(S)							
		4R-9951	4R-9953 4R-5461	4R-9955 4R-5462	4R-9957 4R-5464	4R-9959 4R-5464	4R-2111	4R-2110	
Frame Part Number (total frames total canisters)		C-4R9951		C-4R9953	C-4R9955	C-4R9957	C-4R9959	C-4R2111	C-4R2110
Tank Model Number		C-4R9951		C-4R9953	C-4R9955	C-4R9957	C-4R9959	C-4R2111	C-4R2110
3K		F-C-4R9951-5(25/125)		F-C-4R9953-3(20/60)	F-C-4R9955-3(14/42)	F-C-4R9957-3(8/24)	F-C-4R9959-3(2/6)		
10K		F-C-4R9951-5(56/280) F-C-4R9951-7(56/392)	F-C-4R9953-3(43/129) F-C-4R9953-4(43/172) F-C-4R9953-5(43/215)	F-C-4R9955-3(32/96) F-C-4R9955-4(32/128) F-C-4R9955-5(32/160)	F-C-4R9957-3(30/90) F-C-4R9957-4(30/120) F-C-4R9957-5(30/150)	F-C-4R9959-3(21/63) F-C-4R9959-4(21/84) F-C-4R9959-5(21/105)	F-C-4R2111-3(20/60) F-C-4R2111-4(20/80)	F-C-4R2110-3(18/54)	
Platform Divider Part No.(For Taylor Wharton 10K)		P20-4R9951		P20-4R9953	P20-4R9955	P20-4R9957	P20-4R9959	P20-4R2111	P20-4R2110
24K		F-C-4R9951-5(120/600) F-C-4R9951-7(120/840)	F-C-4R9953-3(108/324) F-C-4R9953-4(108/432) F-C-4R9953-5(108/540)	F-C-4R9955-3(88/264) F-C-4R9955-4(88/352) F-C-4R9955-5(88/440)	F-C-4R9957-3(80/240) F-C-4R9957-4(80/320) F-C-4R9957-5(80/400)	F-C-4R9959-3(57/171) F-C-4R9959-4(57/228) F-C-4R9959-5(57/285)	F-C-4R2111-3(56/168) F-C-4R2111-4(56/224)	F-C-4R2110-3(52/156)	
Platform Divider Part No.(For Taylor Wharton 24K)		P31-4R9951		P31-4R9953	P31-4R9955	P31-4R9957	P31-4R9959	P31-4R2111	P31-4R2110
38K		F-C-4R9951-5(210/1050) F-C-4R9951-7(210/1470)	F-C-4R9953-3(156/468) F-C-4R9953-4(156/624) F-C-4R9953-5(156/780)	F-C-4R9955-3(144/432) F-C-4R9955-4(144/578) F-C-4R9955-5(144/720)	F-C-4R9957-3(136/408) F-C-4R9957-4(136/544) F-C-4R9957-5(136/680)	F-C-4R9959-3(106/318) F-C-4R9959-4(106/424) F-C-4R9959-5(106/530)	F-C-4R2111-3(92/276) F-C-4R2111-4(92/368)	F-C-4R2110-3(81/243)	
Platform Divider Part No.(For Taylor Wharton 38K)		P39-4R9951		P39-4R9953	P39-4R9955	P39-4R9957	P39-4R9959	P39-4R2111	P39-4R2110

Bag Manufacturer		CHARTERMED				GAMBRO		STERICON	MED SEP
		7120-2	7005-2 7303-2 7039-2	7400-2 7403-2 7410-2 7025-2		DF-200	DF-700	RCM-91D	791-05 791-05U 8-0346-0
Frame Part Number (total frames total canisters)		C-7120		C-7400	C-DF-200	C-DF-700	C-91D	CBBC-V1	
Tank Model Number		C-7120		C-7400	C-DF-200	C-DF-700	C-91D	CBBC-V1	
3K		F-C-7120-4(36/144)		F-C-7400-2(18/36)	F-C-DF200-4(11/88)		F-C-91D-3(8/24)		
10K		F-C-7120-4(58/232) F-C-7120-5(58/290) F-C-7120-6(58/348)	F-C-7005-2(34/68) F-C-7005-3(34/102)	F-C-7400-2(34/68)	F-C-DF200-4(20/80) F-C-DF200-6(20/120)	F-C-DF700-4(16/64) F-C-DF700-6(16/96)	F-C-91D-3(30/90) F-C-91D-4(30/120) F-C-91D-5(30/150)	F-C-CBBC-V1-6(90/540) F-C-CBBC-V1-7(90/630) CF-9-3-T(14/882)	
Platform Divider Part No.(For Taylor Wharton 10K)		P20-7120		P20-7400	P20-DF200	P20-DF700	P20-91D	NONE	
24K		F-C-7120-4(152/608) F-C-7120-5(152/760) F-C-7120-6(152/912)	F-C-7005-2(96/192) F-C-7005-3(96/288)	F-C-7400-2(96/192)	F-C-DF200-4(56/224) F-C-DF200-6(56/336)	F-C-DF700-4(40/160) F-C-DF700-6(40/240)	F-C-91D-3(70/210) F-C-91D-4(70/280) F-C-91D-5(70/350)	F-C-CBBC-V1-6(216/1296) F-C-CBBC-V1-7(216/1512) CF-9-3-T(33/2079)	
Platform Divider Part No.(For Taylor Wharton 24K)		P31-7120		P31-7400	P31-DF200	P31-DF700	P31-91D	NONE	
38K		F-C-7120-4(256/1024) F-C-7120-5(256/1280) F-C-7120-6(256/1536)	F-C-7005-2(148/296) F-C-7005-3(148/444)	F-C-7400-2(148/296)	F-C-DF200-4(92/368) F-C-DF200-6(92/552)	F-C-DF700-4(68/272) F-C-DF700-6(68/408)	F-C-91D-3(136/408) F-C-91D-4(136/544) F-C-91D-5(136/680)	F-C-CBBC-V1-6(368/2208) F-C-CBBC-V1-7(368/2576) CF-9-3-T(53/3339)	
Platform Divider Part No.(For Taylor Wharton 38K)		P39-7120		P39-7400	P39-DF200	P39-DF700	P39-91D	NONE	