



Upright Freezer Drawer Racks

For Standard 3" Boxes



[Check page 108 for Boxes]

- ◆ Drawer Racks offer a more organized system for storage and retrieval of samples.
- ◆ Constructed of corrosion resistant stainless steel.
- ◆ Security lock devices are available for each rack.



Rack Only



Rack w/Cardboard Box



Rack w/Metal Box



w/DD Lock Device Option

Part Number	Description	Dimension L x H x W	Compartment Size l x h x w	Capacity (Box)	Lock Device Option**	Shipping Weight lb(Kg)
UFD-223	Rack Only	11 ³ / ₈ x 6 ⁷ / ₁₆ x 5 ¹ / ₂ inch 289 x 164 x 140 mm	5 ³ / ₈ x 3 x 5 ¹ / ₄ inch 137.4 x 75.9 x 133 mm	2 x 2 =4	DD-06	4.1(1.84)
UFD-233	Rack Only	11 ³ / ₈ x 10 ¹ / ₁₆ x 5 ¹ / ₂ inch 289 x 256 x 140 mm	5 ³ / ₈ x 3 ¹ / ₈ x 5 ¹ / ₄ inch 137.4 x 79.4 x 133 mm	2 x 3 =6	DD-09	5.7(2.59)
UFD-243	Rack Only	11 ³ / ₈ x 13 ⁷ / ₁₆ x 5 ¹ / ₂ inch 289 x 341 x 140 mm	5 ³ / ₈ x 3 ¹ / ₈ x 5 ¹ / ₄ inch 137.4 x 79.4 x 133 mm	2 x 4 =8	DD-13	7.4(3.35)
UFD-323	Rack Only	16 ³ / ₄ x 6 ⁷ / ₁₆ x 5 ¹ / ₂ inch 426 x 164 x 140 mm	5 ³ / ₈ x 3 x 5 ¹ / ₄ inch 137.4 x 75.9 x 133 mm	3 x 2 =6	DD-06	5.4(2.46)
UFD-333	Rack Only	16 ³ / ₄ x 10 ¹ / ₁₆ x 5 ¹ / ₂ inch 426 x 256 x 140 mm	5 ³ / ₈ x 3 ¹ / ₈ x 5 ¹ / ₄ inch 137.4 x 79.4 x 133 mm	3 x 3 =9	DD-09	7.9(3.58)
UFD-343	Rack Only	16 ³ / ₄ x 13 ⁷ / ₁₆ x 5 ¹ / ₂ inch 426 x 341 x 140 mm	5 ³ / ₈ x 3 ¹ / ₈ x 5 ¹ / ₄ inch 137.4 x 79.4 x 133 mm	3 x 4 =12	DD-13	10.0(4.54)
UFD-423	Rack Only	22 ⁵ / ₁₆ x 6 ⁷ / ₁₆ x 5 ¹ / ₂ inch 566 x 164 x 140 mm	5 ³ / ₈ x 3 x 5 ¹ / ₄ inch 137.4 x 75.9 x 133 mm	4 x 2 =8	DD-06	6.7(3.05)
UFD-433	Rack Only	22 ⁵ / ₁₆ x 10 ¹ / ₁₆ x 5 ¹ / ₂ inch 566 x 256 x 140 mm	5 ³ / ₈ x 3 ¹ / ₈ x 5 ¹ / ₄ inch 137.4 x 79.4 x 133 mm	4 x 3 =12	DD-09	9.7(4.39)
UFD-433-W-XX*	Rack w/3" box and divider					9.7(4.39)+12ω***
UFD-443	Rack Only	22 ⁵ / ₁₆ x 13 ¹ / ₂ x 5 ¹ / ₂ inch 566 x 341 x 140 mm	5 ³ / ₈ x 3 ¹ / ₈ x 5 ¹ / ₄ inch 137.4 x 79.4 x 133 mm	4 x 4 =16	DD-13	12.5(5.69)
UFD-453	Rack Only	22 ⁵ / ₁₆ x 16 ⁹ / ₁₆ x 5 ¹ / ₂ inch 566 x 420 x 140 mm	5 ³ / ₈ x 3 ¹ / ₈ x 5 ¹ / ₄ inch 137.4 x 79.4 x 133 mm	4 x 5 =20	DD-16	15.4(7.00)
UFD-523	Rack Only	27 ¹ / ₁₆ x 6 ⁷ / ₁₆ x 5 ¹ / ₂ inch 687 x 164 x 140 mm	5 ⁵ / ₁₆ x 3 x 5 ¹ / ₄ inch 134.7 x 75.9 x 133 mm	5 x 2 =10	DD-06	8.1(3.69)
UFD-533	Rack Only	27 ¹ / ₁₆ x 10 ¹ / ₁₆ x 5 ¹ / ₂ inch 687 x 256 x 140 mm	5 ⁵ / ₁₆ x 3 ¹ / ₈ x 5 ¹ / ₄ inch 134.7 x 79.4 x 133 mm	5 x 3 =15	DD-09	11.6(5.25)
UFD-543	Rack Only	27 ¹ / ₁₆ x 13 ⁷ / ₁₆ x 5 ¹ / ₂ inch 687 x 341 x 140 mm	5 ⁵ / ₁₆ x 3 ¹ / ₈ x 5 ¹ / ₄ inch 134.7 x 79.4 x 133 mm	5 x 4 =20	DD-13	15.0(6.81)


* To order, replace XX with either 81 or 100 cell divider.

** Lock device is ordered separately.

*** ω¹=0.22lb(0.1kg) for cardboard box, 0.37lb(0.17kg) for plastic box, 0.84lb(0.38kg) for stainless steel box, 0.53lb(0.24kg) for aluminum box.


Upright Sliding Tray Racks

For Standard 3" Boxes



[Check page 108 for Boxes]

- ◆ Designed for best cost performance. Function is about the same as our upright freezer drawer racks.
- ◆ Both rack frame and sliding tray are constructed of corrosion resistant stainless steel.
- ◆ Offer a more organized system for storage and retrieval of samples.



Part Number	Description	Dimension L x H x W	Compartment Size l x h x w	Capacity (Box)	Shipping Weight lb(Kg)
UFS-333	Rack Only	16 ⁷ / ₈ x 9 ¹ / ₂ x 5 ¹ / ₂ inch 428 x 241 x 140 mm	5 ³ / ₈ x 3 ¹ / ₁₆ x 5 ¹ / ₄ inch 137.4 x 76.7 x 134 mm	3 x 3 =9	4.4(2.70)
UFS-343	Rack Only	16 ⁷ / ₈ x 12 ¹¹ / ₁₆ x 5 ¹ / ₂ inch 428 x 322 x 140 mm	5 ³ / ₈ x 3 ¹ / ₁₆ x 5 ¹ / ₄ inch 137.4 x 76.7 x 134 mm	3 x 4 =12	7.9(3.60)
UFS-423	Rack Only	22 ³ / ₈ x 6 ⁵ / ₁₆ x 5 ¹ / ₂ inch 568 x 161 x 140 mm	5 ³ / ₈ x 3 ¹ / ₁₆ x 5 ¹ / ₄ inch 137.4 x 76.7 x 134 mm	4 x 2 =8	4.9(2.22)
UFS-433	Rack Only	22 ³ / ₈ x 9 ¹ / ₂ x 5 ¹ / ₂ inch 568 x 241 x 140 mm	5 ³ / ₈ x 3 ¹ / ₁₆ x 5 ¹ / ₄ inch 137.4 x 76.7 x 134 mm	4 x 3 =12	7.4(3.37)
UFS-443	Rack Only	22 ³ / ₈ x 12 ¹¹ / ₁₆ x 5 ¹ / ₂ inch 568 x 322 x 140 mm	5 ³ / ₈ x 3 ¹ / ₁₆ x 5 ¹ / ₄ inch 137.4 x 76.7 x 134 mm	4 x 4 =16	9.9(4.51)
UFS-453	Rack Only	22 ³ / ₈ x 15 ¹³ / ₁₆ x 5 ¹ / ₂ inch 568 x 402 x 140 mm	5 ³ / ₈ x 3 ¹ / ₁₆ x 5 ¹ / ₄ inch 137.4 x 76.7 x 134 mm	4 x 5 =20	12.5(5.65)
UFS-523	Rack Only	27 ¹ / ₁₆ x 6 ⁵ / ₁₆ x 5 ¹ / ₂ inch 688 x 161 x 140 mm	5 ⁵ / ₁₆ x 3 ¹ / ₁₆ x 5 ¹ / ₄ inch 134.7 x 76.7 x 134 mm	5 x 2 =10	5.9(2.68)
UFS-533	Rack Only	27 ¹ / ₁₆ x 9 ¹ / ₂ x 5 ¹ / ₂ inch 688 x 241 x 140 mm	5 ⁵ / ₁₆ x 3 ¹ / ₁₆ x 5 ¹ / ₄ inch 134.7 x 76.7 x 134 mm	5 x 3 =15	8.9(4.04)
UFS-543	Rack Only	27 ¹ / ₁₆ x 12 ¹¹ / ₁₆ x 5 ¹ / ₂ inch 688 x 322 x 140 mm	5 ⁵ / ₁₆ x 3 ¹ / ₁₆ x 5 ¹ / ₄ inch 134.7 x 76.7 x 134 mm	5 x 4 =20	11.9(5.40)