



Upright Freezer Racks

For 15ml & 50ml Centrifuge Tubes Boxes



- ◆ Racks are for standard cardboard 15ml& 50ml centrifuge tubes boxes.
- ◆ Made of corrosion resistant stainless steel.
- ◆ Security lock devices are available for each rack.



Chest Freezer Rack w/Cardboard Box



w/DR Lock Device Option


Part Number	Description	Dimension L x H x W	Compartment Size l x h x w	Capacity (Box)	Lock Device Option*	Shipping Weight (lb/Kg)
UFLB-22	Rack Only	12 ³ / ₄ x 10 ³ / ₁₆ x 5 ⁷ / ₈ inch 324 x 259 x 149 mm	6 x 5 ¹ / ₁₆ x 5 ¹³ / ₁₆ inch 152 x 128 x 147.5 mm	2 x 2 =4	DP-1210	3.7(1.67)
UFLB-32	Rack Only	18 ¹³ / ₁₆ x 10 ³ / ₁₆ x 5 ⁷ / ₈ inch 478 x 259 x 149 mm	6 x 5 ¹ / ₁₆ x 5 ¹³ / ₁₆ inch 152 x 128 x 147.5 mm	3 x 2 =6	DP-1810	5.2(2.36)
UFLB-42	Rack Only	24 ¹³ / ₁₆ x 10 ³ / ₁₆ x 5 ⁷ / ₈ inch 631 x 259 x 149 mm	6 x 5 ¹ / ₁₆ x 5 ¹³ / ₁₆ inch 152 x 128 x 147.5 mm	4 x 2 =8	DP-2410	6.6(2.98)

* Lock device ordered separately.


Chest Freezer Rack

(Locking Rod Included)

For 15ml & 50ml Centrifuge Tubes Boxes



- ◆ Racks are for standard cardboard 15ml&50ml centrifuge tubes boxes.
- ◆ Made of corrosion resistant stainless steel.
- ◆ Security lock devices are available for each rack.



Rack w/DR Lock
Device Option



Part Number	Description	Dimension H x D x W	Compartment Size h x d x w	Capacity (Box)	Lock Device Option*	Shipping Weight (lb/Kg)
CFLB-4	Rack Only	20 ³ / ₈ x 6 x 6 ³ / ₁₆ inch 517 x 153 x 157 mm	5 x 5 ⁷ / ₈ x 5 ⁷ / ₈ inch 126.8 x 149 x 149 mm	4	DR-19	2.7(1.22)
CFLB-5	Rack Only	25 ³ / ₈ x 6 x 6 ³ / ₁₆ inch 644 x 153 x 157 mm	5 x 5 ⁷ / ₈ x 5 ⁷ / ₈ inch 126.8 x 149 x 149 mm	5	DR-23	3.2(1.45)
CFLB-6	Rack Only	30 ³ / ₈ x 6 x 6 ³ / ₁₆ inch 771 x 153 x 157 mm	5 x 5 ⁷ / ₈ x 5 ⁷ / ₈ inch 126.8 x 149 x 149 mm	6	DR-29	3.7(1.68)

* Lock device ordered separately.

Upright Freezer Drawer Racks

For Centrifuge Tubes

- ◆ Racks are specially designed to provide an ideal alternative storing large quantities of samples in 15ml and 50ml centrifuge tubes.
- ◆ Made of corrosion resistant stainless steel.
- ◆ Drawer Racks facilitate both access to and storage of samples without having to remove the whole racks from the freezer.
- ◆ Security lock devices are available for each rack.



Content to Store	Rack Length	Rack Height	Number of Drawers	Part Number	Description	Dimension L x H x W	Capacity	Lock Device Option *	Shipping Weight (lb/kg)
15ml 	16 ¹ / ₂ inch 424mm	4 ¹⁵ / ₁₆ inch 126 mm	1	UFD-T15-1	Rack with Drawer for 15ml Centrifuge Tubes	16 ¹ / ₁₆ x 4 ¹⁵ / ₁₆ x 5 ¹ / ₂ inch 424 x 126 x 140 mm	60 Tubes (15 ml)	DD-04	4.0 / 1.83
	22 inch 564mm			UFD-LT15-1		22 ³ / ₁₆ x 4 ¹⁵ / ₁₆ x 5 ¹ / ₂ inch 564 x 126 x 140 mm	80 Tubes (15 ml)		5.2 / 2.35
	26 ³ / ₄ inch 685mm			UFD-XLT15-1		26 ¹⁵ / ₁₆ x 4 ¹⁵ / ₁₆ x 5 ¹ / ₂ inch 685 x 126 x 140 mm	100 Tubes (15 ml)		6.1 / 2.75
	16 ¹ / ₂ inch 424mm	9 ¹⁵ / ₁₆ inch 252 mm	2	UFD-T15-2		16 ¹ / ₁₆ x 9 ¹⁵ / ₁₆ x 5 ¹ / ₂ inch 424 x 252 x 140 mm	120 Tubes (15 ml)	DD-09	7.4 / 3.34
	22 inch 564mm			UFD-LT15-2		22 ³ / ₁₆ x 9 ¹⁵ / ₁₆ x 5 ¹ / ₂ inch 564 x 252 x 140 mm	160 Tubes (15 ml)		9.3 / 4.22
	26 ³ / ₄ inch 684mm			UFD-XLT15-2		26 ¹⁵ / ₁₆ x 9 ¹⁵ / ₁₆ x 5 ¹ / ₂ inch 684 x 252 x 140 mm	200 Tubes (15 ml)		10.8 / 4.90
50ml 	16 ¹ / ₂ inch 424mm	4 ¹⁵ / ₁₆ inch 126 mm	1	UFD-T50-1	Rack with Drawer for 50ml Centrifuge Tubes	16 ¹ / ₁₆ x 4 ¹⁵ / ₁₆ x 5 ¹ / ₂ inch 424 x 126 x 140 mm	30 Tubes (50 ml)	DD-04	3.9 / 1.78
	22 inch 564mm			UFD-LT50-1		22 ³ / ₁₆ x 4 ¹⁵ / ₁₆ x 5 ¹ / ₂ inch 564 x 126 x 140 mm	39 Tubes (50 ml)		5.1 / 2.30
	26 ³ / ₄ inch 685mm			UFD-XLT50-1		26 ¹⁵ / ₁₆ x 4 ¹⁵ / ₁₆ x 5 ¹ / ₂ inch 685 x 126 x 140 mm	51 Tubes (50 ml)		5.9 / 2.68
	16 ¹ / ₂ inch 424mm	9 ¹⁵ / ₁₆ inch 252 mm	2	UFD-T50-2		16 ¹ / ₁₆ x 9 ¹⁵ / ₁₆ x 5 ¹ / ₂ inch 424 x 252 x 140 mm	60 Tubes (50 ml)	DD-09	7.2 / 3.25
	22 inch 564mm			UFD-LT50-2		22 ³ / ₁₆ x 9 ¹⁵ / ₁₆ x 5 ¹ / ₂ inch 564 x 252 x 140 mm	78 Tubes (50 ml)		9.1 / 4.12
	26 ³ / ₄ inch 684mm			UFD-XLT50-2		26 ¹⁵ / ₁₆ x 9 ¹⁵ / ₁₆ x 5 ¹ / ₂ inch 684 x 252 x 140 mm	102 Tubes (50 ml)		10.4 / 4.76
15/50ml 	16 ¹ / ₂ inch 424mm	9 ¹⁵ / ₁₆ inch 252 mm	2	UFD-T15/50	Rack with Drawer for 15/50ml Centrifuge Tubes	16 ¹ / ₁₆ x 9 ¹⁵ / ₁₆ x 5 ¹ / ₂ inch 424 x 252 x 140mm	60 Tubes(15ml) + 30 Tubes (50 ml)	DD-04	4.6 / 2.10
	22 inch 564mm			UFD-LT15/50		22 ³ / ₁₆ x 9 ¹⁵ / ₁₆ x 5 ¹ / ₂ inch 564 x 252 x 140mm	80 Tubes(15ml) + 39 Tubes (50 ml)		5.5 / 2.50
	26 ³ / ₄ inch 685mm			UFD-XLT15/50		26 ¹⁵ / ₁₆ x 9 ¹⁵ / ₁₆ x 5 ¹ / ₂ inch 685 x 252 x 140mm	100 Tubes(15ml) + 51 Tubes (50 ml)		5.5 / 2.50

* Lock device is ordered separately.