


Security Lock Devices for Chest Freezer & Liquid Nitrogen Tank Racks

Security lock devices for chest freezer and liquid nitrogen tank racks are made of high quality stainless steel and secured with a lock included with each system. Select your lock device by identifying your specific chest freezer or liquid nitrogen tank rack and then the security lock device for your particular freezer rack.



Lock Device	Lock Dimension	Applicable to Rack	Rack Height
	H		
DR-06	6 ⁷ / ₁₆ inch 160 mm	CF-3-2 CFA-3-2 MCF-3-2	6 ⁷ / ₈ inch 168 mm
DR-06S	6 ⁷ / ₁₆ inch 160 mm	CFM-305	6 ³ / ₄ inch 171 mm
		CFM-315	6 ¹⁵ / ₁₆ inch 176 mm
DR-09	8 ¹¹ / ₁₆ inch 220 mm	CF-4-2 CFA-4-2 MCF-4-2	8 ¹³ / ₁₆ inch 224 mm
DR-09S	8 ¹¹ / ₁₆ inch 220 mm	CFM-405	9 inch 228 mm
		CFM-415	9 ¹ / ₄ inch 235 mm
DR-11	10 ¹³ / ₁₆ inch 275 mm	CF-5-2 CFA-5-2 MCF-5-2	11 inch 280 mm
DR-11S	10 ¹³ / ₁₆ inch 275 mm	CFM-505	11 ¹ / ₄ inch 285 mm
		CFM-515	11 ⁹ / ₁₆ inch 293 mm
DR-12	12 ⁵ / ₈ inch 320 mm	CF-6-2 CFA-6-2 MCF-6-2	13 ¹ / ₄ inch 336 mm
DR-12S	12 ² / ₈ inch 320 mm	CF-4-3 CFA-4-3 MCF-4-3	12 ¹¹ / ₁₆ inch 323 mm
		CFM-605	13 ⁷ / ₁₆ inch 341 mm
		CFM-615	13 ⁷ / ₈ inch 352 mm
DR-15	15 ³ / ₈ inch 390 mm	CF-7-2 CFA-7-2 MCF-7-2	15 ⁷ / ₁₆ inch 392 mm
		CF-5-3 CFA-5-3 MCF-5-3	15 ⁷ / ₈ inch 404 mm
		CFDP-9	16 ⁷ / ₁₆ inch 417 mm
DR-15S	15 ³ / ₈ inch 390 mm	CFMP-5	16 ⁹ / ₁₆ inch 421 mm
		CFM-705	15 ¹¹ / ₁₆ inch 398 mm
DR-17	17 ¹ / ₂ inch 445 mm	CFM-715	16 ¹ / ₈ inch 410 mm
		CF-8-2 CFA-8-2 MCF-8-2	17 ³ / ₈ inch 447 mm
		CF-6-3 CFA-6-3 MCF-6-3	19 ¹ / ₈ inch 485 mm
		MCF-6-3	19 ¹ / ₈ inch 485 mm
		CFHT-8	18 ⁷ / ₁₆ inch 469 mm
		CFTB-8	17 ¹¹ / ₁₆ inch 449 mm
DR-17S	17 ¹ / ₂ inch 445 mm	CFDP-10	18 ¹ / ₄ inch 463 mm
		CFM-805	17 ¹⁵ / ₁₆ inch 455 mm
DR-19	19 ¹ / ₂ inch 496 mm	CFM-815	18 ⁷ / ₁₆ inch 469 mm
		CF-9-2 CFA-9-2	19 ¹³ / ₁₆ inch 503 mm
		MCF-9-2 MCF-9-2-S MCFA-9-2	19 ¹³ / ₁₆ inch 503 mm
		CFHT-9	20 ³ / ₄ inch 527 mm
		CFLB-4	20 ¹ / ₁₆ inch 510 mm
		CFTB-9	19 ¹⁵ / ₁₆ inch 506 mm
		CFDP-11	20 ¹ / ₁₆ inch 509 mm
DR-19S	19 ¹ / ₂ inch 496 mm	CFMP-6	19 ⁷ / ₈ inch 505 mm
		CFM-905	20 ³ / ₁₆ inch 512 mm
		CFM-915	20 ³ / ₄ inch 527 mm

*Protect your samples from unauthorized access using these security lock devices.

*No tools are required for installing lock devices.

*Locking devices are made of durable stainless steel.

(Continued on next page)

Security Lock Devices for Chest Freezer & Liquid Nitrogen Tank Racks

Lock device	Lock Dimension	Applicable to Rack	Rack Height
	H		
DR-22	21 ¹¹ / ₁₆ inch 550 mm	CF-10-2 CF-10-2-S CFA-10-2	22 inch 559 mm
		MCF-10-2 MCF-10-2-S MCFA-10-2	22 inch 559 mm
		CF-7-3 CF-7-3-S CFA-7-3	22 ¹ / ₄ inch 565 mm
		MCF-7-3 MCFA-7-3	22 ¹ / ₄ inch 565 mm
		CFHT-10	23 ¹ / ₁₆ inch 586 mm
		CFTB-10	22 ¹ / ₈ inch 562 mm
		CFDP-12	21 ⁷ / ₈ inch 556 mm
		CFMP-7	23 ³ / ₁₆ inch 589 mm
DR-22S	21 ¹¹ / ₁₆ inch 550 mm	CFM-1005	22 ³ / ₈ inch 569 mm
		CFM-1015	23 ¹ / ₁₆ inch 586 mm
DR-23	23 ³ / ₈ inch 595 mm	CF-11-2 CF-11-2-S CFA-11-2	24 ¹ / ₁₆ inch 615 mm
		MCF-11-2 MCF-11-2-S MCFA-11-2	24 ³ / ₁₆ inch 615 mm
		CFTB-11	24 ⁵ / ₁₆ inch 618 mm
		CFDP-13	23 ¹¹ / ₁₆ inch 602 mm
		CFLB-5	25 ¹ / ₁₆ inch 637 mm
DR-23S	23 ³ / ₈ inch 595 mm	CFM-1105	24 ³ / ₈ inch 625 mm
DR-25	25 ¹ / ₈ inch 638 mm	CF-8-3 CF-8-3-S CFA-8-3	25 ⁷ / ₁₆ inch 646 mm
		MCF-8-3 MCFA-8-3	25 ⁷ / ₁₆ inch 646 mm
		CFHT-11	25 ³ / ₈ inch 644 mm
		CFDP-14	25 ¹ / ₂ inch 648 mm
		CFMP-8	26 ¹ / ₂ inch 673 mm
DR-25S	25 ¹ / ₈ inch 638 mm	CFM-1115	25 ³ / ₈ inch 644 mm
DR-26	26 inch 660 mm	CF-12-2 CF-12-2-S CFA-12-2	26 ³ / ₈ inch 670 mm
		MCF-12-2 MCF-12-2-S MCFA-12-2	26 ³ / ₈ inch 670 mm
		CFTB-12	26 ⁷ / ₁₆ inch 674mm
DR-26S	26 inch 660 mm	CFM-1205	26 ⁷ / ₈ inch 682 mm
DR-27	27 ³ / ₁₆ inch 690 mm	CF-13-2 CF-13-2-S CFA-13-2	28 ⁹ / ₁₆ inch 726 mm
		MCF-13-2 MCF-13-2-S MCFA-13-2	28 ⁹ / ₁₆ inch 726 mm
		CF-9-3 CF-9-3-S CFA-9-3	28 ⁹ / ₁₆ inch 726 mm
		MCF-9-3 MCFA-9-3	28 ⁹ / ₁₆ inch 726 mm
		CFHT-12	27 ¹¹ / ₁₆ inch 703 mm
		CFDP-15	27 ⁵ / ₁₆ inch 694 mm
DR-27S	27 ³ / ₁₆ inch 690 mm	CFM-1215	27 ¹¹ / ₁₆ inch 703 mm
DR-29	28 ⁵ / ₈ inch 726 mm	CFHT-13	29 ¹⁵ / ₁₆ inch 761 mm
		CFLB-6	30 ¹ / ₁₆ inch 764 mm
		CFDP-16	29 ³ / ₁₆ inch 741 mm
		CFMP-9	29 ¹³ / ₁₆ inch 757 mm
		CFTB-13	28 ³ / ₄ inch 730 mm
DR-29S	28 ⁵ / ₈ inch 726 mm	CFM-1305	29 ¹ / ₈ inch 739 mm
		CFM-1315	29 ¹⁵ / ₁₆ inch 761 mm
DR-31	31 inch 776 mm	CF-14-2 CF-14-2-S CFA-14-2	30 ¹³ / ₁₆ inch 782 mm
		MCF-14-2 MCF-14-2-S MCFA-14-2	30 ¹³ / ₁₆ inch 782 mm
		CF-10-3 CF-10-3-S CAF-10-3	31 ³ / ₄ inch 807 mm
		MCF-10-3 MCFA-10-3	31 ³ / ₄ inch 807 mm
DR-32	32 ⁵ / ₁₆ inch 820 mm	CF-15-2 CF-15-2-S CFA-15-2	33 inch 838 mm
		MCF-15-2 MCF-15-2-S MCFA-15-2	33 inch 838 mm

*Protect your samples from unauthorized access using these security lock devices.


*No tools are required for installing lock devices.

*Locking devices are made of durable stainless steel.

(Continued on next page)

Security Lock Devices for Upright Freezer Racks

Security lock devices for upright freezer racks are made of high quality stainless steel and secured with a lock included with each system. Select your lock device by identifying your specific upright freezer rack and then the security lock device for your particular freezer rack.



Lock Device	Lock Dimension		Applicable to Rack	Rack Height
	L	H		
DP-0507	4 ³ / ₄ inch 120 mm	6 ³ / ₄ -7 ¹ / ₈ inch 171-181 mm	UFM-1305	6 ⁵ / ₈ inch 169 mm
DP-0509	4 ³ / ₄ inch 120 mm	8 ⁷ / ₈ -9 ⁵ / ₁₆ inch 226-236 mm	UFM-1405	8 ¹³ / ₁₆ inch 224 mm
DP-0511	4 ³ / ₄ inch 120 mm	11 ¹ / ₈ -11 ¹ / ₂ inch 282-292 mm	UFM-1505	11 inch 280 mm
DP-0607	5 ¹ / ₂ inch 140 mm	7 - 7 ⁷ / ₁₆ inch 178-188 mm	UFM-1315	6 ¹⁵ / ₁₆ inch 176 mm
DP-0609	5 ¹ / ₂ inch 140 mm	9 ⁵ / ₁₆ -9 ³ / ₄ inch 237-247 mm	UFM-1415	9 ¹ / ₄ inch 235 mm
DP-0612	5 ¹ / ₂ inch 140 mm	11 ⁵ / ₈ -12 inch 295-305 mm	UFM-1515	11 ⁹ / ₁₆ inch 293 mm
DP-1107	10 ³ / ₁₆ inch 259 mm	6 ¹¹ / ₁₆ -7 ¹ / ₈ inch 171-181 mm	UFM-2305	6 ⁵ / ₈ inch 169 mm
DP-1109	10 ³ / ₁₆ inch 259 mm	8 ⁷ / ₈ -9 ⁵ / ₁₆ inch 226-236 mm	UFM-2405	8 ¹³ / ₁₆ inch 224 mm
DP-1111	10 ³ / ₁₆ inch 259 mm	11 ¹ / ₈ -11 ¹ / ₂ inch 282-292 mm	UFM-2505	11 inch 280 mm
DP-1210	11 ⁷ / ₁₆ inch 290 mm	9 ¹⁵ / ₁₆ -10 ³ / ₈ inch 253-263 mm	UFLB-22	10 ³ / ₁₆ inch 259 mm
DP-1307	11 ³ / ₄ inch 298 mm	7 - 7 ⁷ / ₁₆ inch 178-188 mm	UFHT-23	6 ¹⁵ / ₁₆ inch 176 mm
			UFM-2315	6 ¹⁵ / ₁₆ inch 176 mm
DP-1309	11 ³ / ₄ inch 298 mm	9 ⁵ / ₁₆ -9 ³ / ₄ inch 237-247 mm	UFHT-24	9 ¹ / ₄ inch 235 mm
			UFM-2415	9 ¹ / ₄ inch 235 mm
DP-1312	11 ³ / ₄ inch 298 mm	11 ⁵ / ₈ -12 inch 295-305 mm	UFHT-25	11 ⁹ / ₁₆ inch 293 mm
			UFM-2515	11 ⁹ / ₁₆ inch 293 mm
DP-1607	15 ⁹ / ₁₆ inch 395 mm	6 ¹¹ / ₁₆ -7 ¹ / ₈ inch 170-180 mm	UF-332	6 ⁵ / ₈ inch 168 mm
			UF-323	6 ⁵ / ₈ inch 168 mm
DP-1609	15 ⁹ / ₁₆ inch 395 mm	8 ¹⁵ / ₁₆ -9 ⁵ / ₁₆ inch 227-237 mm	UF-342	8 ⁷ / ₈ inch 225 mm
DP-1610	15 ⁹ / ₁₆ inch 395 mm	9 ⁷ / ₁₆ -9 ¹³ / ₁₆ inch 239-249 mm	UF-333	9 ³ / ₈ inch 238 mm
DP-1611	15 ⁹ / ₁₆ inch 395 mm	11 ¹ / ₈ -11 ¹ / ₂ inch 282-292 mm	UF-352	11 inch 280 mm
DP-1613	15 ⁹ / ₁₆ inch 395 mm	12 ³ / ₄ -13 ¹ / ₈ inch 323-333 mm	UF-362	12 ³ / ₈ inch 321mm
DP-1707	15 ¹¹ / ₁₆ inch 398 mm	6 ³ / ₄ -7 ¹ / ₈ inch 171-181 mm	UFM-3305	6 ⁵ / ₈ inch 169 mm
DP-1709	15 ¹¹ / ₁₆ inch 398 mm	8 ⁷ / ₈ -9 ⁵ / ₁₆ inch 226-236 mm	UFM-3405	8 ¹³ / ₁₆ inch 224 mm
DP-1711	15 ¹¹ / ₁₆ inch 398 mm	11 ¹ / ₈ -11 ¹ / ₂ inch 282-292 mm	UFM-3505	11 inch 280 mm

*Protect your samples from unauthorized access using these security lock devices.

*No tools are required for installing lock devices.

*Locking devices are made of durable stainless steel.

(Continued on next page)

(Continued from last page)

Security Lock Devices for Upright Freezer Racks

Inventory Control System

Lock Device	Lock Dimension		Applicable to Rack	Rack Height
	L	H		
DP-1809	17 ⁹ / ₁₆ inch 446 mm	8 ⁷ / ₈ -9 ⁵ / ₁₆ inch 226-236 mm	UFTB-54	8 ¹³ / ₁₆ inch 224 mm
		9 ¹ / ₄ -10 ¹ / ₁₆ inch 235-255 mm	UFDP-55	9 ⁷ / ₁₆ inch 233 mm
			UFMP-509	9 ⁷ / ₁₆ inch 233 mm
DP-1810	17 ⁷ / ₁₆ inch 443 mm	9 ¹⁵ / ₁₆ -10 ⁵ / ₁₆ inch 252-262 mm	UFLB-32	10 ³ / ₁₆ inch 259 mm
DP-1907	17 ¹⁵ / ₁₆ inch 456 mm	7-7 ⁷ / ₁₆ inch 178-188 mm	UFHT-33	6 ¹⁵ / ₁₆ inch 176 mm
			UFM-3315	6 ¹⁵ / ₁₆ inch 176 mm
DP-1909	17 ¹⁵ / ₁₆ inch 456 mm	9 ⁵ / ₁₆ -9 ³ / ₄ inch 237-247 mm	UFHT-34	9 ¹ / ₄ inch 235 mm
			UFM-3415	9 ¹ / ₄ inch 235 mm
DP-1911	17 ⁹ / ₁₆ inch 446 mm	11 ¹ / ₈ -11 ¹ / ₂ inch 282-292 mm	UFDP-56	11 inch 280 mm
			UFTB-55	11 inch 280 mm
			UFMP-511	11 inch 280 mm
DP-1912	17 ¹⁵ / ₁₆ inch 456 mm	11 ⁵ / ₈ -12 inch 295-305 mm	UFHT-35	11 ⁹ / ₁₆ inch 293 mm
			UFM-3515	11 ⁹ / ₁₆ inch 293 mm
DP-2107	20 ⁷ / ₈ inch 530 mm	6 ¹¹ / ₁₆ -7 ¹ / ₈ inch 170-180 mm	UF-432	6 ⁵ / ₈ inch 168 mm
			UF-423	6 ⁵ / ₈ inch 168 mm
DP-2109	20 ⁷ / ₈ inch 530 mm	8 ¹⁵ / ₁₆ -9 ⁵ / ₁₆ inch 227-237 mm	UF-442	8 ⁷ / ₈ inch 225 mm
DP-2110	20 ⁷ / ₈ inch 530 mm	9 ⁷ / ₁₆ -9 ¹³ / ₁₆ inch 239-249 mm	UF-433	9 ³ / ₈ inch 238 mm
DP-2111	20 ⁷ / ₈ inch 530 mm	11 ¹ / ₈ -11 ¹ / ₂ inch 282-292 mm	UF-452	11 inch 280 mm
DP-2112	21 ³ / ₈ inch 542 mm	12 ³ / ₄ -13 ¹ / ₈ inch 323-333 mm	UF-443	12 ⁵ / ₈ inch 321 mm
DP-2113	20 ⁷ / ₈ inch 530 mm	12 ³ / ₄ -13 ¹ / ₈ inch 323-333 mm	UF-462	12 ⁵ / ₈ inch 321 mm
DP-2207	21 ¹ / ₈ inch 536 mm	6 ³ / ₄ -7 ¹ / ₈ inch 171-181 mm	UFM-4305	6 ⁵ / ₈ inch 169 mm
DP-2209	21 ³ / ₁₆ inch 538 mm	9 ¹ / ₄ -10 ¹ / ₁₆ inch 235-255 mm	UFDP-65	9 ⁷ / ₁₆ inch 233 mm
			UFMP-609	9 ⁷ / ₁₆ inch 233 mm
		8 ⁷ / ₈ -9 ⁵ / ₁₆ inch 226-236 mm	UFTB-64	8 ¹³ / ₁₆ inch 224 mm
			UFM-4405	8 ¹³ / ₁₆ inch 224 mm
DP-2211	21 ³ / ₁₆ inch 538 mm	11 ¹ / ₈ -11 ¹ / ₂ inch 282-292 mm	UFDP-66	11 inch 280 mm
			UFTB-65	11 inch 280 mm
			UFMP-611	11 inch 280 mm
			UFM-4505	11 inch 280 mm
DP-2410	23 ⁵ / ₈ inch 600 mm	9 ¹⁵ / ₁₆ -10 ³ / ₈ inch 253-263 mm	UFLB-42	10 ³ / ₁₆ inch 259 mm
DP-2609	24 ³ / ₁₆ inch 631 mm	9 ¹ / ₄ -9 ³ / ₄ inch 235-245 mm	UFMP-709	9 ⁷ / ₁₆ inch 233 mm
DP-2611	24 ¹³ / ₁₆ inch 631 mm	11 ¹ / ₈ -11 ¹ / ₂ inch 282-292 mm	UFMP-711	11 inch 280 mm
DP-2707	25 ¹⁵ / ₁₆ inch 659 mm	6 ¹¹ / ₁₆ -7 ¹ / ₈ inch 170-180 mm	UF-532	6 ⁵ / ₈ inch 168 mm
DP-2709	25 ¹⁵ / ₁₆ inch 659 mm	8 ¹⁵ / ₁₆ -9 ⁵ / ₁₆ inch 227-237 mm	UF-542	8 ⁷ / ₈ inch 225 mm
DP-2710	25 ¹⁵ / ₁₆ inch 659 mm	9 ⁹ / ₁₆ -9 ¹⁵ / ₁₆ inch 243-253 mm	UF-533	9 ¹ / ₂ inch 241 mm
DP-2711	25 ¹⁵ / ₁₆ inch 659 mm	11 ¹ / ₈ -11 ¹ / ₂ inch 282-292 mm	UF-552	11 inch 280 mm
DP-2713	25 ¹⁵ / ₁₆ inch 659 mm	12 ³ / ₄ -13 ¹ / ₈ inch 323-333 mm	UF-543	12 ⁵ / ₈ inch 321 mm

*Protect your samples from unauthorized access using these security lock devices.


*No tools are required for installing lock devices.

*Locking devices are made of durable stainless steel.

(Continued on next page)

Security Lock Devices for Upright Freezer Drawer Racks

Security lock devices for upright freezer drawer racks are made of high quality stainless steel and secured with a lock included with each system. Select your lock device by identifying your specific upright freezer drawer rack and then the security lock device for your particular freezer rack.



Lock Device	Lock Dimension		Applicable to Rack	Rack Height
	D	H		
DD-03	5 ¹¹ / ₁₆ inch 144 mm	3 ¹¹ / ₁₆ inch 94 mm	UFD-B03-2	6 ¹ / ₁₆ inch 154 mm
			UFD-T03-2	6 ¹ / ₁₆ inch 154 mm
			UFD-LB03-2	6 ¹ / ₁₆ inch 154 mm
			UFD-LT03-2	6 ¹ / ₁₆ inch 154 mm
			UFD-T01-2	6 ¹ / ₁₆ inch 154 mm
			UFD-LT01-2	6 ¹ / ₁₆ inch 154 mm
DD-04	5 ¹¹ / ₁₆ inch 144 mm	4 ³ / ₁₆ inch 106 mm	UFD-B10-1	4 ¹ / ₂ inch 114 mm
			UFD-T10-1	4 ¹ / ₂ inch 114 mm
			UFD-LB10-1	4 ¹ / ₂ inch 114 mm
			UFD-LT10-1	4 ¹ / ₂ inch 114 mm
			UFD-T06-1	4 ¹ / ₂ inch 114 mm
			UFD-LT06-1	4 ¹ / ₂ inch 114 mm
			UFD-B-1	4 ¹⁵ / ₁₆ inch 126 mm
			UFD-T15-1	4 ¹⁵ / ₁₆ inch 126 mm
			UFD-T50-1	4 ¹⁵ / ₁₆ inch 126 mm
			UFD-LB-1	4 ¹⁵ / ₁₆ inch 126 mm
			UFD-LT15-1	4 ¹⁵ / ₁₆ inch 126 mm
			UFD-LT50-1	4 ¹⁵ / ₁₆ inch 126 mm
			UFD-XLB-1	4 ¹⁵ / ₁₆ inch 126 mm
			UFD-XLT15-1	4 ¹⁵ / ₁₆ inch 126 mm
UFD-XLT50-1	4 ¹⁵ / ₁₆ inch 126 mm			
DD-05	5 ¹¹ / ₁₆ inch 144 mm	5 ¹ / ₂ inch 140 mm	UFDDP-33	6 ¹ / ₈ inch 154 mm
			UFDDP-43	6 ¹ / ₈ inch 154 mm
			UFDDP-53	6 ¹ / ₈ inch 154 mm
			UFDDP-63	6 ¹ / ₈ inch 154 mm
			UFDDP-73	6 ¹ / ₈ inch 154 mm
DD-06	5 ¹¹ / ₁₆ inch 144 mm	6 ³ / ₈ inch 161 mm	UFD-332	7 ¹ / ₈ inch 181 mm
			UFD-432	7 ¹ / ₈ inch 181 mm
			UFD-532	7 ¹ / ₈ inch 181 mm
			UFD-323	6 ⁷ / ₁₆ inch 164 mm
			UFD-423	6 ⁷ / ₁₆ inch 164 mm
			UFD-B10-2	9 inch 228 mm
			UFD-T10-2	9 inch 228 mm
			UFD-LB10-2	9 inch 228 mm
			UFD-LT10-2	9 inch 228 mm
			UFD-T06-2	9 inch 228 mm
			UFD-LT06-2	9 inch 228 mm
			UFD-B03-3	9 ¹ / ₈ inch 231 mm
			UFD-T03-3	9 ¹ / ₈ inch 231 mm
			UFD-LB03-3	9 ¹ / ₈ inch 231 mm
			UFD-LT03-3	9 ¹ / ₈ inch 231 mm
			UFD-T01-3	9 ¹ / ₈ inch 231 mm
UFD-LT01-3	9 ¹ / ₈ inch 231 mm			

*Protect your samples from unauthorized access using these security lock devices.

*No tools are required for installing lock devices.

*Locking devices are made of durable stainless steel.

(Continued on next page)

Security Lock Devices for Upright Freezer Drawer Racks

Lock Device	Lock Dimension		Applicable to Rack	Rack Height
	D	H		
DD-06	5 ¹¹ / ₁₆ inch 144 mm	6 ³ / ₈ inch 161 mm	UFDMP-33	7 ⁵ / ₁₆ inch 186 mm
			UFDMP-43	7 ⁵ / ₁₆ inch 186 mm
			UFDMP-53	7 ⁵ / ₁₆ inch 186 mm
			UFDMP-63	7 ⁵ / ₁₆ inch 186 mm
			UFDMP-73	7 ⁵ / ₁₆ inch 186 mm
DD-07	5 ¹¹ / ₁₆ inch 144 mm	7 ⁵ / ₁₆ inch 186 mm	UFDDP-34	8 ¹ / ₁₆ inch 205 mm
			UFDDP-44	8 ¹ / ₁₆ inch 205 mm
			UFDDP-54	8 ¹ / ₁₆ inch 205 mm
			UFDDP-64	8 ¹ / ₁₆ inch 205 mm
			UFDDP-74	8 ¹ / ₁₆ inch 205 mm
DD-09	5 ¹¹ / ₁₆ inch 144 mm	8 ¹¹ / ₁₆ inch 220 mm	UFD-342	9 ⁷ / ₁₆ inch 240 mm
			UFD-442	9 ⁷ / ₁₆ inch 240 mm
			UFD-542	9 ⁷ / ₁₆ inch 240 mm
			UFD-333	10 ¹ / ₁₆ inch 256 mm
			UFD-433	10 ¹ / ₁₆ inch 256 mm
			UFD-533	10 ¹ / ₁₆ inch 256 mm
			UFD-B-2	9 ¹⁵ / ₁₆ inch 252 mm
			UFD-T15-2	9 ¹⁵ / ₁₆ inch 252 mm
			UFD-T50-2	9 ¹⁵ / ₁₆ inch 252 mm
			UFD-LB-2	9 ¹⁵ / ₁₆ inch 252 mm
			UFD-LT15-2	9 ¹⁵ / ₁₆ inch 252 mm
			UFD-LT50-2	9 ¹⁵ / ₁₆ inch 252 mm
			UFD-XLB-2	9 ¹⁵ / ₁₆ inch 252 mm
			UFD-XLT15-2	9 ¹⁵ / ₁₆ inch 252 mm
			UFD-XLT50-2	9 ¹¹ / ₁₆ inch 246 mm
			UFD-B-222	9 ¹¹ / ₁₆ inch 246 mm
			UFD-B-322	9 ¹¹ / ₁₆ inch 246 mm
			UFD-LB-422	9 ¹¹ / ₁₆ inch 246 mm
			UFD-T15-222	9 ¹¹ / ₁₆ inch 246 mm
			UFD-T15-322	9 ¹¹ / ₁₆ inch 246 mm
			UFD-LT15-422	9 ¹¹ / ₁₆ inch 246 mm
			UFD-T50-222	9 ¹¹ / ₁₆ inch 246 mm
			UFD-T50-322	9 ¹¹ / ₁₆ inch 246 mm
			UFD-LT50-422	9 ¹¹ / ₁₆ inch 246 mm
			UFDMP-34	9 ³ / ₄ inch 247 mm
			UFDMP-44	9 ³ / ₄ inch 247 mm
			UFDMP-54	9 ³ / ₄ inch 247 mm
UFDMP-64	9 ³ / ₄ inch 247 mm			
UFDMP-74	9 ³ / ₄ inch 247 mm			
DD-10	5 ¹¹ / ₁₆ inch 144 mm	9 ¹ / ₈ inch 232 mm	UFDDP-35	10 ¹ / ₁₆ inch 256 mm
			UFDDP-45	10 ¹ / ₁₆ inch 256 mm
			UFDDP-55	10 ¹ / ₁₆ inch 256 mm
			UFDDP-65	10 ¹ / ₁₆ inch 256 mm
			UFDDP-75	10 ¹ / ₁₆ inch 256 mm
DD-11	5 ¹¹ / ₁₆ inch 144 mm	11 inch 279 mm	UFD-352	11 ¹³ / ₁₆ inch 300mm
			UFD-452	11 ¹³ / ₁₆ inch 300mm
			UFDMP-35	12 ¹ / ₈ inch 308mm
			UFDMP-45	12 ¹ / ₈ inch 308mm
			UFDMP-55	12 ¹ / ₈ inch 308mm
			UFDMP-65	12 ¹ / ₈ inch 308mm
DDHT-06	6 ³ / ₈ inch 168 mm	6 ³ / ₈ inch 161 mm	UFDHT-23	7 ⁵ / ₁₆ inch 186 mm
			UFDHT-33	7 ⁵ / ₁₆ inch 186 mm
			UFDHT-43	7 ⁵ / ₁₆ inch 186 mm
DDHT-09	6 ³ / ₈ inch 168 mm	8 ¹¹ / ₁₆ inch 220 mm	UFDHT-24	9 ³ / ₄ inch 247 mm
			UFDHT-34	9 ³ / ₄ inch 247 mm
			UFDHT-44	9 ³ / ₄ inch 247 mm
DDHT-11	6 ³ / ₈ inch 168 mm	11 inch 279 mm	UFDHT-25	12 ¹ / ₈ inch 308mm
			UFDHT-35	12 ¹ / ₈ inch 308mm

*Protect your samples from unauthorized access using these security lock devices.

*No tools are required for installing lock devices.

*Locking devices are made of durable stainless steel.

(Continued on next page)