




Upright Freezer Racks

For 25-Place Slide Boxes

- ◆ Made of corrosion resistant stainless steel.
- ◆ Security lock devices are available for each rack.



L1 x W1 x H1: 5.5" x 3.6" x 1.4" / 140 x 91 x 35mm




L1 x W1 x H1: 5.5" x 3.3" x 1.4" / 140 x 83 x 34mm

Part Number	Dimension L x H x W	Compartment Size l x h x w	Capacity	Lock Device Option*	Shipping Weight(lb/Kg)
UFS-25-34	11 ⁵ / ₈ x 6 ³ / ₁₆ x 5 ¹¹ / ₁₆ inch 295.5 x 156.5 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	12	DP25-1006	4.0 / 1.83
UFS-25-44	15 ⁵ / ₁₆ x 6 ³ / ₁₆ x 5 ¹¹ / ₁₆ inch 389 x 156.5 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	16	DP25-1406	5.2 / 2.35
UFS-25-54	19 x 6 ³ / ₁₆ x 5 ¹¹ / ₁₆ inch 483 x 156.5 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	20	DP25-1806	6.4 / 2.92
UFS-25-64	22 ¹¹ / ₁₆ x 6 ³ / ₁₆ x 5 ¹¹ / ₁₆ inch 576.5 x 156.5 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	24	DP25-2106	7.7 / 3.47
UFS-25-74	26 ³ / ₈ x 6 ³ / ₁₆ x 5 ¹¹ / ₁₆ inch 670 x 156.5 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	28	DP25-2506	8.8 / 4.00
UFS-25-35	11 ⁵ / ₈ x 7 ¹¹ / ₁₆ x 5 ¹¹ / ₁₆ inch 295.5 x 195 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	15	DP25-1008	4.8 / 2.16
UFS-25-45	15 ⁵ / ₁₆ x 7 ¹¹ / ₁₆ x 5 ¹¹ / ₁₆ inch 389 x 195 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	20	DP25-1408	6.2 / 2.82
UFS-25-55	19 x 7 ¹¹ / ₁₆ x 5 ¹¹ / ₁₆ inch 483 x 195 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	25	DP25-1808	7.7 / 3.50
UFS-25-65	22 ¹¹ / ₁₆ x 7 ¹¹ / ₁₆ x 5 ¹¹ / ₁₆ inch 576.5 x 195 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	30	DP25-2108	9.1 / 4.12
UFS-25-75	26 ³ / ₈ x 7 ¹¹ / ₁₆ x 5 ¹¹ / ₁₆ inch 670 x 195 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	35	DP25-2508	10.4 / 4.74
UFS-25-36	11 ⁵ / ₈ x 9 ³ / ₁₆ x 5 ¹¹ / ₁₆ inch 295.5 x 233.5 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	18	DP25-1009	5.6 / 2.52
UFS-25-46	15 ⁵ / ₁₆ x 9 ³ / ₁₆ x 5 ¹¹ / ₁₆ inch 389 x 233.5 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	24	DP25-1409	7.2 / 3.28
UFS-25-56	19 x 9 ³ / ₁₆ x 5 ¹¹ / ₁₆ inch 483 x 233.5 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	30	DP25-1809	8.9 / 4.03
UFS-25-66	22 ¹¹ / ₁₆ x 9 ³ / ₁₆ x 5 ¹¹ / ₁₆ inch 576.5 x 233.5 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	36	DP25-2109	10.6 / 4.80
UFS-25-76	26 ³ / ₈ x 9 ³ / ₁₆ x 5 ¹¹ / ₁₆ inch 670 x 233.5 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	42	DP25-2509	12.2 / 5.53
UFS-25-37	11 ⁵ / ₈ x 10 ¹¹ / ₁₆ x 5 ¹¹ / ₁₆ inch 295.5 x 272 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	21	DP25-1011	6.3 / 2.87
UFS-25-47	15 ⁵ / ₁₆ x 10 ¹¹ / ₁₆ x 5 ¹¹ / ₁₆ inch 389 x 272 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	28	DP25-1411	8.2 / 3.72
UFS-25-57	19 x 10 ¹¹ / ₁₆ x 5 ¹¹ / ₁₆ inch 483 x 272 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	35	DP25-1811	10.1 / 4.57
UFS-25-67	22 ¹¹ / ₁₆ x 10 ¹¹ / ₁₆ x 5 ¹¹ / ₁₆ inch 576.5 x 272 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	42	DP25-2111	12.0 / 5.43
UFS-25-77	26 ³ / ₈ x 10 ¹¹ / ₁₆ x 5 ¹¹ / ₁₆ inch 670 x 272 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	49	DP25-2511	13.9 / 6.29
UFS-25-38	11 ⁵ / ₈ x 12 ¹ / ₄ x 5 ¹¹ / ₁₆ inch 295.5 x 310.5 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	24	DP25-1012	7.1 / 3.24
UFS-25-48	15 ⁵ / ₁₆ x 12 ¹ / ₄ x 5 ¹¹ / ₁₆ inch 389 x 310.5 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	32	DP25-1412	9.3 / 4.20
UFS-25-58	19 x 12 ¹ / ₄ x 5 ¹¹ / ₁₆ inch 483 x 310.5 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	40	DP25-1812	11.4 / 5.17
UFS-25-68	22 ¹¹ / ₁₆ x 12 ¹ / ₄ x 5 ¹¹ / ₁₆ inch 576.5 x 310.5 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	48	DP25-2112	13.5 / 6.14
UFS-25-78	26 ³ / ₈ x 12 ¹ / ₄ x 5 ¹¹ / ₁₆ inch 670 x 310.5 x 145 mm	3 ⁵ / ₈ x 1 ¹ / ₂ x 5 ⁵ / ₈ inch 92 x 38 x 143 mm	56	DP25-2512	15.7 / 7.12

*Security lock device is ordered separately.


Upright Freezer Drawer Rack




L1 x W1 x H1: 5.5" x 3.6" x 1.4"
140 x 91 x 35mm

For 25-Place Slide Boxes

- ◆ Made of corrosion resistant stainless steel.
- ◆ Drawer racks allow easy access to and storage of samples without having to remove the whole rack from the freezer.
- ◆ Security lock devices are available for each rack.



L1 x W1 x H1: 5.5" x 3.3" x 1.4"
140 x 83 x 34mm



Part Number	Dimension L x H x W	Compartment Size l x h x w	Capacity	Lock Device Option*	Shipping Weight(lb/Kg)
UFDS-25-34	11 ¹ / ₂ x 6 ⁷ / ₈ x 6 inch 291.5 x 175 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	12	DD25-06	6.6 / 3.00
UFDS-25-44	15 ³ / ₁₆ x 6 ⁷ / ₈ x 6 inch 385 x 175 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	16	DD25-06	8.4 / 3.81
UFDS-25-54	18 ⁷ / ₈ x 6 ⁷ / ₈ x 6 inch 479 x 175 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	20	DD25-06	10.2 / 4.61
UFDS-25-64	22 ⁹ / ₁₆ x 6 ⁷ / ₈ x 6 inch 572.5 x 175 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	24	DD25-06	12.1 / 5.51
UFDS-25-74	26 ¹ / ₄ x 6 ⁷ / ₈ x 6 inch 666 x 175 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	28	DD25-06	13.9 / 6.30
UFDS-25-35	11 ¹ / ₂ x 8 ⁹ / ₁₆ x 6 inch 291.5 x 218 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	15	DD25-08	7.9 / 3.60
UFDS-25-45	15 ³ / ₁₆ x 8 ⁹ / ₁₆ x 6 inch 385 x 218 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	20	DD25-08	10.4 / 4.70
UFDS-25-55	18 ⁷ / ₈ x 8 ⁹ / ₁₆ x 6 inch 479 x 218 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	25	DD25-08	12.6 / 5.70
UFDS-25-65	22 ⁹ / ₁₆ x 8 ⁹ / ₁₆ x 6 inch 572.5 x 218 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	30	DD25-08	14.9 / 6.76
UFDS-25-75	26 ¹ / ₄ x 8 ⁹ / ₁₆ x 6 inch 666 x 218 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	35	DD25-08	17.1 / 7.74
UFDS-25-36	11 ¹ / ₂ x 10 ³ / ₁₆ x 6 inch 291.5 x 259 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	18	DD25-09	9.4 / 4.27
UFDS-25-46	15 ³ / ₁₆ x 10 ³ / ₁₆ x 6 inch 385 x 259 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	24	DD25-09	12.2 / 5.54
UFDS-25-56	18 ⁷ / ₈ x 10 ³ / ₁₆ x 6 inch 479 x 259 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	30	DD25-09	14.9 / 6.74
UFDS-25-66	22 ⁹ / ₁₆ x 10 ³ / ₁₆ x 6 inch 572.5 x 259 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	36	DD25-09	17.6 / 8.00
UFDS-25-76	26 ¹ / ₄ x 10 ³ / ₁₆ x 6 inch 666 x 259 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	42	DD25-09	20.2 / 9.16
UFDS-25-37	11 ¹ / ₂ x 12 x 6 inch 291.5 x 304.5 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	21	DD25-11	11.2 / 5.06
UFDS-25-47	15 ³ / ₁₆ x 12 x 6 inch 385 x 304.5 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	28	DD25-11	14.1 / 6.40
UFDS-25-57	18 ⁷ / ₈ x 12 x 6 inch 479 x 304.5 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	35	DD25-11	17.2 / 7.80
UFDS-25-67	22 ⁹ / ₁₆ x 12 x 6 inch 572.5 x 304.5 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	42	DD25-11	20.4 / 9.24
UFDS-25-77	26 ¹ / ₄ x 12 x 6 inch 666 x 304.5 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	49	DD25-11	23.4 / 10.60
UFDS-25-38	11 ¹ / ₂ x 13 ¹¹ / ₁₆ x 6 inch 291.5 x 347.5 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	24	DD25-13	12.7 / 5.74
UFDS-25-48	15 ³ / ₁₆ x 13 ¹¹ / ₁₆ x 6 inch 385 x 347.5 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	32	DD25-13	16.0 / 7.24
UFDS-25-58	18 ⁷ / ₈ x 13 ¹¹ / ₁₆ x 6 inch 479 x 347.5 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	40	DD25-13	19.5 / 8.83
UFDS-25-68	22 ⁹ / ₁₆ x 13 ¹¹ / ₁₆ x 6 inch 572.5 x 347.5 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	48	DD25-13	23.1 / 10.50
UFDS-25-78	26 ¹ / ₄ x 13 ¹¹ / ₁₆ x 6 inch 666 x 347.5 x 152 mm	3 ¹¹ / ₁₆ x 1 ¹ / ₂ x 5 ⁹ / ₁₆ inch 92.5 x 38 x 141 mm	56	DD25-13	26.5 / 12.00

*Security lock device is ordered separately.