



Upright Freezer Racks

For 96-Well & 384-Well

Microtiter Plates



- ◆ Made of corrosion resistant stainless steel.
- ◆ Locking rods are included.
- ◆ Security lock devices are available for each rack.




*UFML-311L / UFML-411L / UFML-511L / UFML-611L / UFML-711L (k:138mm) UFML-312L / UFML-412L / UFML-512L / UFML-612L / UFML-712L (k:151mm)
 UFML-316L / UFML-416L / UFML-516L / UFML-616L / UFML-716L (k:202mm) UFML-319L / UFML-419L / UFML-519L / UFML-619L / UFML-719L (k:239mm)

Part Number	Dimension L x H x W	Compartment Size l x h x w	Capacity (Plate)	Lock Device Option*	Shipping Weight (lb/Kg)
UFMP-309L	11 1/2 x 9 3/16 x 5 1/2 inch 292 x 233 x 140 mm	3 9/16 x 4 1/2 x 5 3/16 inch 91 x 115 x 130.5 mm	42 Plates w/lids 48 Plates w/o lids	DP-1109	3.8 / 1.71
UFMP-311L	11 1/2 x 11 x 5 1/2 inch 292 x 280 x 140 mm	3 9/16 x 5 7/16 x 5 3/16 inch 91 x 138 x 130.5 mm	51 Plates w/lids 54 Plates w/o lids	DP-1111	4.2 / 1.89
UFMP-312L	11 1/2 x 12 x 5 1/2 inch 292 x 305 x 140 mm	3 9/16 x 5 15/16 x 5 3/16 inch 91 x 151 x 130.5 mm	54 Plates w/lids 60 Plates w/o lids	DP-1112	4.4 / 1.99
UFMP-316L	11 1/2 x 16 x 5 1/2 inch 292 x 407 x 140 mm	3 9/16 x 7 15/16 x 5 3/16 inch 91 x 202 x 130.5 mm	72 Plates w/lids 84 Plates w/o lids	DP-1116	5.3 / 2.40
UFMP-319L	11 1/2 x 18 15/16 x 5 1/2 inch 292 x 481 x 140 mm	3 9/16 x 9 7/16 x 5 3/16 inch 91 x 239 x 130.5 mm	84 Plates w/lids 96 Plates w/o lids	DP-1119	5.9 / 2.69
UFMP-409L	15 3/16 x 9 3/16 x 5 1/2 inch 386 x 233 x 140 mm	3 9/16 x 4 1/2 x 5 3/16 inch 91 x 115 x 130.5 mm	56 Plates w/lids 64 Plates w/o lids	DP-1409	4.9 / 2.20
UFMP-411L	15 3/16 x 11 x 5 1/2 inch 386 x 280 x 140 mm	3 9/16 x 5 7/16 x 5 3/16 inch 91 x 138 x 130.5 mm	64 Plates w/lids 80 Plates w/o lids	DP-1411	5.4 / 2.43
UFMP-412L	15 3/16 x 12 x 5 1/2 inch 386 x 305 x 140 mm	3 9/16 x 5 15/16 x 5 3/16 inch 91 x 151 x 130.5 mm	72 Plates w/lids 84 Plates w/o lids	DP-1412	5.6 / 2.55
UFMP-416L	15 3/16 x 16 x 5 1/2 inch 386 x 407 x 140 mm	3 9/16 x 7 15/16 x 5 3/16 inch 91 x 202 x 130.5 mm	96 Plates w/lids 112 Plates w/o lids	DP-1416	6.7 / 3.06
UFMP-419L	15 3/16 x 18 15/16 x 5 1/2 inch 386 x 481 x 140 mm	3 9/16 x 9 7/16 x 5 3/16 inch 91 x 239 x 130.5 mm	116 Plates w/lids 136 Plates w/o lids	DP-1419	7.5 / 3.42
UFMP-509L	18 13/16 x 9 3/16 x 5 1/2 inch 478 x 233 x 140 mm	3 9/16 x 4 1/2 x 5 3/16 inch 91 x 115 x 130.5 mm	70 Plates w/lids 80 Plates w/o lids	DP-1809	6.1 / 2.77
UFMP-511L	18 13/16 x 11 x 5 1/2 inch 478 x 280 x 140 mm	3 9/16 x 5 7/16 x 5 3/16 inch 91 x 138 x 130.5 mm	85 Plates w/lids 95 Plates w/o lids	DP-1811	6.7 / 3.04
UFMP-512L	18 13/16 x 12 x 5 1/2 inch 478 x 305 x 140 mm	3 9/16 x 5 15/16 x 5 3/16 inch 91 x 151 x 130.5 mm	90 Plates w/lids 100 Plates w/o lids	DP-1812	6.9 / 3.11
UFMP-516L	18 13/16 x 16 x 5 1/2 inch 478 x 407 x 140 mm	3 9/16 x 7 15/16 x 5 3/16 inch 91 x 202 x 130.5 mm	120 Plates w/lids 140 Plates w/o lids	DP-1816	8.2 / 3.72
UFMP-519L	18 13/16 x 18 15/16 x 5 1/2 inch 478 x 481 x 140 mm	3 9/16 x 9 7/16 x 5 3/16 inch 91 x 239 x 130.5 mm	140 Plates w/lids 170 Plates w/o lids	DP-1819	9.2 / 4.16
UFMP-609L	22 1/2 x 9 3/16 x 5 1/2 inch 571 x 233 x 140 mm	3 9/16 x 4 1/2 x 5 3/16 inch 91 x 115 x 130.5 mm	84 Plates w/lids 96 Plates w/o lids	DP-2209	7.2 / 3.27
UFMP-611L	22 1/2 x 11 x 5 1/2 inch 571 x 280 x 140 mm	3 9/16 x 5 7/16 x 5 3/16 inch 91 x 138 x 130.5 mm	96 Plates w/lids 120 Plates w/o lids	DP-2211	7.9 / 3.59
UFMP-612L	22 1/2 x 12 x 5 1/2 inch 571 x 305 x 140 mm	3 9/16 x 5 15/16 x 5 3/16 inch 91 x 151 x 130.5 mm	108 Plates w/lids 132 Plates w/o lids	DP-2212	8.1 / 3.67
UFMP-616L	22 1/2 x 16 x 5 1/2 inch 571 x 407 x 140 mm	3 9/16 x 7 15/16 x 5 3/16 inch 91 x 202 x 130.5 mm	144 Plates w/lids 168 Plates w/o lids	DP-2216	9.7 / 4.38
UFMP-619L	22 1/2 x 18 15/16 x 5 1/2 inch 571 x 481 x 140 mm	3 9/16 x 9 7/16 x 5 3/16 inch 91 x 239 x 130.5 mm	168 Plates w/lids 204 Plates w/o lids	DP-2219	10.8 / 4.90
UFMP-709L	26 1/8 x 9 3/16 x 5 1/2 inch 664 x 233 x 140 mm	3 9/16 x 4 1/2 x 5 3/16 inch 91 x 115 x 130.5 mm	98 Plates w/lids 112 Plates w/o lids	DP-2609	8.3 / 3.77
UFMP-711L	26 1/8 x 11 x 5 1/2 inch 664 x 280 x 140 mm	3 9/16 x 5 7/16 x 5 3/16 inch 91 x 138 x 130.5 mm	112 Plates w/lids 140 Plates w/o lids	DP-2611	9.1 / 4.13
UFMP-712L	26 1/8 x 12 x 5 1/2 inch 664 x 305 x 140 mm	3 9/16 x 5 15/16 x 5 3/16 inch 91 x 151 x 130.5 mm	126 Plates w/lids 154 Plates w/o lids	DP-2612	9.3 / 4.22
UFMP-716L	26 1/8 x 16 x 5 1/2 inch 664 x 407 x 140 mm	3 9/16 x 7 15/16 x 5 3/16 inch 91 x 202 x 130.5 mm	168 Plates w/lids 196 Plates w/o lids	DP-2616	11.1 / 5.05
UFMP-719L	26 1/8 x 18 15/16 x 5 1/2 inch 664 x 481 x 140 mm	3 9/16 x 9 7/16 x 5 3/16 inch 91 x 239 x 130.5 mm	196 Plates w/lids 238 Plates w/o lids	DP-2619	12.4 / 5.64


* Lock device is ordered separately.

Upright Freezer Drawer Racks

For 96-Well & 384-Well Microtiter Plates



- ◆ Made of corrosion resistant stainless steel.
- ◆ Drawer racks allow easy access to fewer microtiter plates at a time thus protecting other samples in the rack from warming up.
- ◆ Security lock devices are available for each rack.




Part Number	Description	Dimension L x H x W	Compartment Size l x h x w	Capacity (Plate)	Lock Device Option*	Shipping Weight (lb/Kg)
UFDMP-33	Rack w/3 drawers (3 compartments per drawer) for Microtiter Plates	11 ¹ / ₈ x 7 ⁷ / ₁₆ x 5 ⁷ / ₁₆ inch 282 x 186 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	27 w/lids 36w/o lids	DD-06	5.4 / 2.45
UFDMP-34	Rack w/4 drawers (3 compartments per drawer) for Microtiter Plates	11 ¹ / ₈ x 9 ³ / ₄ x 5 ⁷ / ₁₆ inch 282 x 247 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	36 w/lids 48w/o lids	DD-09	6.9 / 3.13
UFDMP-35	Rack w/5 drawers (3 compartments per drawer) for Microtiter Plates	11 ¹ / ₈ x 12 ¹ / ₈ x 5 ⁷ / ₁₆ inch 282 x 308 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	45 w/lids 60w/o lids	DD-11	8.6 / 3.90
UFDMP-36	Rack w/6 drawers (3 compartments per drawer) for Microtiter Plates	11 ¹ / ₈ x 14 ¹ / ₂ x 5 ⁷ / ₁₆ inch 282 x 368 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	54 w/lids 72 w/o lids	DD-14	10.2 / 4.64
UFDMP-43	Rack w/3 drawers (4 compartments per drawer) for Microtiter Plates	14 ⁵ / ₈ x 7 ⁷ / ₁₆ x 5 ⁷ / ₁₆ inch 372 x 186 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	36 w/lids 48w/o lids	DD-06	6.8 / 3.10
UFDMP-44	Rack w/4 drawers (4 compartments per drawer) for Microtiter Plates	14 ⁵ / ₈ x 9 ³ / ₄ x 5 ⁷ / ₁₆ inch 372 x 247 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	48 w/lids 64 w/o lids	DD-09	8.9 / 4.02
UFDMP-45	Rack w/5 drawers (4 compartments per drawer) for Microtiter Plates	14 ⁵ / ₈ x 12 ¹ / ₈ x 5 ⁷ / ₁₆ inch 372 x 308 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	60 w/lids 80 w/o lids	DD-11	10.9 / 4.94
UFDMP-46	Rack w/6 drawers (4 compartments per drawer) for Microtiter Plates	14 ⁵ / ₈ x 14 ¹ / ₂ x 5 ⁷ / ₁₆ inch 372 x 368 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	72 w/lids 96 w/o lids	DD-14	12.9 / 5.86
UFDMP-53	Rack w/3 drawers (5 compartments per drawer) for Microtiter Plates	18 ³ / ₁₆ x 7 ⁷ / ₁₆ x 5 ⁷ / ₁₆ inch 462 x 168 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	45 w/lids 60 w/o lids	DD-06	8.2 / 3.70
UFDMP-54	Rack w/4 drawers (5 compartments per drawer) for Microtiter Plates	18 ³ / ₁₆ x 9 ³ / ₄ x 5 ⁷ / ₁₆ inch 462 x 247 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	60 w/lids 80 w/o lids	DD-09	10.6 / 4.79
UFDMP-55	Rack w/5 drawers (5 compartments per drawer) for Microtiter Plates	18 ³ / ₁₆ x 12 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 462 x 308 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	75 w/lids 100 w/o lids	DD-11	13.0 / 5.90
UFDMP-56	Rack w/6 drawers (5 compartments per drawer) for Microtiter Plates	18 ³ / ₁₆ x 14 ¹ / ₂ x 5 ⁷ / ₁₆ inch 462 x 368 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	90 w/lids 120 w/o lids	DD-14	15.4 / 7.00
UFDMP-63	Rack w/3 drawers (6 compartments per drawer) for Microtiter Plates	21 ³ / ₄ x 7 ⁷ / ₁₆ x 5 ⁷ / ₁₆ inch 552 x 186 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	54 w/lids 72 w/o lids	DD-06	9.7 / 4.42
UFDMP-64	Rack w/4 drawers (6 compartments per drawer) for Microtiter Plates	21 ³ / ₄ x 9 ³ / ₄ x 5 ⁷ / ₁₆ inch 552 x 247 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	72 w/lids 96 w/o lids	DD-09	12.6 / 5.73
UFDMP-65	Rack w/5 drawers (6 compartments per drawer) for Microtiter Plates	21 ³ / ₄ x 12 ¹ / ₈ x 5 ⁷ / ₁₆ inch 552 x 308 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	90 w/lids 120 w/o lids	DD-11	15.5 / 7.05
UFDMP-66	Rack w/6 drawers (6 compartments per drawer) for Microtiter Plates	21 ³ / ₄ x 14 ¹ / ₂ x 5 ⁷ / ₁₆ inch 552 x 368 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	108 w/lids 144 w/o lids	DD-14	17.9 / 8.14
UFDMP-73	Rack w/3 drawers (7 compartments per drawer) for Microtiter Plates	25 ¹ / ₄ x 7 ⁷ / ₁₆ x 5 ⁷ / ₁₆ inch 642 x 186 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	63 w/lids 84 w/o lids	DD-06	11.0 / 4.97
UFDMP-74	Rack w/4 drawers (7 compartments per drawer) for Microtiter Plates	25 ¹ / ₄ x 9 ³ / ₄ x 5 ⁷ / ₁₆ inch 642 x 247 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	84 w/lids 112 w/o lids	DD-09	14.2 / 6.43
UFDMP-75	Rack w/5 drawers (7 compartments per drawer) for Microtiter Plates	25 ¹ / ₄ x 12 ¹ / ₈ x 5 ⁷ / ₁₆ inch 642 x 308 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	105 w/lids 140 w/o lids	DD-11	17.4 / 7.91
UFDMP-76	Rack w/6 drawers (7 compartments per drawer) for Microtiter Plates	25 ¹ / ₄ x 14 ¹ / ₂ x 5 ⁷ / ₁₆ inch 642 x 368 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	126 w/lids 168 w/o lids	DD-14	20.7 / 9.38


* Lock device is ordered separately.

Upright Freezer Racks

For 96 Deep-Well Microtiter Plates



- ◆ Made of corrosion resistant stainless steel.
- ◆ Locking rods are included.
- ◆ Security lock devices are available for each rack.




Part Number	Dimension L x H x W	Compartment Size l x h x w	Capacity (Plate)	Lock Device Option*	Shipping Weight(lb/Kg)
UFDP-34L	11 1/2 x 7 1/2 x 5 1/2 inch 292.5 x 191 x 140 mm	3 9/16 x 1 13/16 x 5 3/16 inch 91 x 45.7 x 130.5 mm	3 x 4 = 12		DP-1107 4.2 / 1.91
UFDP-35L	11 1/2 x 9 3/8 x 5 1/2 inch 292.5 x 237.5 x 140 mm	3 9/16 x 1 13/16 x 5 3/16 inch 91 x 45.7 x 130.5 mm	3 x 5 = 15		DP-1109 5.0 / 2.28
UFDP-36L	11 1/2 x 11 3/16 x 5 1/2 inch 292.5 x 283.5 x 140 mm	3 9/16 x 1 13/16 x 5 3/16 inch 91 x 45.7 x 130.5 mm	3 x 6 = 18		DP-1111 5.8 / 2.64
UFDP-44L	15 3/16 x 7 7/2 x 5 1/2 inch 385.3 x 191 x 140 mm	3 9/16 x 1 13/16 x 5 3/16 inch 91 x 45.7 x 130.5 mm	4 x 4 = 16		DP-1507 5.4 / 2.47
UFDP-45L	15 3/16 x 9 3/8 x 5 1/2 inch 385.3 x 237.5 x 140mm	3 9/16 x 1 13/16 x 5 3/16 inch 91 x 45.7 x 130.5 mm	4 x 5 = 20		DP-1509 6.5 / 2.95
UFDP-46L	15 3/16 x 11 3/16 x 5 1/2 inch 385.3 x 283.5 x 140 mm	3 9/16 x 1 13/16 x 5 3/16 inch 91 x 45.7 x 130.5 mm	4 x 6 = 24		DP-1511 7.5 / 3.42
UFDP-54L	18 13/16 x 7 1/2 x 5 3/16 inch 478 x 191 x 140 mm	3 9/16 x 1 13/16 x 5 3/16 inch 91 x 45.7 x 130.5 mm	5 x 4 = 20		DP-1807 6.7 / 3.03
UFDP-55L	18 13/16 x 9 3/16 x 5 1/2 inch 478 x 237 x 140 mm	3 9/16 x 1 13/16 x 5 3/16 inch 91 x 45.7 x 130.5 mm	5 x 5 = 25		DP-1809 8.2 / 3.72
UFDP-56L	18 13/16 x 11 3/16 x 5 1/2 inch 478 x 284 x 140mm	3 9/16 x 1 13/16 x 5 3/16 inch 91 x 45.7 x 130.5 mm	5 x 6 = 30		DP-1811 9.5 / 4.32
UFDP-64L	22 1/2 x 7 1/2 x 5 1/2 inch 571 x 191 x 140 mm	3 9/16 x 1 13/16 x 5 3/16 inch 91 x 45.7 x 130.5 mm	6 x 4 = 24		DP-2207 8.0 / 3.62
UFDP-65L	22 1/2 x 9 3/16 x 5 1/2 inch 571 x 237 x 140 mm	3 9/16 x 1 13/16 x 5 3/16 inch 91 x 45.7 x 130.5 mm	6 x 5 = 30		DP-2209 9.8 / 4.43
UFDP-66L	22 1/2 x 11 3/16 x 5 1/2 inch 571 x 284 x 140 mm	3 9/16 x 1 13/16 x 5 3/16 inch 91 x 45.7 x 130.5 mm	6 x 6 = 36		DP-2211 11.3 / 5.13
UFDP-74L	26 1/8 x 7 1/2 x 5 1/2 inch 663.7 x 191 x 140 mm	3 9/16 x 1 13/16 x 5 3/16 inch 91 x 45.7 x 130.5 mm	7 x 4 = 28		DP-2507 9.2 / 4.17
UFDP-75L	26 1/8 x 9 3/8 x 5 1/2 inch 663.7 x 237.5 x 140 mm	3 9/16 x 1 13/16 x 5 3/16 inch 91 x 45.7 x 130.5 mm	7 x 5 = 35		DP-2509 11.3 / 5.11
UFDP-76L	26 1/8 x 11 3/16 x 5 1/2 inch 663.7 x 283.5 x 140 mm	3 9/16 x 1 13/16 x 5 3/16 inch 91 x 45.7 x 130.5 mm	7 x 6 = 42		DP-2511 13.1 / 5.93


* Lock device is ordered separately.

Upright Freezer Drawer Racks

For 96 Deep-Well Microtiter Plates



- ◆ Made of corrosion resistant stainless steel.
- ◆ Drawer racks allow easy access to fewer microtiter plates at a time thus protecting other samples in the rack from warming up.
- ◆ Security lock devices are available for each rack.



Part Number	Description	Dimension L x H x W	Compartment Size l x h x w	Capacity	Lock Device Option*	Shipping Weight(lb/Kg)
UFDDP-33	Rack w/3 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	11 ¹ / ₈ x 6 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 282 x 154 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	3x3=9	DD-05	8
UFDDP-34	Rack w/4 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	11 ¹ / ₈ x 8 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 282 x 205 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	3x4=12	DD-07	6.8 / 3.09
UFDDP-35	Rack w/5 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	11 ¹ / ₈ x 10 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 282 x 256 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	3x5=15	DR-10	8.4 / 3.79
UFDDP-36	Rack w/6 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	11 ¹ / ₈ x 12 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 282 x 307 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	3x6=18	DD-12	9.9 / 4.50
UFDDP-37	Rack w/7 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	11 ¹ / ₈ x 14 ¹ / ₈ x 5 ⁷ / ₁₆ inch 282 x 358 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	3x7=21	DD-14	11.7 / 5.30
UFDDP-43	Rack w/3 drawers (4 compartments per drawer) for 96Deep Well Microtiter Plates	14 ⁵ / ₈ x 6 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 372 x 154 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	4x3=12	DD-05	6.6 / 3.00
UFDDP-44	Rack w/4 drawers (4 compartments per drawer) for 96Deep Well Microtiter Plates	14 ⁵ / ₈ x 8 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 372 x 205 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	4x4=16	DD-07	8.6 / 3.89
UFDDP-45	Rack w/5 drawers (4 compartments per drawer) for 96Deep Well Microtiter Plates	14 ⁵ / ₈ x 10 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 372 x 256 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	4x5=20	DR-10	10.5 / 4.77
UFDDP-46	Rack w/6 drawers (4 compartments per drawer) for 96Deep Well Microtiter Plates	14 ⁵ / ₈ x 12 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 372 x 307 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	4x6=24	DD-12	12.5 / 5.66
UFDDP-47	Rack w/7 drawers (4 compartments per drawer) for 96Deep Well Microtiter Plates	14 ⁵ / ₈ x 14 ¹ / ₈ x 5 ⁷ / ₁₆ inch 372 x 358 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	4x7=28	DD-14	14.6 / 6.60
UFDDP-53	Rack w/3 drawers (5 compartments per drawer) for 96Deep Well Microtiter Plates	18 ³ / ₁₆ x 6 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 462 x 154 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	5x3=15	DD-05	8.0 / 3.61
UFDDP-54	Rack w/4 drawers (5 compartments per drawer) for 96Deep Well Microtiter Plates	18 ³ / ₁₆ x 8 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 462 x 205 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	5x4=20	DD-07	10.3 / 4.68
UFDDP-55	Rack w/5 drawers (5 compartments per drawer) for 96Deep Well Microtiter Plates	18 ³ / ₁₆ x 10 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 462 x 256 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	5x5=25	DR-10	12.7 / 5.75
UFDDP-56	Rack w/6 drawers (5 compartments per drawer) for 96Deep Well Microtiter Plates	18 ³ / ₁₆ x 12 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 462 x 307 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	5x6=30	DD-12	15.0 / 6.81
UFDDP-57	Rack w/7 drawers (5 compartments per drawer) for 96Deep Well Microtiter Plates	18 ³ / ₁₆ x 14 ¹ / ₈ x 5 ⁷ / ₁₆ inch 462 x 358 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	5x7=35	DD-14	17.4 / 7.90
UFDDP-63	Rack w/3 drawers (6 compartments per drawer) for 96Deep Well Microtiter Plates	21 ³ / ₄ x 6 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 552 x 154 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	6x3=18	DD-05	9.3 / 4.23
UFDDP-64	Rack w/4 drawers (6 compartments per drawer) for 96Deep Well Microtiter Plates	21 ³ / ₄ x 8 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 552 x 205 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	6x4=24	DD-07	12.1 / 5.48
UFDDP-65	Rack w/5 drawers (6 compartments per drawer) for 96Deep Well Microtiter Plates	21 ³ / ₄ x 10 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 552 x 256 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	6x5=30	DR-10	14.8 / 6.73
UFDDP-66	Rack w/6 drawers (6 compartments per drawer) for 96Deep Well Microtiter Plates	21 ³ / ₄ x 12 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 552 x 307 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	6x6=36	DD-12	17.6 / 7.97
UFDDP-67	Rack w/7 drawers (6 compartments per drawer) for 96Deep Well Microtiter Plates	21 ³ / ₄ x 14 ¹ / ₈ x 5 ⁷ / ₁₆ inch 552 x 358 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	6x7=42	DD-14	20.3 / 9.20
UFDDP-73	Rack w/3 drawers (7 compartments per drawer) for 96Deep Well Microtiter Plates	25 ³ / ₁₆ x 6 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 640 x 154 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	7x3=21	DD-05	10.7 / 4.84
UFDDP-74	Rack w/4 drawers (7 compartments per drawer) for 96 Deep Well Microtiter Plates	25 ³ / ₁₆ x 8 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 640 x 205 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	7x4=28	DD-07	13.8 / 6.27
UFDDP-75	Rack w/5 drawers (7 compartments per drawer) for 96 Deep Well Microtiter Plates	25 ³ / ₁₆ x 10 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 640 x 256 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	7x5=35	DR-10	17.0 / 7.70
UFDDP-76	Rack w/6 drawers (7 compartments per drawer) for 96 Deep Well Microtiter Plates	25 ³ / ₁₆ x 12 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 640 x 307 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	7x6=42	DD-12	20.1 / 9.12
UFDDP-77	Rack w/7 drawers (7 compartments per drawer) for 96 Deep Well Microtiter Plates	25 ³ / ₁₆ x 14 ¹ / ₈ x 5 ⁷ / ₁₆ inch 640 x 358 x 138 mm	3 ¹ / ₂ x 1 ³ / ₄ x 5 ¹ / ₈ inch 89.4 x 45.4 x 130 mm	7x7=49	DD-14	24.1 / 10.63

* Lock device is ordered separately.