

Upright Freezer Drawer Racks

For Blood Sample Tubes

- ◆ Racks are specially designed for storing blood sample tubes.
- ◆ Made of corrosion resistant stainless steel.
- ◆ Drawer Racks facilitate both access to and storage of samples without having to remove the whole rack from the freezer.
- ◆ Security lock devices are available for each rack.



Content to Store	Rack Length	Rack Height	Number of Drawers	Part Number	Description	Dimension L x H x W	Capacity	Lock Device Option*	Shipping Weight (lb/Kg)
8ml - 10ml	16 1/2 inch 419 mm	4 1/2 inch 114 mm	1	UFD-T10-1	Rack with 1 Drawer for 8 ml, 8.5ml, 10ml Tubes (Φ ≤ 5/8 inch H ≤ 4 3/16 inch)	16 11/16 x 4 1/2 x 5 1/2 inch 424 x 114 x 140 mm	60 Tubes of 8ml, 8.5ml or 10ml	DD-04	3.8 / 1.74
	22 inch 559 mm			UFD-LT10-1		22 3/16 x 4 1/2 x 5 1/2 inch 564 x 114 x 140 mm	80 Tubes of 8ml, 8.5ml or 10ml		4.7 / 2.14
	16 1/2 inch 419 mm	9 inch 228 mm	2	UFD-T10-2	Rack with 2 Drawers for 8 ml, 8.5ml, 10ml Tubes (Φ ≤ 5/8 inch H ≤ 4 3/16 inch)	16 11/16 x 9 x 5 1/2 inch 424 x 228 x 140 mm	120 Tubes of 8ml, 8.5ml or 10ml	DD-06	6.8 / 3.10
	22 inch 559 mm			UFD-LT10-2		22 3/16 x 9 x 5 1/2 inch 564 x 228 x 140 mm	160 Tubes of 8ml, 8.5ml or 10ml		8.4 / 3.79
4ml - 7ml	16 1/2 inch 419 mm	4 1/2 inch 114 mm	1	UFD-T06-1	Rack with 1 Drawer for 4 ml, 4.5ml, 6ml, 7ml Tubes (Φ ≤ 1/2 inch H ≤ 4 3/16 inch)	16 11/16 x 4 1/2 x 5 1/2 inch 424 x 114 x 140 mm	85 Tubes of 4ml, 4.5ml or 7ml	DD-04	3.6 / 1.64
	22 inch 559 mm			UFD-LT06-1		22 3/16 x 4 1/2 x 5 1/2 inch 564 x 114 x 140 mm	115 Tubes of 4ml, 4.5ml, 6ml or 7ml		4.6 / 2.09
	16 1/2 inch 419 mm	9 inch 228 mm	2	UFD-T06-2	Rack with 2 Drawers for 4 ml, 4.5ml, 6ml, 7ml Tubes (Φ ≤ 1/2 inch H ≤ 4 3/16 inch)	16 11/16 x 9 x 5 1/2 inch 424 x 228 x 140 mm	170 Tubes of 4ml, 4.5ml, 6ml or 7ml	DD-06	6.5 / 2.94
	22 inch 559 mm			UFD-LT06-2		22 3/16 x 9 x 5 1/2 inch 564 x 228 x 140 mm	230 Tubes of 4ml, 4.5ml, 6ml or 7ml		8.1 / 3.68
3ml	16 1/2 inch 419 mm	6 1/16 inch 154 mm	2	UFD-T03-2	Rack with 2 Drawers for 3ml Tubes (Φ ≤ 7/16 inch H ≤ 2 3/4 inch)	16 11/16 x 6 1/16 x 5 1/2 inch 424 x 154 x 140 mm	180 Tubes of 3ml	DD-03	6.2 / 2.83
	22 inch 559 mm			UFD-LT03-2		22 3/16 x 6 1/16 x 5 1/2 inch 564 x 154 x 140 mm	250 Tubes of 3ml		7.9 / 3.59
	16 1/2 inch 419 mm	9 1/8 inch 231 mm	3	UFD-T03-3	Rack with 3 Drawers for 3ml Tubes (Φ ≤ 7/16 inch H ≤ 2 3/4 inch)	16 11/16 x 9 1/8 x 5 1/2 inch 424 x 231 x 140 mm	270 Tubes of 3ml	DD-06	8.9 / 4.05
	22 inch 559 mm			UFD-LT03-3		22 3/16 x 9 1/8 x 5 1/2 inch 564 x 231 x 140 mm	375 Tubes of 3ml		11.4 / 5.16
1ml	16 1/2 inch 419 mm	6 1/16 inch 154 mm	2	UFD-T01-2	Rack with 2 Drawers for 1ml Tubes (Φ ≤ 5/16 inch H ≤ 2 3/4 inch)	16 11/16 x 6 1/16 x 5 1/2 inch 424 x 154 x 140 mm	240 Tubes of 1ml	DD-03	6.3 / 2.85
	22 inch 559 mm			UFD-LT01-2		22 3/16 x 6 1/16 x 5 1/2 inch 564 x 154 x 140 mm	324 Tubes of 1ml		7.9 / 3.58
	16 1/2 inch 419 mm	9 1/8 inch 231 mm	3	UFD-T01-3	Rack with 3 Drawers for 1ml Tubes (Φ ≤ 5/16 inch H ≤ 2 3/4 inch)	16 11/16 x 9 1/8 x 5 1/2 inch 424 x 231 x 140 mm	360 Tubes of 1ml	DD-06	9.0 / 4.08
	22 inch 559 mm			UFD-LT01-3		22 3/16 x 9 1/8 x 5 1/2 inch 564 x 231 x 140 mm	486 Tubes of 1ml		11.3 / 5.13

* Lock device is ordered separately.