



Upright Freezer Rack



For Matrix

500µl (0.5ml) Latch Rack (H=1 1/16 inch / 27.48mm)

For Micronic


Loborack-96 with Low Cover (H= 15/16 inch / 23.8mm)

For Fluidx

Low Profile 96-Format Rack for all 0.5ml and 0.75ml Tube with Low Profile Lid (H=1 inch / 26.0mm)

For Abgene


0.3ml 2D Barcoded Storage Rack (H=1 1/8 inch / 28.45mm)



- ◆ Made of corrosion resistant stainless steel.
- ◆ Locking rods are included.

Part Number	Dimension L x H x W	Compartment Size l x h x w	Capacity	Shipping Weight (lb/Kg)
UFMM-43	14 11/16 x 3 7/8 x 5 7/16 inch 373 x 98 x 138 mm	3 7/16 x 1 3/16 x 5 1/8 inch 88 x 30 x 130 mm	4 x 3 = 12	3.7 / 1.66
UFMM-44	14 11/16 x 5 1/16 x 5 7/16 inch 373 x 129 x 138 mm	3 7/16 x 1 3/16 x 5 1/8 inch 88 x 30 x 130 mm	4 x 4 = 16	4.7 / 2.15
UFMM-45	14 11/16 x 6 5/16 x 5 7/16 inch 373 x 160 x 138 mm	3 7/16 x 1 3/16 x 5 1/8 inch 88 x 30 x 130 mm	4 x 5 = 20	6.0 / 2.70
UFMM-46	14 11/16 x 7 1/2 x 5 7/16 inch 373 x 191 x 138 mm	3 7/16 x 1 3/16 x 5 1/8 inch 88 x 30 x 130 mm	4 x 6 = 24	7.3 / 3.29
UFMM-53	18 1/4 x 3 7/8 x 5 7/16 inch 463 x 98 x 138 mm	3 7/16 x 1 3/16 x 5 1/8 inch 88 x 30 x 130 mm	5 x 3 = 15	4.3 / 1.97
UFMM-54	18 1/4 x 5 1/16 x 5 7/16 inch 463 x 129 x 138 mm	3 7/16 x 1 3/16 x 5 1/8 inch 88 x 30 x 130 mm	5 x 4 = 20	5.6 / 2.52
UFMM-55	18 1/4 x 6 5/16 x 5 7/16 inch 463 x 160 x 138 mm	3 7/16 x 1 3/16 x 5 1/8 inch 88 x 30 x 130 mm	5 x 5 = 25	6.9 / 3.14
UFMM-56	18 1/4 x 7 1/2 x 5 7/16 inch 463 x 191 x 138 mm	3 7/16 x 1 3/16 x 5 1/8 inch 88 x 30 x 130 mm	5 x 6 = 30	8.4 / 3.79
UFMM-63	21 3/4 x 3 7/8 x 5 7/16 inch 552 x 98 x 138 mm	3 7/16 x 1 3/16 x 5 1/8 inch 88 x 30 x 130 mm	6 x 3 = 18	5.0 / 2.27
UFMM-64	21 3/4 x 5 1/16 x 5 7/16 inch 552 x 129 x 138 mm	3 7/16 x 1 3/16 x 5 1/8 inch 88 x 30 x 130 mm	6 x 4 = 24	6.4 / 2.89
UFMM-65	21 3/4 x 6 5/16 x 5 7/16 inch 552 x 160 x 138 mm	3 7/16 x 1 3/16 x 5 1/8 inch 88 x 30 x 130 mm	6 x 5 = 30	7.9 / 3.58
UFMM-66	21 3/4 x 7 1/2 x 5 7/16 inch 552 x 191 x 138 mm	3 7/16 x 1 3/16 x 5 1/8 inch 88 x 30 x 130 mm	6 x 6 = 36	9.5 / 4.30

Upright Freezer Rack



H: 1 1/4" - 1 1/2" (31.7-38.7mm)

For Micronic

Loborack-96 with Standard Cover (H=1 1/4 inch / 31.7mm)

Loborack-96 with High Cover (H=1 1/2 inch / 38.7mm)

For Fluidx

Low Profile 96-Format Rack for all 0.5ml and 0.75ml Tube with Standard Profile Lid (H=1 5/16 inch / 32.9mm)


96-Format Rack for Tube with Protective Jacket with Low Cover (H=1 1/2 inch / 38mm)

For Abgene

0.5ml 2D Barcoded Lock Storage Tube Rack (H=1 1/2 inch / 38.5mm)


0.6ml 2D Barcoded Storage Tube Rack (H=1 1/2 inch / 38.5mm)

- ◆ Made of corrosion resistant stainless steel.
- ◆ Locking rods are included.



Part Number	Dimension L x H x W	Compartment Size l x h x w	Capacity	Shipping Weight(lb/Kg)
UFML-43L	14 ¹¹ / ₁₆ x 5 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 373 x 128 x 138 mm	3 ⁷ / ₁₆ x 1 ⁹ / ₁₆ x 5 ¹ / ₈ inch 88 x 40 x 130 mm	4 x 3 = 12	3.8 / 1.74
UFML-44L	14 ¹¹ / ₁₆ x 6 ⁵ / ₈ x 5 ¹ / ₁₆ inch 373 x 168 x 138 mm	3 ⁷ / ₁₆ x 1 ⁹ / ₁₆ x 5 ¹ / ₈ inch 88 x 40 x 130 mm	4 x 4 = 16	5.0 / 2.28
UFML-45L	14 ¹¹ / ₁₆ x 8 ¹ / ₄ x 5 ⁷ / ₁₆ inch 373 x 209 x 138 mm	3 ⁷ / ₁₆ x 1 ⁹ / ₁₆ x 5 ¹ / ₈ inch 88 x 40 x 130 mm	4 x 5 = 20	6.3 / 2.86
UFML-46L	14 ¹¹ / ₁₆ x 9 ¹³ / ₁₆ x 5 ⁷ / ₁₆ inch 373 x 249 x 138 mm	3 ⁷ / ₁₆ x 1 ⁹ / ₁₆ x 5 ¹ / ₈ inch 88 x 40 x 130 mm	4 x 6 = 24	7.8 / 3.52
UFML-53L	18 ¹ / ₄ x 5 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 463 x 128 x 138 mm	3 ⁷ / ₁₆ x 1 ⁹ / ₁₆ x 5 ¹ / ₈ inch 88 x 40 x 130 mm	5 x 3 = 15	4.5 / 2.05
UFML-54L	18 ¹ / ₄ x 6 ⁵ / ₈ x 5 ⁷ / ₁₆ inch 463 x 168 x 138 mm	3 ⁷ / ₁₆ x 1 ⁹ / ₁₆ x 5 ¹ / ₈ inch 88 x 40 x 130 mm	5 x 4 = 20	5.9 / 2.66
UFML-55L	18 ¹ / ₄ x 8 ¹ / ₄ x 5 ⁷ / ₁₆ inch 463 x 209 x 138 mm	3 ⁷ / ₁₆ x 1 ⁹ / ₁₆ x 5 ¹ / ₈ inch 88 x 40 x 130 mm	5 x 5 = 25	7.3 / 3.31
UFML-56L	18 ¹ / ₄ x 9 ¹³ / ₁₆ x 5 ⁷ / ₁₆ inch 463 x 249 x 138 mm	3 ⁷ / ₁₆ x 1 ⁹ / ₁₆ x 5 ¹ / ₈ inch 88 x 40 x 130 mm	5 x 6 = 30	8.9 / 4.04
UFML-63L	21 ³ / ₄ x 5 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 552 x 128 x 138 mm	3 ⁷ / ₁₆ x 1 ⁹ / ₁₆ x 5 ¹ / ₈ inch 88 x 40 x 130 mm	6 x 3 = 18	5.2 / 2.35
UFML-64L	21 ³ / ₄ x 6 ⁵ / ₈ x 5 ⁷ / ₁₆ inch 552 x 168 x 138 mm	3 ⁷ / ₁₆ x 1 ⁹ / ₁₆ x 5 ¹ / ₈ inch 88 x 40 x 130 mm	6 x 4 = 24	6.6 / 2.99
UFML-65L	21 ³ / ₄ x 8 ¹ / ₄ x 5 ⁷ / ₁₆ inch 552 x 209 x 138 mm	3 ⁷ / ₁₆ x 1 ⁹ / ₁₆ x 5 ¹ / ₈ inch 88 x 40 x 130 mm	6 x 5 = 30	8.3 / 3.76
UFML-66L	21 ³ / ₄ x 9 ¹³ / ₁₆ x 5 ⁷ / ₁₆ inch 552 x 249 x 138 mm	3 ⁷ / ₁₆ x 1 ⁹ / ₁₆ x 5 ¹ / ₈ inch 88 x 40 x 130 mm	6 x 6 = 36	10.1 / 4.56

Upright Freezer Rack



For Matrix

- 500µl Screw-Top Latch Rack (H= 1³/₄inch / 44.32mm)
- 0.75ml Latch Rack (H= 1¹/₂inch / 39.37mm)




For Fluidx

- 384-Format Rack for Square Tube (H= 9/16 inch / 14.4mm, each compartment holds 3 racks.)
- Low Profile 96-Format Rack for all 0.5ml and 0.75ml Tube with Tall Profile Lid (H=1³/₄inch / 43.87mm)
- 96-Format Rack for 1.0ml Tube with Low Profile Lid (H= 1³/₄inch / 43.9mm)
- 96-Format Rack for Tube with Protective Jacket with Standard Profile Lid (H=1³/₄inch / 44.9mm)

For Micronic

- Roborack-96 with Low Cover (H=1⁵/₈inch / 41.2mm)




Rack w/ Locking Rod

- ◆ Made of corrosion resistant stainless steel.
- ◆ Locking rods are included.

Part Number	Dimension L x H x W	Compartment Size l x h x w	Capacity	Shipping Weight(lb/Kg)
UFDP-34L	11 ¹ / ₂ x 7 ¹ / ₂ x 5 ¹ / ₂ inch 292.5 x 191.1 x 140 mm	3 ⁹ / ₁₆ x 1 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 45.7 x 130.5 mm	3 x 4 = 12	4.2 / 1.91
UFDP-35L	11 ¹ / ₂ x 9 ³ / ₈ x 5 ¹ / ₂ inch 292.5 x 237.5 x 140 mm	3 ⁹ / ₁₆ x 1 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 45.7 x 130.5 mm	3 x 5 = 15	5.0 / 2.28
UFDP-36L	11 ¹ / ₂ x 11 ³ / ₁₆ x 5 ¹ / ₂ inch 292.5 x 283.5 x 140 mm	3 ⁹ / ₁₆ x 1 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 45.7 x 130.5 mm	3 x 6 = 18	5.8 / 2.64
UFDP-44L	15 ³ / ₁₆ x 7 ¹ / ₂ x 5 ¹ / ₂ inch 385.3 x 191.1 x 140 mm	3 ⁹ / ₁₆ x 1 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 45.7 x 130.5 mm	4 x 4 = 16	5.4 / 2.47
UFDP-45L	15 ³ / ₁₆ x 9 ³ / ₈ x 5 ¹ / ₂ inch 385.3 x 237.5 x 140 mm	3 ⁹ / ₁₆ x 1 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 45.7 x 130.5 mm	4 x 5 = 20	6.5 / 2.95
UFDP-46L	15 ³ / ₁₆ x 11 ³ / ₁₆ x 5 ¹ / ₂ inch 385.3 x 283.5 x 140 mm	3 ⁹ / ₁₆ x 1 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 45.7 x 130.5 mm	4 x 6 = 24	7.5 / 3.42
UFDP-54L	18 ¹³ / ₁₆ x 7 ¹ / ₂ x 5 ¹ / ₂ inch 478.1 x 191.1 x 140 mm	3 ⁹ / ₁₆ x 1 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 45.7 x 130.5 mm	5 x 4 = 20	6.7 / 3.03
UFDP-55L	18 ¹³ / ₁₆ x 9 ³ / ₈ x 5 ¹ / ₂ inch 478.1 x 237 x 140 mm	3 ⁹ / ₁₆ x 1 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 45.7 x 130.5 mm	5 x 5 = 25	8.2 / 3.72
UFDP-56L	18 ¹³ / ₁₆ x 11 ³ / ₁₆ x 5 ¹ / ₂ inch 478.1 x 284 x 140 mm	3 ⁹ / ₁₆ x 1 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 45.7 x 130.5 mm	5 x 6 = 30	9.5 / 4.32
UFDP-64L	22 ¹ / ₂ x 7 ¹ / ₂ x 5 ¹ / ₂ inch 571 x 191.1 x 140 mm	3 ⁹ / ₁₆ x 1 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 45.7 x 130.5 mm	6 x 4 = 24	8.0 / 3.62
UFDP-65L	22 ¹ / ₂ x 9 ⁵ / ₁₆ x 5 ¹ / ₂ inch 571 x 237 x 140 mm	3 ⁹ / ₁₆ x 1 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 45.7 x 130.5 mm	6 x 5 = 30	9.8 / 4.43
UFDP-66L	22 ¹ / ₂ x 11 ¹ / ₁₆ x 5 ¹ / ₂ inch 571 x 284 x 140 mm	3 ⁹ / ₁₆ x 1 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 45.7 x 130.5 mm	6 x 6 = 36	11.3 / 5.13
UFDP-74L	26 ¹ / ₈ x 7 ¹ / ₂ x 5 ¹ / ₂ inch 663.7 x 191.1 x 140 mm	3 ⁹ / ₁₆ x 1 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 45.7 x 130.5 mm	7 x 4 = 28	9.2 / 4.17
UFDP-75L	26 ¹ / ₈ x 9 ³ / ₈ x 5 ¹ / ₂ inch 663.7 x 237.5 x 140 mm	3 ⁹ / ₁₆ x 1 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 45.7 x 130.5 mm	7 x 5 = 35	11.3 / 5.11
UFDP-76L	26 ¹ / ₈ x 11 ³ / ₁₆ x 5 ¹ / ₂ inch 663.7 x 283.5 x 140 mm	3 ⁹ / ₁₆ x 1 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 45.7 x 130.5 mm	7 x 6 = 42	13.1 / 5.93

Upright Freezer Rack



For Matrix

1.4ml Latch Rack (H=1¹⁵/₁₆ inch / 49.84mm)

H: 5¹/₈" - 2¹/₈" (17.3-53.5mm)

For Micronic

Roborack-96 with Standard Cover (H=1¹⁵/₁₆ inch / 49.1mm)

Comorack-96 (H=2 inch / 52mm)

Comorack-24 (H=2 inch / 58mm)

For Fluidx

384-Format Rack for Round Tube (H=5⁷/₈ inch / 17.1mm, each compartment holds 3 racks.)

96-Format Rack for 1.0ml and 1.4ml Tube with Standard Profile Lid (H=2 inch / 50.8mm)

48-Format Rack for Cryo Tube with Standard Profile Lid (H=2¹/₁₆ inch / 52.8mm)


48-Format Rack for Glass Tube (H=2 inch / 51.6mm)

For Abgene

2.0ml 2D Barcoded Screw Cap Storage Tube Rack (H=2¹/₁₆ inch / 52.5mm)


1.0ml 2D Barcoded Twist Lock Storage Tube Rack (H=2¹/₈ inch / 53.5mm)

- ◆ Made of corrosion resistant stainless steel.
- ◆ Locking rods are included.



Part Number	Dimension L x H x W	Compartment Size l x h x w	Capacity	Lock Device Option*	Shipping Weight(Ib/Kg)
UFTB-33L	11 ⁹ / ₁₆ x 6 ¹³ / ₁₆ x 5 ¹ / ₂ inch 293 x 173 x 140 mm	3 ⁹ / ₁₆ x 2 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	3 x 3 = 9	DP-1107	3.7 / 1.70
UFTB-34L	11 ⁹ / ₁₆ x 9 x 5 ¹ / ₂ inch 293 x 228 x 140 mm	3 ⁹ / ₁₆ x 2 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	3 x 4 = 12	DP-1109	4.7 / 2.12
UFTB-35L	11 ⁹ / ₁₆ x 11 ³ / ₁₆ x 5 ¹ / ₂ inch 293 x 284 x 140 mm	3 ⁹ / ₁₆ x 2 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	3 x 5 = 15	DP-1111	5.6 / 2.53
UFTB-36L	11 ⁹ / ₁₆ x 13 ³ / ₈ x 5 ¹ / ₂ inch 293 x 339 x 140 mm	3 ⁹ / ₁₆ x 2 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	3 x 6 = 18	DP-1113	6.5 / 2.95
UFTB-43L	15 ³ / ₁₆ x 6 ¹³ / ₁₆ x 5 ¹ / ₂ inch 385 x 173 x 140 mm	3 ⁹ / ₁₆ x 2 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	4 x 3 = 12	DP-1407	4.9 / 2.20
UFTB-44L	15 ³ / ₁₆ x 9 x 5 ¹ / ₂ inch 385 x 228 x 140 mm	3 ⁹ / ₁₆ x 2 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	4 x 4 = 16	DP-1409	6.0 / 2.74
UFTB-45L	15 ³ / ₁₆ x 11 ³ / ₁₆ x 5 ¹ / ₂ inch 385 x 284 x 140 mm	3 ⁹ / ₁₆ x 2 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	4 x 5 = 20	DP-1411	7.2 / 3.27
UFTB-46L	15 ³ / ₁₆ x 13 ³ / ₈ x 5 ¹ / ₂ inch 385 x 339 x 140 mm	3 ⁹ / ₁₆ x 2 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	4 x 6 = 24	DP-1413	8.4 / 3.81
UFTB-53L	18 ¹³ / ₁₆ x 6 ¹³ / ₁₆ x 5 ¹ / ₂ inch 478 x 173 x 140 mm	3 ⁹ / ₁₆ x 2 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	5 x 3 = 15	DP-1807	5.9 / 2.69
UFTB-54L	18 ¹³ / ₁₆ x 9 x 5 ¹ / ₂ inch 478 x 228 x 140 mm	3 ⁹ / ₁₆ x 2 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	5 x 4 = 20	DP-1809	7.4 / 3.36
UFTB-55L	18 ¹³ / ₁₆ x 11 ³ / ₁₆ x 5 ¹ / ₂ inch 478 x 284 x 140 mm	3 ⁹ / ₁₆ x 2 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	5 x 5 = 25	DP-1811	8.9 / 4.02
UFTB-56L	18 ¹³ / ₁₆ x 13 ³ / ₈ x 5 ¹ / ₂ inch 478 x 339 x 140 mm	3 ⁹ / ₁₆ x 2 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	5 x 6 = 30	DP-1813	10.3 / 4.68
UFTB-63L	22 ¹ / ₂ x 6 ¹³ / ₁₆ x 5 ¹ / ₂ inch 571 x 173 x 140 mm	3 ⁹ / ₁₆ x 2 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	6 x 3 = 18	DP-2207	7.0 / 3.17
UFTB-64L	22 ¹ / ₂ x 9 x 5 ¹ / ₂ inch 571 x 228 x 140 mm	3 ⁹ / ₁₆ x 2 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	6 x 4 = 24	DP-2209	8.8 / 3.97
UFTB-65L	22 ¹ / ₂ x 11 ³ / ₁₆ x 5 ¹ / ₂ inch 571 x 284 x 140 mm	3 ⁹ / ₁₆ x 2 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	6 x 5 = 30	DP-2211	10.5 / 4.76
UFTB-66L	22 ¹ / ₂ x 13 ³ / ₈ x 5 ¹ / ₂ inch 571 x 339 x 140 mm	3 ⁹ / ₁₆ x 2 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	6 x 6 = 36	DP-2213	12.2 / 5.54
UFTB-73L	26 ¹ / ₈ x 6 ¹³ / ₁₆ x 5 ¹ / ₂ inch 664 x 173 x 140 mm	3 ⁹ / ₁₆ x 2 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	7 x 3 = 21	DP-2607	8.1 / 3.68
UFTB-74L	26 ¹ / ₈ x 9 x 5 ¹ / ₂ inch 664 x 228 x 140 mm	3 ⁹ / ₁₆ x 2 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	7 x 4 = 28	DP-2609	10.1 / 4.59
UFTB-75L	26 ¹ / ₈ x 11 ³ / ₁₆ x 5 ¹ / ₂ inch 664 x 284 x 140 mm	3 ⁹ / ₁₆ x 2 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	7 x 5 = 35	DP-2611	11.9 / 5.42
UFTB-76L	26 ¹ / ₈ x 13 ³ / ₈ x 5 ¹ / ₂ inch 664 x 339 x 140 mm	3 ⁹ / ₁₆ x 2 ¹³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	7 x 6 = 42	DP-2613	14.2 / 6.42

Upright Freezer Rack



For Matrix

1.0ml Screw-Top Latch Rack (H=2 1/8 inch / 54.8mm)

H: 2 5/16" - 2 7/16" (54.8-61.8mm)

For Micronic

Roborack-96 with High Cover (H=2 1/4 inch / 56.1mm)

For Fluidx

96-Format Rack for Tube with Protective Jacket with Tall Profile Lid (H=2 3/16 inch / 55mm)


96-Format Rack Tack for 1.0ml and 1.4ml Tube with Tall Profile Lid (H=2 7/16 inch / 61.8mm)

24-Format Rack (H=2 1/4 inch / 56mm)

For Abgene

1.2ml 2D Barcoded Tube in Stackable Rack (H=2 1/4 inch / 56.13mm)


1.2ml 2D Barcoded Storage Tube Rack (H=2 1/4 inch / 56.13mm)



- ◆ Made of corrosion resistant stainless steel.
- ◆ Locking rods are included.

Part Number	Dimension L x H x W	Compartment Size l x h x w	Capacity	Shipping Weight(lb/Kg)
UFMX-332	11 1/8 x 7 11/16 x 5 7/16 inch 283 x 195 x 138 mm	3 7/16 x 2 1/2 x 5 1/8 inch 88 x 62.5 x 130 mm	3 x 3 = 9	3.9 / 1.75
UFMX-342	11 1/8 x 10 3/16 x 5 7/16 inch 283 x 258 x 138 mm	3 7/16 x 2 1/2 x 5 1/8 inch 88 x 62.5 x 130 mm	3 x 4 = 12	4.8 / 2.19
UFMX-352	11 1/8 x 12 5/8 x 5 7/16 inch 283 x 321 x 138 mm	3 7/16 x 2 1/2 x 5 1/8 inch 88 x 62.5 x 130 mm	3 x 5 = 15	5.8 / 2.62
UFMX-432	14 11/16 x 7 11/16 x 5 7/16 inch 373 x 195 x 138 mm	3 7/16 x 2 1/2 x 5 1/8 inch 88 x 62.5 x 130 mm	4 x 3 = 12	5.0 / 2.25
UFMX-442	14 11/16 x 10 3/16 x 5 7/16 inch 373 x 258 x 138 mm	3 7/16 x 2 1/2 x 5 1/8 inch 88 x 62.5 x 130 mm	4 x 4 = 16	6.2 / 2.82
UFMX-452	14 11/16 x 12 5/8 x 5 7/16 inch 373 x 321 x 138 mm	3 7/16 x 2 1/2 x 5 1/8 inch 88 x 62.5 x 130 mm	4 x 5 = 20	7.5 / 3.38
UFMX-532	18 1/4 x 7 11/16 x 5 7/16 inch 463 x 195 x 138 mm	3 7/16 x 2 1/2 x 5 1/8 inch 88 x 62.5 x 130 mm	5 x 3 = 15	6.1 / 2.75
UFMX-542	18 1/4 x 10 3/16 x 5 7/16 inch 463 x 258 x 138 mm	3 7/16 x 2 1/2 x 5 1/8 inch 88 x 62.5 x 130 mm	5 x 4 = 20	7.6 / 3.45
UFMX-552	18 1/4 x 12 5/8 x 5 7/16 inch 463 x 321 x 138 mm	3 7/16 x 2 1/2 x 5 1/8 inch 88 x 62.5 x 130 mm	5 x 5 = 25	9.1 / 4.14
UFMX-562	18 1/4 x 15 1/8 x 5 7/16 inch 463 x 384 x 138 mm	3 7/16 x 2 1/2 x 5 1/8 inch 88 x 62.5 x 130 mm	5 x 6 = 30	10.5 / 4.75
UFMX-632	21 3/4 x 7 11/16 x 5 7/16 inch 552 x 195 x 138 mm	3 7/16 x 2 1/2 x 5 1/8 inch 88 x 62.5 x 130 mm	6 x 3 = 18	7.2 / 3.26
UFMX-642	21 3/4 x 10 3/16 x 5 7/16 inch 552 x 258 x 138 mm	3 7/16 x 2 1/2 x 5 1/8 inch 88 x 62.5 x 130 mm	6 x 4 = 24	9.0 / 4.08
UFMX-652	21 3/4 x 12 5/8 x 5 7/16 inch 552 x 321 x 138 mm	3 7/16 x 2 1/2 x 5 1/8 inch 88 x 62.5 x 130 mm	6 x 5 = 30	10.8 / 4.89
UFMX-662	21 3/4 x 15 1/8 x 5 7/16 inch 552 x 384 x 138 mm	3 7/16 x 2 1/2 x 5 1/8 inch 88 x 62.5 x 130 mm	6 x 6 = 36	12.4 / 5.61

Upright Freezer Rack





For Micronic

Roborack-96 with Extra High Cover (H=3 1/2 inch / 88.5mm)

For Fluidx

48-Format Rack for Cryo Tube with Tall Profile Lid (H=3 1/2 inch / 88.5mm)



- ◆ Made of corrosion resistant stainless steel.
- ◆ Locking rods are included.

Part Number	Dimension L x H x W	Compartment Size l x h x w	Capacity	Shipping Weight(lb/Kg)
UFMH-42L	14 ¹¹ / ₁₆ x 7 ⁷ / ₁₆ x 5 ³ / ₈ inch 373 x 189 x 138 mm	3 ⁷ / ₁₆ x 3 ⁷ / ₁₆ x 5 ¹ / ₈ inch 88 x 91 x 130 mm	8	4.0 / 1.82
UFMH-43L	14 ¹¹ / ₁₆ x 11 ¹ / ₁₆ x 5 ³ / ₈ inch 373 x 281 x 138 mm	3 ⁷ / ₁₆ x 3 ⁹ / ₁₆ x 5 ¹ / ₈ inch 88 x 91 x 130 mm	12	5.4 / 2.47
UFMH-52L	18 ³ / ₁₆ x 7 ⁷ / ₁₆ x 5 ³ / ₈ inch 462 x 189 x 138 mm	3 ⁷ / ₁₆ x 3 ⁹ / ₁₆ x 5 ¹ / ₈ inch 88 x 91 x 130 mm	10	4.9 / 2.23
UFMH-53L	18 ³ / ₁₆ x 11 ¹ / ₁₆ x 5 ³ / ₈ inch 462 x 281 x 138 mm	3 ⁷ / ₁₆ x 3 ⁹ / ₁₆ x 5 ¹ / ₈ inch 88 x 91 x 130 mm	15	6.5 / 2.96



For Matrix

12ml Screw-Top Latch Rack Boxes (H=4 1/4 inch / 107.82mm)

For Abgene


5.0ml 2D Barcoded Screw Cap Storage Tube Rack (H=4 inch / 101.5mm)

- ◆ Made of corrosion resistant stainless steel.
- ◆ Locking rods are included.



Part Number	Dimension L x H x W	Compartment Size l x h x w	Capacity	Shipping Weight(lb/Kg)
UFMR-42	14 ¹¹ / ₁₆ x 9 x 5 ⁷ / ₁₆ inch 373 x 229 x 138 mm	3 ⁷ / ₁₆ x 4 ³ / ₈ x 5 ¹ / ₈ inch 88 x 110.5 x 130 mm	4 x 2 = 8	4.8 / 2.17
UFMR-43	14 ¹¹ / ₁₆ x 13 ³ / ₈ x 5 ⁷ / ₁₆ inch 373 x 340 x 138 mm	3 ⁷ / ₁₆ x 4 ³ / ₈ x 5 ¹ / ₈ inch 88 x 110.5 x 130 mm	4 x 3 = 12	6.6 / 2.98
UFMR-52	18 ¹ / ₄ x 9 x 5 ⁷ / ₁₆ inch 463 x 229 x 138 mm	3 ⁷ / ₁₆ x 4 ³ / ₈ x 5 ¹ / ₈ inch 88 x 110.5 x 130 mm	5 x 2 = 10	5.8 / 2.64
UFMR-53	18 ¹ / ₄ x 13 ³ / ₈ x 5 ⁷ / ₁₆ inch 463 x 340 x 138 mm	3 ⁷ / ₁₆ x 4 ³ / ₈ x 5 ¹ / ₈ inch 88 x 110.5 x 130 mm	5 x 3 = 15	8.0 / 3.63
UFMR-62	21 ³ / ₄ x 9 x 5 ⁷ / ₁₆ inch 552 x 229 x 138 mm	3 ⁷ / ₁₆ x 4 ³ / ₈ x 5 ¹ / ₈ inch 88 x 110.5 x 130 mm	6 x 2 = 12	6.9 / 3.13
UFMR-63	21 ³ / ₄ x 13 ³ / ₈ x 5 ⁷ / ₁₆ inch 552 x 340 x 138 mm	3 ⁷ / ₁₆ x 4 ³ / ₈ x 5 ¹ / ₈ inch 88 x 110.5 x 130 mm	6 x 3 = 18	9.5 / 4.29
UFMR-72	25 ¹ / ₄ x 9 x 5 ⁷ / ₁₆ inch 642 x 229 x 138 mm	3 ⁷ / ₁₆ x 4 ³ / ₈ x 5 ¹ / ₈ inch 88 x 110.5 x 130 mm	7 x 2 = 14	8.0 / 3.62
UFMR-73	25 ¹ / ₄ x 13 ³ / ₈ x 5 ⁷ / ₁₆ inch 642 x 340 x 138 mm	3 ⁷ / ₁₆ x 4 ³ / ₈ x 5 ¹ / ₈ inch 88 x 110.5 x 130 mm	7 x 3 = 21	10.7 / 4.85

Upright Drawer Freezer Rack



For Matrix

500µl (0.5ml) Latch Rack (H=1 1/16 inch / 27.48mm)

For Micronic


Loborack-96 with Low Cover (H=1 5/16 inch / 23.8mm)

For Fluidx


Low Profile 96-Format Rack for all 0.5ml and 0.75ml Tube with Low Profile Lid (H=1 inch / 26.0mm)

For Abgene

0.3ml 2D Barcoded Storage Rack (H=1 1/8 inch / 28.45mm)




- ◆ Made of corrosion resistant stainless steel.
- ◆ Drawer racks allow easy access to and storage of samples without having to remove the whole rack from the freezer.



Part Number	Dimension L x H x W	Compartment Size l x h x w	Capacity	Shipping Weight(lb/Kg)
UFDMM-32	11 ³ / ₁₆ x 5 ³ / ₈ x 5 ⁷ / ₁₆ inch 283.4 x 136 x 138 mm	3 ¹ / ₂ x 2 ³ / ₈ x 5 ¹ / ₈ inch 89 x 61 x 130 mm	3 x 2 = 6	4.0 / 1.83
UFDMM-33	11 ³ / ₁₆ x 8 x 5 ⁷ / ₁₆ inch 283.4 x 203.4 x 138 mm	3 ¹ / ₂ x 2 ³ / ₈ x 5 ¹ / ₈ inch 89 x 61 x 130 mm	3 x 3 = 9	5.7 / 2.60
UFDMM-34	11 ³ / ₁₆ x 10 ¹¹ / ₁₆ x 5 ⁷ / ₁₆ inch 283.4 x 271 x 138 mm	3 ¹ / ₂ x 2 ³ / ₈ x 5 ¹ / ₈ inch 89 x 61 x 130 mm	3 x 4 = 12	7.5 / 3.38
UFDMM-42	14 ⁹ / ₁₆ x 5 ³ / ₈ x 5 ⁷ / ₁₆ inch 370 x 136 x 138 mm	3 ¹ / ₂ x 2 ³ / ₈ x 5 ¹ / ₈ inch 89 x 61 x 130 mm	4 x 2 = 8	5.0 / 2.27
UFDMM-43	14 ⁹ / ₁₆ x 8 x 5 ⁷ / ₁₆ inch 370 x 203.4 x 138 mm	3 ¹ / ₂ x 2 ³ / ₈ x 5 ¹ / ₈ inch 89 x 61 x 130 mm	4 x 3 = 12	7.1 / 3.24
UFDMM-44	14 ⁹ / ₁₆ x 10 ¹¹ / ₁₆ x 5 ⁷ / ₁₆ inch 370 x 271 x 138 mm	3 ¹ / ₂ x 2 ³ / ₈ x 5 ¹ / ₈ inch 89 x 61 x 130 mm	4 x 4 = 16	9.3 / 4.22
UFDMM-52	18 ¹ / ₈ x 5 ³ / ₈ x 5 ⁷ / ₁₆ inch 459.6 x 136 x 138 mm	3 ¹ / ₂ x 2 ³ / ₈ x 5 ¹ / ₈ inch 89 x 61 x 130 mm	5 x 2 = 10	6.0 / 2.74
UFDMM-53	18 ¹ / ₈ x 8 x 5 ⁷ / ₁₆ inch 459.6 x 203.4 x 138 mm	3 ¹ / ₂ x 2 ³ / ₈ x 5 ¹ / ₈ inch 89 x 61 x 130 mm	5 x 3 = 15	8.6 / 3.91
UFDMM-54	18 ¹ / ₈ x 10 ¹¹ / ₁₆ x 5 ⁷ / ₁₆ inch 459.6 x 271 x 138 mm	3 ¹ / ₂ x 2 ³ / ₈ x 5 ¹ / ₈ inch 89 x 61 x 130 mm	5 x 4 = 20	11.2 / 5.08
UFDMM-62	21 ⁵ / ₈ x 5 ³ / ₈ x 5 ⁷ / ₁₆ inch 549.2 x 136 x 138 mm	3 ¹ / ₂ x 2 ³ / ₈ x 5 ¹ / ₈ inch 89 x 61 x 130 mm	6 x 2 = 12	7.0 / 3.19
UFDMM-63	21 ⁵ / ₈ x 8 x 5 ⁷ / ₁₆ inch 549.2 x 203.4 x 138 mm	3 ¹ / ₂ x 2 ³ / ₈ x 5 ¹ / ₈ inch 89 x 61 x 130 mm	6 x 3 = 18	10.0 / 4.54
UFDMM-64	21 ⁵ / ₈ x 10 ¹¹ / ₁₆ x 5 ⁷ / ₁₆ inch 549.2 x 271 x 138 mm	3 ¹ / ₂ x 2 ³ / ₈ x 5 ¹ / ₈ inch 89 x 61 x 130 mm	6 x 4 = 24	13.0 / 5.90

Upright Freezer Drawer Rack



For Micronic

Loborack-96 with Standard Cover (H=1 1/4 inch / 31.7mm)

Loborack-96 with High Cover (H=1 1/2 inch / 38.7mm)

For Fluidx

Low Profile 96-Format Rack for all 0.5ml and 0.75ml Tube with Standard Profile Lid (H=1 5/16 inch / 32.9mm)


96-Format Rack for Tube with Protective Jacket with Low Cover (H=1 1/2 inch / 38mm)

For Abgene

0.5ml 2D Barcoded Lock Storage Tube Rack (H=1 1/2 inch / 38.5mm)


0.6ml 2D Barcoded Storage Tube Rack (H=1 1/2 inch / 38.5mm)

- ◆ Made of corrosion resistant stainless steel.
- ◆ Drawer racks allow easy access to and storage of samples without having to remove the whole rack from the freezer.



Part Number	Dimension L x H x W	Compartment Size l x h x w	Capacity	Shipping Weight (lb/Kg)
UFDML-43	14 ⁷ / ₁₆ x 5 ⁹ / ₁₆ x 5 ⁷ / ₁₆ inch 366 x 142 x 138 mm	3 ⁷ / ₁₆ x 1 ⁵ / ₈ x 5 ¹ / ₈ inch 88 x 41.4 x 130 mm	4 x 3 = 12	6.0 / 2.73
UFDML-44	14 ⁷ / ₁₆ x 7 ³ / ₈ x 5 ⁷ / ₁₆ inch 366 x 188 x 138 mm	3 ⁷ / ₁₆ x 1 ⁵ / ₈ x 5 ¹ / ₈ inch 88 x 41.4 x 130 mm	4 x 4 = 16	7.8 / 3.53
UFDML-45	14 ⁷ / ₁₆ x 9 ¹ / ₄ x 5 ⁷ / ₁₆ inch 366 x 235 x 138 mm	3 ⁷ / ₁₆ x 1 ⁵ / ₈ x 5 ¹ / ₈ inch 88 x 41.4 x 130 mm	4 x 5 = 20	9.5 / 4.33
UFDML-46	14 ⁷ / ₁₆ x 11 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 366 x 281 x 138 mm	3 ⁷ / ₁₆ x 1 ⁵ / ₈ x 5 ¹ / ₈ inch 88 x 41.4 x 130 mm	4 x 6 = 24	11.3 / 5.13
UFDML-53	17 ⁷ / ₈ x 5 ⁹ / ₁₆ x 5 ⁷ / ₁₆ inch 454 x 142 x 138 mm	3 ⁷ / ₁₆ x 1 ⁵ / ₈ x 5 ¹ / ₈ inch 88 x 41.4 x 130 mm	5 x 3 = 15	7.3 / 3.31
UFDML-54	17 ⁷ / ₈ x 7 ³ / ₈ x 5 ⁷ / ₁₆ inch 454 x 188 x 138 mm	3 ⁷ / ₁₆ x 1 ⁵ / ₈ x 5 ¹ / ₈ inch 88 x 41.4 x 130 mm	5 x 4 = 20	9.4 / 4.27
UFDML-55	17 ⁷ / ₈ x 9 ¹ / ₄ x 5 ⁷ / ₁₆ inch 454 x 235 x 138 mm	3 ⁷ / ₁₆ x 1 ⁵ / ₈ x 5 ¹ / ₈ inch 88 x 41.4 x 130 mm	5 x 5 = 25	11.6 / 5.24
UFDML-56	17 ⁷ / ₈ x 11 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 454 x 281 x 138 mm	3 ⁷ / ₁₆ x 1 ⁵ / ₈ x 5 ¹ / ₈ inch 88 x 41.4 x 130 mm	5 x 6 = 30	13.7 / 6.21
UFDML-63	21 ³ / ₈ x 5 ⁹ / ₁₆ x 5 ⁷ / ₁₆ inch 543 x 142 x 138 mm	3 ⁷ / ₁₆ x 1 ⁵ / ₈ x 5 ¹ / ₈ inch 88 x 41.4 x 130 mm	6 x 3 = 18	8.6 / 3.88
UFDML-64	21 ³ / ₈ x 7 ³ / ₈ x 5 ⁷ / ₁₆ inch 543 x 188 x 138 mm	3 ⁷ / ₁₆ x 1 ⁵ / ₈ x 5 ¹ / ₈ inch 88 x 41.4 x 130 mm	6 x 4 = 24	11.1 / 5.02
UFDML-65	21 ³ / ₈ x 9 ¹ / ₄ x 5 ⁷ / ₁₆ inch 543 x 235 x 138 mm	3 ⁷ / ₁₆ x 1 ⁵ / ₈ x 5 ¹ / ₈ inch 88 x 41.4 x 130 mm	6 x 5 = 30	15.1 / 6.86
UFDML-66	21 ³ / ₈ x 11 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 543 x 281 x 138 mm	3 ⁷ / ₁₆ x 1 ⁵ / ₈ x 5 ¹ / ₈ inch 88 x 41.4 x 130 mm	6 x 6 = 36	16.1 / 7.29

Upright Freezer Drawer Rack



For Matrix

500µl Screw-Top Latch Rack (H=1 3/4 inch / 44.32mm)

0.75ml Latch Rack (H=1 1/2 inch / 39.37mm)

H: 9 1/16" - 1 3/4" (14.4-44.9mm)

For Micronic

Roborack-96 with Low Cover (H=1 5/8 inch / 41.2mm)


For Fluidx

384-Format Rack for Square Tube (H=9/16 inch / 14.4mm, each compartment holds 3 racks.)

Low Profile 96-Format Rack for all 0.5ml and 0.75ml Tube with Tall Profile Lid (H=1 3/4 inch / 43.87mm)

96-Format Rack for 1.0ml Tube with Low Profile Lid (H=1 3/4 inch / 43.9mm)


96-Format Rack for Tube with Protective Jacket with Standard Profile Lid (H=1 3/4 inch / 44.9mm)



- ◆ Made of corrosion resistant stainless steel.
- ◆ Drawer racks allow easy access to and storage of samples without having to remove the whole rack from the freezer.

Part Number	Dimension L x H x W	Compartment Size l x h x w	Capacity	Shipping Weight(lb/Kg)
UFDDP-33	11 1/8 x 6 1/16 x 5 7/16 inch 282 x 154 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	3 x 3 = 9	5.2 / 2.38
UFDDP-34	11 1/8 x 8 1/16 x 5 7/16 inch 282 x 205 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	3 x 4 = 12	6.8 / 3.09
UFDDP-35	11 1/8 x 10 1/16 x 5 7/16 inch 282 x 256 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	3 x 5 = 15	8.4 / 3.79
UFDDP-36	11 1/8 x 12 1/8 x 5 7/16 inch 282 x 307.2 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	3 x 6 = 18	9.9 / 4.50
UFDDP-37	11 1/8 x 14 1/8 x 5 7/16 inch 282 x 358 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	3 x 7 = 21	11.7 / 5.30
UFDDP-43	14 3/8 x 6 1/16 x 5 7/16 inch 372 x 154 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	4 x 3 = 12	6.6 / 3.00
UFDDP-44	14 5/8 x 8 1/16 x 5 7/16 inch 372 x 205 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	4 x 4 = 16	8.6 / 3.89
UFDDP-45	14 3/8 x 10 1/16 x 5 7/16 inch 372 x 256 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	4 x 5 = 20	10.5 / 4.77
UFDDP-46	14 3/8 x 12 1/8 x 5 7/16 inch 372 x 307.2 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	4 x 6 = 24	12.5 / 5.66
UFDDP-47	14 5/8 x 14 1/8 x 5 7/16 inch 372 x 358 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	4 x 7 = 28	14.6 / 6.60
UFDDP-53	18 3/16 x 6 1/16 x 5 7/16 inch 462 x 154 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	5 x 3 = 15	8.0 / 3.61
UFDDP-54	18 3/16 x 8 1/16 x 5 7/16 inch 462 x 205 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	5 x 4 = 20	10.3 / 4.68
UFDDP-55	18 3/16 x 10 1/16 x 5 7/16 inch 462 x 256 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	5 x 5 = 25	12.7 / 5.75
UFDDP-56	18 3/16 x 12 1/8 x 5 7/16 inch 462 x 307.2 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	5 x 6 = 30	15.0 / 6.81
UFDDP-57	18 3/16 x 14 1/8 x 5 7/16 inch 462 x 358 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	5 x 7 = 35	17.4 / 7.90
UFDDP-63	21 3/4 x 6 1/16 x 5 7/16 inch 552 x 154 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	6 x 3 = 18	9.3 / 4.23
UFDDP-64	21 1/4 x 8 1/16 x 5 7/16 inch 552 x 205 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	6 x 4 = 24	12.1 / 5.48
UFDDP-65	21 1/4 x 10 1/16 x 5 7/16 inch 552 x 256 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	6 x 5 = 30	14.8 / 6.73
UFDDP-66	21 3/4 x 12 1/8 x 5 7/16 inch 552 x 307.2 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	6 x 6 = 36	17.6 / 7.97
UFDDP-67	21 1/4 x 14 1/8 x 5 7/16 inch 552 x 358 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	6 x 7 = 42	20.3 / 9.20
UFDDP-73	25 3/16 x 6 1/16 x 5 7/16 inch 640 x 154 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	7 x 3 = 21	10.7 / 4.84
UFDDP-74	25 3/16 x 8 1/16 x 5 7/16 inch 640 x 205 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	7 x 4 = 28	13.8 / 6.27
UFDDP-75	25 1/4 x 10 1/16 x 5 7/16 inch 642 x 256 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	7 x 5 = 35	17.0 / 7.70
UFDDP-76	25 1/4 x 12 1/8 x 5 7/16 inch 642 x 307 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	7 x 6 = 42	20.1 / 9.12
UFDDP-77	25 1/4 x 14 1/8 x 5 7/16 inch 642 x 358 x 138 mm	3 1/2 x 1 3/4 x 5 1/8 inch 89.4 x 45.4 x 130 mm	7 x 7 = 49	23.4 / 10.63

Upright Freezer Drawer Rack



For Matrix

1.4ml Latch Rack (H=1¹⁵/₁₆ inch / 49.84mm)

For Micronic

Roborack-96 with Standard Cover (H=1¹⁵/₁₆ inch / 49.1mm)

Comorack-96 (H=2 inch / 52mm)

Comorack-24 (H=2 inch / 52mm)

For Fluidx

384-Format Rack for Round Tube (H=1¹¹/₁₆ inch / 17.1mm, each compartment holds 3 racks.)

96-Format Rack for 1.0ml and 1.4ml Tube with Standard Profile Lid (H=2 inch / 50.8mm)

48-Format Rack for Cryo Tube with Standard Profile Lid (H=2¹/₁₆ inch / 52.8mm)

48-Format Rack for Glass Tube (H=2 inch / 51.6mm)

For Abgene


2.0ml 2D Barcoded Screw Cap Storage Tube Rack

(H=2 inch / 52.5mm)

1.0ml 2D Barcoded Twist Lock Storage Tube Rack


(H=2¹/₈ inch / 53.5mm)

- ◆ Made of corrosion resistant stainless steel.
- ◆ Drawer racks allow easy access to and storage of samples without having to remove the whole rack from the freezer.



Part Number	Dimension L x H x W	Compartment Size l x h x w	Capacity	Shipping Weight (lb/Kg)
UFDMP-33	11 ¹ / ₈ x 7 ⁷ / ₁₆ x 5 ⁷ / ₁₆ inch 282 x 186 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	3 x 3 = 9	5.4 / 2.45
UFDMP-34	11 ¹ / ₈ x 9 ³ / ₄ x 5 ⁷ / ₁₆ inch 282 x 247 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	3 x 4 = 12	6.9 / 3.13
UFDMP-35	11 ¹ / ₈ x 12 ¹ / ₈ x 5 ⁷ / ₁₆ inch 282 x 308 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	3 x 5 = 15	8.6 / 3.90
UFDMP-36	11 ¹ / ₈ x 14 ¹ / ₂ x 5 ⁷ / ₁₆ inch 282 x 369 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	3 x 6 = 18	10.2 / 4.64
UFDMP-43	14 ³ / ₈ x 7 ⁷ / ₁₆ x 5 ⁷ / ₁₆ inch 372 x 186 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	4 x 3 = 12	6.8 / 3.10
UFDMP-44	14 ³ / ₈ x 9 ³ / ₄ x 5 ⁷ / ₁₆ inch 372 x 247 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	4 x 4 = 16	8.9 / 4.02
UFDMP-45	14 ³ / ₈ x 12 ¹ / ₈ x 5 ⁷ / ₁₆ inch 372 x 308 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	4 x 5 = 20	10.9 / 4.94
UFDMP-46	14 ³ / ₈ x 14 ¹ / ₂ x 5 ⁷ / ₁₆ inch 372 x 369 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	4 x 6 = 24	12.9 / 5.86
UFDMP-53	18 ³ / ₁₆ x 7 ⁵ / ₁₆ x 5 ⁷ / ₁₆ inch 462 x 186 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	5 x 3 = 15	8.2 / 3.70
UFDMP-54	18 ³ / ₁₆ x 9 ³ / ₄ x 5 ⁷ / ₁₆ inch 462 x 247 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	5 x 4 = 20	10.6 / 4.79
UFDMP-55	18 ³ / ₁₆ x 12 ¹ / ₈ x 5 ⁷ / ₁₆ inch 462 x 308 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	5 x 5 = 25	13.0 / 5.90
UFDMP-56	18 ³ / ₁₆ x 14 ¹ / ₂ x 5 ⁷ / ₁₆ inch 462 x 369 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	5 x 6 = 30	15.4 / 7.00
UFDMP-63	21 ³ / ₄ x 7 ⁷ / ₁₆ x 5 ⁷ / ₁₆ inch 552 x 186 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	6 x 3 = 18	9.7 / 4.42
UFDMP-64	21 ³ / ₄ x 9 ³ / ₄ x 5 ⁷ / ₁₆ inch 552 x 247 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	6 x 4 = 24	12.6 / 5.73
UFDMP-65	21 ³ / ₄ x 12 ¹ / ₈ x 5 ⁷ / ₁₆ inch 552 x 308 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	6 x 5 = 30	15.5 / 7.05
UFDMP-66	21 ³ / ₄ x 14 ¹ / ₂ x 5 ⁷ / ₁₆ inch 552 x 369 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	6 x 6 = 36	17.9 / 8.14
UFDMP-73	25 ¹ / ₄ x 7 ⁷ / ₁₆ x 5 ⁷ / ₁₆ inch 642 x 186 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	7 x 3 = 21	11.0 / 4.97
UFDMP-74	25 ¹ / ₄ x 9 ³ / ₄ x 5 ⁷ / ₁₆ inch 642 x 247 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	7 x 4 = 28	14.2 / 6.43
UFDMP-75	25 ¹ / ₄ x 12 ¹ / ₈ x 5 ⁷ / ₁₆ inch 642 x 308 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	7 x 5 = 35	17.4 / 7.91
UFDMP-76	25 ¹ / ₄ x 14 ¹ / ₂ x 5 ⁷ / ₁₆ inch 642 x 369 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	7 x 6 = 42	20.7 / 9.38

Upright Freezer Drawer Rack



For Matrix

1.0ml Screw-Top Latch Rack (H=2¹/₈inch / 54.8mm)



For Micronic

Roborack-96 with High Cover (H=2¹/₄inch / 56.1mm)

For Fluidx

96-Format Rack for Tube with Protective Jacket with Tall Profile Lid (H=2³/₁₆inch / 55mm)

24-Format Rack (H=2¹/₄inch / 56mm)


96-Format Rack for 1.0ml and 1.4ml Tube with Tall Profile Lid (H=2⁷/₁₆inch / 61.8mm)













For Abgene

1.2ml 2D Barcoded Tube in Stackable Rack (H=2¹/₄inch / 56.13mm)

1.2ml 2D Barcoded Storage Tube Rack (H=2¹/₄inch / 56.13mm)

- ◆ Made of corrosion resistant stainless steel.
- ◆ Drawer racks allow easy access to and storage of samples without having to remove the whole rack from the freezer.



Part Number	Dimension L x H x W	Compartment Size l x h x w	Capacity	Shipping Weight(Ib/Kg)
UFDMX-332	10 ⁷ / ₈ x 7 ¹⁵ / ₁₆ x 5 ⁷ / ₁₆ inch 276.5 x 202 x 138 mm	3 ⁷ / ₁₆ x 2 ⁷ / ₁₆ x 5 ¹ / ₁₆ inch 88 x 62.4 x 129 mm	3 x 3 = 9 	5.4 / 2.45
UFDMX-342	10 ⁷ / ₈ x 10 ⁹ / ₁₆ x 5 ⁷ / ₁₆ inch 276.5 x 268.5 x 138 mm	3 ⁷ / ₁₆ x 2 ⁷ / ₁₆ x 5 ¹ / ₁₆ inch 88 x 62.4 x 129 mm	3 x 4 = 12 	7.0 / 3.17
UFDMX-352	10 ⁷ / ₈ x 13 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 276.5 x 335 x 138 mm	3 ⁷ / ₁₆ x 2 ⁷ / ₁₆ x 5 ¹ / ₁₆ inch 88 x 62.4 x 129 mm	3 x 5 = 15 	8.6 / 3.90
UFDMX-432	14 ³ / ₈ x 7 ¹⁵ / ₁₆ x 5 ⁷ / ₁₆ inch 365 x 202 x 138 mm	3 ⁷ / ₁₆ x 2 ⁷ / ₁₆ x 5 ¹ / ₁₆ inch 88 x 62.4 x 129 mm	4 x 3 = 12 	6.7 / 3.05
UFDMX-442	14 ³ / ₈ x 10 ⁹ / ₁₆ x 5 ⁷ / ₁₆ inch 365 x 268.5 x 138 mm	3 ⁷ / ₁₆ x 2 ⁷ / ₁₆ x 5 ¹ / ₁₆ inch 88 x 62.4 x 129 mm	4 x 4 = 16 	8.7 / 3.95
UFDMX-452	14 ³ / ₈ x 13 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 365 x 335 x 138 mm	3 ⁷ / ₁₆ x 2 ⁷ / ₁₆ x 5 ¹ / ₁₆ inch 88 x 62.4 x 129 mm	4 x 5 = 20 	10.7 / 4.86
UFDMX-532	17 ⁷ / ₈ x 7 ¹⁵ / ₁₆ x 5 ⁷ / ₁₆ inch 454 x 202 x 138 mm	3 ⁷ / ₁₆ x 2 ⁷ / ₁₆ x 5 ¹ / ₁₆ inch 88 x 62.4 x 129 mm	5 x 3 = 15 	8.1 / 3.68
UFDMX-542	17 ⁷ / ₈ x 10 ⁹ / ₁₆ x 5 ⁷ / ₁₆ inch 454 x 268.5 x 138 mm	3 ⁷ / ₁₆ x 2 ⁷ / ₁₆ x 5 ¹ / ₁₆ inch 88 x 62.4 x 129 mm	5 x 4 = 20 	10.5 / 4.78
UFDMX-552	17 ⁷ / ₈ x 13 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 454 x 335 x 138 mm	3 ⁷ / ₁₆ x 2 ⁷ / ₁₆ x 5 ¹ / ₁₆ inch 88 x 62.4 x 129 mm	5 x 5 = 25 	12.7 / 5.78
UFDMX-632	21 ³ / ₈ x 7 ¹⁵ / ₁₆ x 5 ⁷ / ₁₆ inch 542.5x202x138 mm	3 ⁷ / ₁₆ x 2 ⁷ / ₁₆ x 5 ¹ / ₁₆ inch 88 x 62.4 x 129 mm	6 x 3 = 18 	9.6 / 4.34
UFDMX-642	21 ³ / ₈ x 10 ⁹ / ₁₆ x 5 ⁷ / ₁₆ inch 542.5x268.5x138 mm	3 ⁷ / ₁₆ x 2 ⁷ / ₁₆ x 5 ¹ / ₁₆ inch 88 x 62.4 x 129 mm	6 x 4 = 24 	12.4 / 5.62
UFDMX-652	21 ³ / ₈ x 13 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 542.5x335x138 mm	3 ⁷ / ₁₆ x 2 ⁷ / ₁₆ x 5 ¹ / ₁₆ inch 88 x 62.4 x 129 mm	6 x 5 = 30 	15.2 / 6.91