








Upright Freezer Drawer Racks

For Storage Bottles

- ◆ Excellent for organizing -20°C and -80°C Freezers & Cold Rooms.
- ◆ Made of corrosion resistant stainless steel.
- ◆ Drawer Racks facilitate both access to and storage of samples without having to remove the whole racks from the freezer.
- ◆ Security lock devices are available for each rack.



Content to Store	Rack Length	Rack Height	Number of Drawers	Part Number	Description	Dimension L x H x W	Capacity	Lock Device Option*	Shipping Weight (Ib/Kg)	
Bottle $H \leq 4\frac{13}{16}$ inch 122 mm 	16 $\frac{1}{2}$ inch 424mm	4 $\frac{15}{16}$ inch 126 mm	1	UFD-B-1	Rack with 1 Drawer for Bottles	16 $\frac{11}{16}$ x 4 $\frac{15}{16}$ x 5 $\frac{1}{2}$ inch 424 x 126 x 140 mm	Drawer Area: 15 x 4 $\frac{1}{2}$ inch (Per Drawer) 383 x 114 mm	DD-04	3.4 / 1.54	
	22 inch 564mm			22 $\frac{7}{16}$ x 4 $\frac{15}{16}$ x 5 $\frac{1}{2}$ inch 564 x 126 x 140 mm		Drawer Area: 20 $\frac{5}{8}$ x 4 $\frac{1}{2}$ inch (Per Drawer) 523 x 114 mm	4.3 / 1.95			
	26 $\frac{3}{4}$ inch 684mm			26 $\frac{15}{16}$ x 4 $\frac{15}{16}$ x 5 $\frac{1}{2}$ inch 684 x 126 x 140 mm		Drawer Area: 25 $\frac{1}{2}$ x 4 $\frac{1}{2}$ inch (Per Drawer) 644 x 114 mm	5.0 / 2.27			
		16 $\frac{1}{2}$ inch 424mm	9 $\frac{15}{16}$ inch 252 mm	2	UFD-B-2	Rack with 2 Drawer for Bottles	16 $\frac{11}{16}$ x 9 $\frac{15}{16}$ x 5 $\frac{1}{2}$ inch 424 x 252 x 140 mm	Drawer Area: 15 x 4 $\frac{1}{2}$ inch (Per Drawer) 383 x 114 mm	DD-09	6.1 / 2.77
		22 inch 564mm			22 $\frac{7}{16}$ x 9 $\frac{15}{16}$ x 5 $\frac{1}{2}$ inch 564 x 252 x 140 mm		Drawer Area: 20 $\frac{5}{8}$ x 4 $\frac{1}{2}$ inch (Per Drawer) 523 x 114 mm	7.5 / 3.42		
		26 $\frac{3}{4}$ inch 684mm			26 $\frac{15}{16}$ x 9 $\frac{15}{16}$ x 5 $\frac{1}{2}$ inch 684 x 252 x 140 mm		Drawer Area: 25 $\frac{1}{2}$ x 4 $\frac{1}{2}$ inch (Per Drawer) 644 x 114 mm	8.7 / 3.94		
Bottle $H \leq 4\frac{7}{16}$ inch 106 mm 	16 $\frac{1}{2}$ inch 424mm	4 $\frac{1}{2}$ inch 114 mm	1	UFD-B10-1	Rack with 1 Drawer for Bottles	16 $\frac{11}{16}$ x 4 $\frac{1}{2}$ x 5 $\frac{1}{2}$ inch 424 x 114 x 140 mm	Drawer Area: 15 x 4 $\frac{1}{2}$ inch (Per Drawer) 383 x 114 mm	DD-04	3.0 / 1.34	
	22 inch 564mm			22 $\frac{7}{16}$ x 4 $\frac{1}{2}$ x 5 $\frac{1}{2}$ inch 564 x 114 x 140 mm		Drawer Area: 20 $\frac{5}{8}$ x 4 $\frac{1}{2}$ inch (Per Drawer) 523 x 114 mm	3.7 / 1.68			
		16 $\frac{1}{2}$ inch 424mm	9 inch 228 mm	2	UFD-B10-2	Rack with 2 Drawer for Bottles	16 $\frac{11}{16}$ x 9 x 5 $\frac{1}{2}$ inch 424 x 228 x 140 mm	Drawer Area: 15 x 4 $\frac{1}{2}$ inch (Per Drawer) 383 x 114 mm	DD-06	5.6 / 2.53
		22 inch 564mm			22 $\frac{7}{16}$ x 9 x 5 $\frac{1}{2}$ inch 564 x 228 x 140 mm		Drawer Area: 20 $\frac{5}{8}$ x 4 $\frac{1}{2}$ inch (Per Drawer) 523 x 114 mm	6.8 / 3.07		
Bottle $H \leq 2\frac{7}{4}$ inch 70 mm 	16 $\frac{1}{2}$ inch 424mm	6 $\frac{1}{16}$ inch 154 mm	2	UFD-B03-2	Rack with 2 Drawer for Bottles	16 $\frac{11}{16}$ x 6 $\frac{1}{16}$ x 5 $\frac{1}{2}$ inch 424 x 154 x 140 mm	Drawer Area: 15 x 4 $\frac{1}{2}$ inch (Per Drawer) 383 x 114 mm	DD-03	4.9 / 2.20	
	22 inch 564mm			22 $\frac{7}{16}$ x 6 $\frac{1}{16}$ x 5 $\frac{1}{2}$ inch 564 x 154 x 140 mm		Drawer Area: 20 $\frac{5}{8}$ x 4 $\frac{1}{2}$ inch (Per Drawer) 523 x 114 mm	6.1 / 2.75			
		16 $\frac{1}{2}$ inch 424mm	9 $\frac{1}{8}$ inch 231 mm	3	UFD-B03-3	Rack with 3 Drawer for Bottles	16 $\frac{11}{16}$ x 9 $\frac{1}{8}$ x 5 $\frac{1}{2}$ inch 424 x 231 x 140 mm	Drawer Area: 15 x 4 $\frac{1}{2}$ inch (Per Drawer) 383 x 114 mm	DD-06	6.9 / 3.12
		22 inch 564mm			22 $\frac{7}{16}$ x 9 $\frac{1}{8}$ x 5 $\frac{1}{2}$ inch 564 x 231 x 140 mm		Drawer Area: 20 $\frac{5}{8}$ x 4 $\frac{1}{2}$ inch (Per Drawer) 523 x 114 mm	8.6 / 3.89		

* Lock device is ordered separately.