




Upright Freezer Racks













For 3.75" High Standard Boxes



Polycarbonate Plastic Box
(Check Page 110 for Boxes)


- ◆ Racks are for 3.75" high standard polycarbonate boxes.
- ◆ Made of corrosion resistant stainless steel.



Part Number	Description	Dimension L x H x W	Compartment Size l x h x w	Capacity (Box)	Shipping Weight lb(Kg)
UF-22-375	Rack Only	11 ³ / ₈ x 8 ¹ / ₁₆ x 5 ¹ / ₂ inch 289 x 205 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 136 x 100 x 137 mm	2 x 2 =4 	2.9(1.33)
UF-23-375	Rack Only	11 ³ / ₈ x 12 ¹ / ₁₆ x 5 ¹ / ₂ inch 289 x 307 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 136 x 100 x 137 mm	2 x 3 =6 	3.9(1.78)
UF-24-375	Rack Only	11 ³ / ₈ x 16 ¹ / ₁₆ x 5 ¹ / ₂ inch 289 x 408 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 136 x 100 x 137 mm	2 x 4 =8 	4.9(2.24)
UF-32-375	Rack Only	16 ¹³ / ₁₆ x 8 ¹ / ₁₆ x 5 ¹ / ₂ inch 427 x 205 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 136 x 100 x 137 mm	3 x 2 =6 	4.1(1.87)
UF-33-375	Rack Only	16 ¹³ / ₁₆ x 12 ¹ / ₁₆ x 5 ¹ / ₂ inch 427 x 307 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 136 x 100 x 137 mm	3 x 3 =9 	5.5(2.51)
UF-34-375	Rack Only	16 ¹³ / ₁₆ x 16 ¹ / ₁₆ x 5 ¹ / ₂ inch 427 x 408 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 136 x 100 x 137 mm	3 x 4 =12 	6.9(3.15)
UF-42-375	Rack Only	22 ¹ / ₈ x 8 ¹ / ₁₆ x 5 ¹ / ₂ inch 562 x 205 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 136 x 100 x 137 mm	4 x 2 =8 	5.3(2.42)
UF-43-375	Rack Only	22 ¹ / ₈ x 12 ¹ / ₁₆ x 5 ¹ / ₂ inch 562 x 307 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 136 x 100 x 137 mm	4 x 3 =12 	7.1(3.24)
UF-44-375	Rack Only	22 ¹ / ₈ x 16 ¹ / ₁₆ x 5 ¹ / ₂ inch 562 x 408 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 136 x 100 x 137 mm	4 x 4 =16 	9.0(4.06)
UF-52-375	Rack Only	27 ³ / ₁₆ x 8 ¹ / ₁₆ x 5 ¹ / ₂ inch 691 x 205 x 140 mm	5 ¹ / ₄ x 3 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 133.5 x 100 x 137 mm	5 x 2 =10 	6.5(2.94)
UF-53-375	Rack Only	27 ³ / ₁₆ x 12 ¹ / ₁₆ x 5 ¹ / ₂ inch 691 x 307 x 140 mm	5 ¹ / ₄ x 3 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 133.5 x 100 x 137 mm	5 x 3 =15 	8.7(3.94)
UF-54-375	Rack Only	27 ³ / ₁₆ x 16 ¹ / ₁₆ x 5 ¹ / ₂ inch 691 x 408 x 140 mm	5 ¹ / ₄ x 3 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 133.5 x 100 x 137 mm	5 x 4 =20 	10.9(4.94)


Upright Freezer Drawer Racks













For 3.75" High Standard Boxes



Polycarbonate Plastic Box
(Check Page 110 for Boxes)

- ◆ Racks are for 3.75" high standard polycarbonate boxes.
- ◆ These drawer racks facilitate both access to and storage of samples without having to remove the whole rack from freezer.
- ◆ Made of corrosion resistant stainless steel.




Part Number	Description	Dimension L x H x W	Compartment Size l x h x w	Capacity (Box)	Shipping Weight lb(Kg)
UFD-22-375	Rack Only	11 ¹ / ₄ x 8 ³ / ₈ x 5 ¹ / ₂ inch 286 x 213 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹ / ₄ inch 137 x 100 x 133 mm	2 x 2 =4	 4.3(1.97)
UFD-23-375	Rack Only	11 ¹ / ₄ x 12 ¹ / ₂ x 5 ¹ / ₂ inch 286 x 318 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹ / ₄ inch 137 x 100 x 133 mm	2 x 3 =6	 7.1(2.81)
UFD-24-375	Rack Only	11 ¹ / ₄ x 16 ¹¹ / ₁₆ x 5 ¹ / ₂ inch 286 x 424 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹ / ₄ inch 137 x 100 x 133 mm	2 x 4 =8	 8.1(3.66)
UFD-32-375	Rack Only	16 ⁵ / ₈ x 8 ³ / ₈ x 5 ¹ / ₂ inch 423 x 213 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹ / ₄ inch 137 x 100 x 133 mm	3 x 2 =6	 6.0(2.65)
UFD-33-375	Rack Only	16 ⁵ / ₈ x 12 ¹ / ₂ x 5 ¹ / ₂ inch 423 x 318 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹ / ₄ inch 137 x 100 x 133 mm	3 x 3 =9	 8.3(3.78)
UFD-34-375	Rack Only	16 ⁵ / ₈ x 16 ¹¹ / ₁₆ x 5 ¹ / ₂ inch 423 x 424 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹ / ₄ inch 137 x 100 x 133 mm	3 x 4 =12	 10.8(4.92)
UFD-42-375	Rack Only	22 ³ / ₁₆ x 8 ³ / ₈ x 5 ¹ / ₂ inch 563 x 213 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹ / ₄ inch 137 x 100 x 133 mm	4 x 2 =8	 7.4(3.34)
UFD-43-375	Rack Only	22 ³ / ₁₆ x 12 ¹ / ₂ x 5 ¹ / ₂ inch 563 x 318 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹ / ₄ inch 137 x 100 x 133 mm	4 x 3 =12	 10.5(4.76)
UFD-44-375	Rack Only	22 ³ / ₁₆ x 16 ¹¹ / ₁₆ x 5 ¹ / ₂ inch 563 x 424 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹ / ₄ inch 137 x 100 x 133 mm	4 x 4 =16	 13.6(6.19)
UFD-52-375	Rack Only	26 ⁷ / ₈ x 8 ³ / ₈ x 5 ¹ / ₂ inch 683 x 213 x 140 mm	5 ⁵ / ₁₆ x 3 ¹⁵ / ₁₆ x 5 ¹ / ₄ inch 134.5 x 100 x 133 mm	5 x 2 =10	 8.8(3.97)
UFD-53-375	Rack Only	26 ⁷ / ₈ x 12 ¹ / ₂ x 5 ¹ / ₂ inch 683 x 318 x 140 mm	5 ⁵ / ₁₆ x 3 ¹⁵ / ₁₆ x 5 ¹ / ₄ inch 134.5 x 100 x 133 mm	5 x 3 =15	 11.9(5.40)
UFD-54-375	Rack Only	26 ⁷ / ₈ x 16 ¹¹ / ₁₆ x 5 ¹ / ₂ inch 683 x 424 x 140 mm	5 ⁵ / ₁₆ x 3 ¹⁵ / ₁₆ x 5 ¹ / ₄ inch 134.5 x 100 x 133 mm	5 x 4 =20	 16.2(7.34)

Vertical Racks


(Locking Rod Included)

For 3.75" High Standard Boxes

- ◆ Racks are for 3.75" high standard polycarbonate boxes.
- ◆ Made of corrosion resistant stainless steel.



Polycarbonate Plastic Box
(Check Page 110 for Boxes)





Part Number	Description	Dimension H x W x D	Compartment Size h x w x d	Capacity (Box)	Shipping Weight lb(Kg)
CF-4-375	Rack Only	16 ¹ / ₄ x 5 ⁵ / ₈ x 5 ¹ / ₂ inch 413 x 143 x 140 mm	3 ¹⁵ / ₁₆ x 5 ¹ / ₄ x 5 ³ / ₈ inch 100 x 134 x 136 mm	4	3.6(1.64)
CF-5-375	Rack Only	20 ¹ / ₄ x 5 ⁵ / ₈ x 5 ¹ / ₂ inch 514 x 143 x 140 mm	3 ¹⁵ / ₁₆ x 5 ¹ / ₄ x 5 ³ / ₈ inch 100 x 134 x 136 mm	5	4.5(2.03)
CF-6-375	Rack Only	24 ¹ / ₄ x 5 ⁵ / ₈ x 5 ¹ / ₂ inch 616 x 143 x 140 mm	3 ¹⁵ / ₁₆ x 5 ¹ / ₄ x 5 ³ / ₈ inch 100 x 134 x 136 mm	6	5.1(2.31)
CF-7-375	Rack Only	28 x 5 ⁵ / ₈ x 5 ¹ / ₂ inch 717 x 143 x 140 mm	3 ¹⁵ / ₁₆ x 5 ¹ / ₄ x 5 ³ / ₈ inch 100 x 134 x 136 mm	7	6.0(2.70)

Upright Freezer Racks

For Large 3.75" High Boxes

- ◆ Made of corrosion resistant stainless steel.





Part Number	Dimension L x H x W	Compartment Size l x h x w	Capacity (Box)	Shipping Weight(lb/Kg)
UF-22-L375	12 ⁷ / ₁₆ x 8 ¹ / ₁₆ x 5 ¹⁵ / ₁₆ inch 315.5 x 205 x 151 mm	5 ⁷ / ₈ x 4 x 5 ¹³ / ₁₆ inch 149 x 101 x 148 mm	2 x 2 =4	3.1 / 1.39
UF-23-L375	12 ⁷ / ₁₆ x 12 ¹ / ₁₆ x 5 ¹⁵ / ₁₆ inch 315.5 x 307 x 151 mm	5 ⁷ / ₈ x 4 x 5 ¹³ / ₁₆ inch 149 x 101 x 148 mm	2 x 3 =6	4.1 / 1.86
UF-24-L375	12 ⁷ / ₁₆ x 16 ¹ / ₁₆ x 5 ¹⁵ / ₁₆ inch 315.5 x 408 x 151 mm	5 ⁷ / ₈ x 4 x 5 ¹³ / ₁₆ inch 149 x 101 x 148 mm	2 x 4 =8	4.7 / 2.14
UF-32-L375	18 ³ / ₈ x 8 ¹ / ₁₆ x 5 ¹⁵ / ₁₆ inch 466 x 205 x 151 mm	5 ⁷ / ₈ x 4 x 5 ¹³ / ₁₆ inch 149 x 101 x 148 mm	3 x 2 =6	4.3 / 1.96
UF-33-L375	18 ³ / ₈ x 12 ¹ / ₁₆ x 5 ¹⁵ / ₁₆ inch 466 x 307 x 151 mm	5 ⁷ / ₈ x 4 x 5 ¹³ / ₁₆ inch 149 x 101 x 148 mm	3 x 3 =9	5.8 / 2.62
UF-34-L375	18 ³ / ₈ x 16 ¹ / ₁₆ x 5 ¹⁵ / ₁₆ inch 466 x 408 x 151 mm	5 ⁷ / ₈ x 4 x 5 ¹³ / ₁₆ inch 149 x 101 x 148 mm	3 x 4 =12	7.2 / 3.28
UF-42-L375	24 ⁵ / ₁₆ x 8 ¹ / ₁₆ x 5 ¹⁵ / ₁₆ inch 617 x 205 x 151 mm	5 ⁷ / ₈ x 4 x 5 ¹³ / ₁₆ inch 149 x 101 x 148 mm	4 x 2 =8	6.0 / 2.71
UF-43-L375	24 ⁵ / ₁₆ x 12 ¹ / ₁₆ x 5 ¹⁵ / ₁₆ inch 617 x 307 x 151 mm	5 ⁷ / ₈ x 4 x 5 ¹³ / ₁₆ inch 149 x 101 x 148 mm	4 x 3 =12	8.0 / 3.62
UF-44-L375	24 ⁵ / ₁₆ x 16 ¹ / ₁₆ x 5 ¹⁵ / ₁₆ inch 617 x 408 x 151 mm	5 ⁷ / ₈ x 4 x 5 ¹³ / ₁₆ inch 149 x 101 x 148 mm	4 x 4 =16	9.9 / 4.48
UF-52-L375	30 ¹ / ₄ x 8 ¹ / ₁₆ x 5 ¹⁵ / ₁₆ inch 768 x 205 x 151 mm	5 ⁷ / ₈ x 4 x 5 ¹³ / ₁₆ inch 149 x 101 x 148 mm	5 x 2 =10	7.3 / 3.30
UF-53-L375	30 ¹ / ₄ x 12 ¹ / ₁₆ x 5 ¹⁵ / ₁₆ inch 768 x 307 x 151 mm	5 ⁷ / ₈ x 4 x 5 ¹³ / ₁₆ inch 149 x 101 x 148 mm	5 x 3 =15	9.8 / 4.45
UF-54-L375	30 ¹ / ₄ x 16 ¹ / ₁₆ x 5 ¹⁵ / ₁₆ inch 768 x 408 x 151 mm	5 ⁷ / ₈ x 4 x 5 ¹³ / ₁₆ inch 149 x 101 x 148 mm	5 x 4 =20	12.1 / 5.50

Upright Freezer Drawer Racks

For Large 3.75" High Boxes

◆ Made of corrosion resistant stainless steel.





Part Number	Description	Dimension L x H x W	Compartment Size l x h x w	Capacity (Box)	Shipping Weight(lb/Kg)
UFD-22-L375	Rack Only	12 ³ / ₁₆ x 8 ³ / ₈ x 6 ¹ / ₈ inch 310 x 213 x 155 mm	5 ⁷ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 149 x 100 x 148 mm	2 x 2 =4	4.9 / 2.22
UFD-23-L375	Rack Only	12 ³ / ₁₆ x 12 ¹ / ₂ x 6 ¹ / ₈ inch 310 x 318 x 155 mm	5 ⁷ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 149 x 100 x 148 mm	2 x 3 =6	7.0 / 3.17
UFD-24-L375	Rack Only	12 ³ / ₁₆ x 16 ¹¹ / ₁₆ x 6 ¹ / ₈ inch 310 x 424 x 155 mm	5 ⁷ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 149 x 100 x 148 mm	2 x 4 =8	9.1 / 4.11
UFD-32-L375	Rack Only	18 ¹ / ₁₆ x 8 ³ / ₈ x 6 ¹ / ₈ inch 459 x 213 x 155 mm	5 ⁷ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 149 x 100 x 148 mm	3 x 2 =6	6.6 / 2.98
UFD-33-L375	Rack Only	18 ¹ / ₁₆ x 12 ¹ / ₂ x 6 ¹ / ₈ inch 459 x 318 x 155 mm	5 ⁷ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 149 x 100 x 148 mm	3 x 3 =9	9.4 / 4.25
UFD-34-L375	Rack Only	18 ¹ / ₁₆ x 16 ¹¹ / ₁₆ x 6 ¹ / ₈ inch 459 x 424 x 155 mm	5 ⁷ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 149 x 100 x 148 mm	3 x 4 =12	12.2 / 5.54
UFD-42-L375	Rack Only	23 ⁷ / ₈ x 8 ³ / ₈ x 6 ¹ / ₈ inch 609 x 213 x 155 mm	5 ⁷ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 149 x 100 x 148 mm	4 x 2 =8	8.3 / 3.76
UFD-43-L375	Rack Only	24 x 12 ¹ / ₂ x 6 ¹ / ₈ inch 609 x 318 x 155 mm	5 ⁷ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 149 x 100 x 148 mm	4 x 3 =12	11.8 / 5.36
UFD-44-L375	Rack Only	24 x 16 ¹¹ / ₁₆ x 6 ¹ / ₈ inch 609 x 424 x 155 mm	5 ⁷ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 149 x 100 x 148 mm	4 x 4 =16	15.3 / 6.92
UFD-52-L375	Rack Only	29 ¹³ / ₁₆ x 8 ³ / ₈ x 6 ¹ / ₈ inch 758 x 213 x 155 mm	5 ⁷ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 149 x 100 x 148 mm	5 x 2 =10	10.0 / 4.53
UFD-53-L375	Rack Only	29 ¹³ / ₁₆ x 12 ¹ / ₂ x 6 ¹ / ₈ inch 758 x 318 x 155 mm	5 ⁷ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 149 x 100 x 148 mm	5 x 3 =15	14.2 / 6.44
UFD-54-L375	Rack Only	29 ¹³ / ₁₆ x 16 ¹¹ / ₁₆ x 6 ¹ / ₈ inch 758 x 424 x 155 mm	5 ⁷ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 149 x 100 x 148 mm	5 x 4 =20	18.4 / 8.35

Vertical Racks

(Locking Rod Included)

For Large 3.75" High Boxes

◆ Made of corrosion resistant stainless steel.



Part Number	Description	Dimension H x W x D	Compartment Size h x w x d	Capacity (Box)	Shipping Weight(lb/Kg)
CF-4-L375	Rack Only	16 ¹ / ₄ x 6 ¹ / ₄ x 6 inch 413 x 158 x 153 mm	4 x 5 ⁷ / ₈ x 5 ⁷ / ₈ inch 101 x 150 x 149 mm	4	2.5 / 1.12
CF-5-L375	Rack Only	20 ¹ / ₄ x 6 ¹ / ₄ x 6 inch 514 x 158 x 153 mm	4 x 5 ⁷ / ₈ x 5 ⁷ / ₈ inch 101 x 150 x 149 mm	5	3.0 / 1.34
CF-6-L375	Rack Only	24 ¹ / ₄ x 6 ¹ / ₄ x 6 inch 616 x 158 x 153 mm	4 x 5 ⁷ / ₈ x 5 ⁷ / ₈ inch 101 x 150 x 149 mm	6	3.4 / 1.55
CF-7-L375	Rack Only	28 ¹ / ₄ x 6 ¹ / ₄ x 6 inch 717 x 158 x 153 mm	4 x 5 ⁷ / ₈ x 5 ⁷ / ₈ inch 101 x 150 x 149 mm	7	3.9 / 1.76