



# Chest Freezer Racks

## For 100 Cell Hinged Lid Storage Boxes



[ Check page 110 for Boxes ]

- ◆ Made of corrosion resistant stainless steel.
- ◆ Locking rods are included.
- ◆ Security lock devices are available for each rack.



Chest Freezer Rack Only



Chest Freezer Rack w/100-Cell Hinged Plastic Storage Box



w/DR Lock Device Option

Part Number	Description	Dimension H x D x W	Compartment Size h x d x w	Capacity (Box)	Lock Device Option*	Shipping Weight lb(Kg)
CFHT-8	Rack Only	18 <sup>3</sup> / <sub>4</sub> x 5 <sup>7</sup> / <sub>8</sub> x 6 <sup>1</sup> / <sub>4</sub> inch 476 x 149 x 159 mm	2 <sup>5</sup> / <sub>16</sub> x 5 <sup>9</sup> / <sub>16</sub> x 6 <sup>1</sup> / <sub>8</sub> inch 58 x 141 x 155 mm	8	DR-17	3.6(1.62)
CFHT-9	Rack Only	21 x 5 <sup>7</sup> / <sub>8</sub> x 6 <sup>1</sup> / <sub>4</sub> inch 534 x 149 x 159 mm	2 <sup>5</sup> / <sub>16</sub> x 5 <sup>9</sup> / <sub>16</sub> x 6 <sup>1</sup> / <sub>8</sub> inch 58 x 141 x 155 mm	9	DR-19	3.9(1.77)
CFHT-10	Rack Only	23 <sup>3</sup> / <sub>8</sub> x 5 <sup>7</sup> / <sub>8</sub> x 6 <sup>1</sup> / <sub>4</sub> inch 593 x 149 x 159 mm	2 <sup>5</sup> / <sub>16</sub> x 5 <sup>9</sup> / <sub>16</sub> x 6 <sup>1</sup> / <sub>8</sub> inch 58 x 141 x 155 mm	10	DR-22	4.3(1.95)
CFHT-11	Rack Only	25 <sup>5</sup> / <sub>8</sub> x 5 <sup>7</sup> / <sub>8</sub> x 6 <sup>1</sup> / <sub>4</sub> inch 651 x 149 x 159 mm	2 <sup>5</sup> / <sub>16</sub> x 5 <sup>9</sup> / <sub>16</sub> x 6 <sup>1</sup> / <sub>8</sub> inch 58 x 141 x 155 mm	11	DR-25	4.7(2.12)
CFHT-12	Rack Only	27 <sup>15</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>8</sub> x 6 <sup>1</sup> / <sub>4</sub> inch 710 x 149 x 159 mm	2 <sup>5</sup> / <sub>16</sub> x 5 <sup>9</sup> / <sub>16</sub> x 6 <sup>1</sup> / <sub>8</sub> inch 58 x 141 x 155 mm	12	DR-27	5.0(2.28)
CFHT-13	Rack Only	30 <sup>1</sup> / <sub>4</sub> x 5 <sup>7</sup> / <sub>8</sub> x 6 <sup>1</sup> / <sub>4</sub> inch 768 x 149 x 159 mm	2 <sup>5</sup> / <sub>16</sub> x 5 <sup>9</sup> / <sub>16</sub> x 6 <sup>1</sup> / <sub>8</sub> inch 58 x 141 x 155 mm	13	DR-29	5.4(2.45)

\* Lock device is ordered separately.