



Chest Freezer Racks

For 25-Place Slide Boxes


- ◆ Made of corrosion resistant stainless steel.
- ◆ Locking rods are included.
- ◆ Security lock devices are available for each rack.



L1 x W1 x H1: 5.5" x 3.6" x 1.4" / 140 x 91 x 35mm



L1 x W1 x H1: 5.5" x 3.3" x 1.4" / 140 x 83 x 34mm



Part Number	Dimension H x D x W	Compartment Size h x d x w	Capacity	Lock Device Option*	Shipping Weight(lb/Kg)
CFS-25-4	6 ³ / ₈ x 3 ¹³ / ₁₆ x 6 inch 162 x 97 x 152 mm	1 ¹ / ₂ x 3 ¹¹ / ₁₆ x 5 ⁵ / ₈ inch 38 x 93 x 143 mm	4	DR25-06	1.4 / 0.62
CFS-25-6	9 ⁷ / ₁₆ x 3 ¹³ / ₁₆ x 6 inch 239 x 97 x 152 mm	1 ¹ / ₂ x 3 ¹¹ / ₁₆ x 5 ⁵ / ₈ inch 38 x 93 x 143 mm	6	DR25-09	1.8 / 0.82
CFS-25-8	12 ⁷ / ₁₆ x 3 ¹³ / ₁₆ x 6 inch 316 x 97 x 152 mm	1 ¹ / ₂ x 3 ¹¹ / ₁₆ x 5 ⁵ / ₈ inch 38 x 93 x 143 mm	8	DR25-12	2.3 / 1.03
CFS-25-10	15 ¹ / ₂ x 3 ¹³ / ₁₆ x 6 inch 393 x 97 x 152 mm	1 ¹ / ₂ x 3 ¹¹ / ₁₆ x 5 ⁵ / ₈ inch 38 x 93 x 143 mm	10	DR25-15	2.7 / 1.23
CFS-25-12	18 ¹ / ₂ x 3 ¹³ / ₁₆ x 6 inch 470 x 97 x 152 mm	1 ¹ / ₂ x 3 ¹¹ / ₁₆ x 5 ⁵ / ₈ inch 38 x 93 x 143 mm	12	DR25-18	3.2 / 1.45
CFS-25-14	21 ³ / ₁₀ x 3 ¹³ / ₁₆ x 6 inch 547 x 97 x 152 mm	1 ¹ / ₂ x 3 ¹¹ / ₁₆ x 5 ⁵ / ₈ inch 38 x 93 x 143 mm	14	DR25-21	3.7 / 1.66
CFS-25-16	24 ⁹ / ₁₆ x 3 ¹³ / ₁₆ x 6 inch 624 x 97 x 152 mm	1 ¹ / ₂ x 3 ¹¹ / ₁₆ x 5 ⁵ / ₈ inch 38 x 93 x 143 mm	16	DR25-24	4.1 / 1.87
CFS-25-18	27 ⁵ / ₈ x 3 ¹³ / ₁₆ x 6 inch 701 x 97 x 152 mm	1 ¹ / ₂ x 3 ¹¹ / ₁₆ x 5 ⁵ / ₈ inch 38 x 93 x 143 mm	18	DR25-27	4.6 / 2.07
CFS-25-20	30 ⁵ / ₈ x 3 ¹³ / ₁₆ x 6 inch 778 x 97 x 152 mm	1 ¹ / ₂ x 3 ¹¹ / ₁₆ x 5 ⁵ / ₈ inch 38 x 93 x 143 mm	20	DR25-30	5.0 / 2.28

*Security lock device is ordered separately.