



# Upright Freezer Racks

## For 100-Place Slide Boxes



L1 x W1 x H1: 8.25" x 7.0" x 1.34"  
210 x 178 x 34mm

- ◆ Made of corrosion resistant stainless steel.
- ◆ Security lock devices are available for each rack.



Part Number	Dimension L x H x W	Compartment Size l x h x w	Capacity	Lock Device Option*	Shipping Weight (lb/Kg)
UFS-100-24	14 <sup>7</sup> / <sub>8</sub> x 6 <sup>1</sup> / <sub>4</sub> x 8 <sup>3</sup> / <sub>8</sub> inch 377.5 x 159 x 213 mm	7 <sup>3</sup> / <sub>16</sub> x 1 <sup>1</sup> / <sub>2</sub> x 8 <sup>5</sup> / <sub>16</sub> inch 182 x 38 x 211 mm	8	DP100-1306	5.5 / 2.50
UFS-100-34	22 <sup>1</sup> / <sub>16</sub> x 6 <sup>1</sup> / <sub>4</sub> x 8 <sup>3</sup> / <sub>8</sub> inch 559.5 x 159 x 213 mm	7 <sup>3</sup> / <sub>16</sub> x 1 <sup>1</sup> / <sub>2</sub> x 8 <sup>5</sup> / <sub>16</sub> inch 182 x 38 x 211 mm	12	DP100-2106	7.9 / 3.60
UFS-100-44	29 <sup>13</sup> / <sub>16</sub> x 6 <sup>1</sup> / <sub>4</sub> x 8 <sup>3</sup> / <sub>8</sub> inch 741 x 159 x 213 mm	7 <sup>3</sup> / <sub>16</sub> x 1 <sup>1</sup> / <sub>2</sub> x 8 <sup>5</sup> / <sub>16</sub> inch 182 x 38 x 211 mm	16	DP100-2806	8.8 / 4.00
UFS-100-25	14 <sup>7</sup> / <sub>8</sub> x 7 <sup>13</sup> / <sub>16</sub> x 8 <sup>3</sup> / <sub>8</sub> inch 377.5 x 198 x 213 mm	7 <sup>3</sup> / <sub>16</sub> x 1 <sup>1</sup> / <sub>2</sub> x 8 <sup>5</sup> / <sub>16</sub> inch 182 x 38 x 211 mm	10	DP100-1308	6.7 / 3.04
UFS-100-35	22 <sup>1</sup> / <sub>16</sub> x 7 <sup>13</sup> / <sub>16</sub> x 8 <sup>3</sup> / <sub>8</sub> inch 559.5 x 198 x 213 mm	7 <sup>3</sup> / <sub>16</sub> x 1 <sup>1</sup> / <sub>2</sub> x 8 <sup>5</sup> / <sub>16</sub> inch 182 x 38 x 211 mm	15	DP100-2108	9.5 / 4.32
UFS-100-45	29 <sup>13</sup> / <sub>16</sub> x 7 <sup>13</sup> / <sub>16</sub> x 8 <sup>3</sup> / <sub>8</sub> inch 741 x 198 x 213 mm	7 <sup>3</sup> / <sub>16</sub> x 1 <sup>1</sup> / <sub>2</sub> x 8 <sup>5</sup> / <sub>16</sub> inch 182 x 38 x 211 mm	20	DP100-2808	12.6 / 5.70
UFS-100-26	14 <sup>7</sup> / <sub>8</sub> x 9 <sup>5</sup> / <sub>16</sub> x 8 <sup>3</sup> / <sub>8</sub> inch 377.5 x 237 x 213 mm	7 <sup>3</sup> / <sub>16</sub> x 1 <sup>1</sup> / <sub>2</sub> x 8 <sup>5</sup> / <sub>16</sub> inch 182 x 38 x 211 mm	12	DP100-1309	7.4 / 3.34
UFS-100-36	22 <sup>1</sup> / <sub>16</sub> x 9 <sup>5</sup> / <sub>16</sub> x 8 <sup>3</sup> / <sub>8</sub> inch 559.5 x 237 x 213mm	7 <sup>3</sup> / <sub>16</sub> x 1 <sup>1</sup> / <sub>2</sub> x 8 <sup>5</sup> / <sub>16</sub> inch 182 x 38 x 211 mm	18	DP100-2109	10.6 / 4.83
UFS-100-46	29 <sup>13</sup> / <sub>16</sub> x 9 <sup>5</sup> / <sub>16</sub> x 8 <sup>3</sup> / <sub>8</sub> inch 741 x 237 x 213 mm	7 <sup>3</sup> / <sub>16</sub> x 1 <sup>1</sup> / <sub>2</sub> x 8 <sup>5</sup> / <sub>16</sub> inch 182 x 38 x 211 mm	24	DP100-2809	13.9 / 6.30
UFS-100-27	14 <sup>7</sup> / <sub>8</sub> x 10 <sup>7</sup> / <sub>8</sub> x 8 <sup>3</sup> / <sub>8</sub> inch 377.5 x 276 x 213 mm	7 <sup>3</sup> / <sub>16</sub> x 1 <sup>1</sup> / <sub>2</sub> x 8 <sup>5</sup> / <sub>16</sub> inch 182 x 38 x 211 mm	14	DP100-1311	8.2 / 3.72
UFS-100-37	22 <sup>1</sup> / <sub>16</sub> x 10 <sup>7</sup> / <sub>8</sub> x 8 <sup>3</sup> / <sub>8</sub> inch 559.5 x 276 x 213 mm	7 <sup>3</sup> / <sub>16</sub> x 1 <sup>1</sup> / <sub>2</sub> x 8 <sup>5</sup> / <sub>16</sub> inch 182 x 38 x 211 mm	21	DP100-2111	11.8 / 5.35
UFS-100-47	29 <sup>13</sup> / <sub>16</sub> x 10 <sup>7</sup> / <sub>8</sub> x 8 <sup>3</sup> / <sub>8</sub> inch 741 x 276 x 213 mm	7 <sup>3</sup> / <sub>16</sub> x 1 <sup>1</sup> / <sub>2</sub> x 8 <sup>5</sup> / <sub>16</sub> inch 182 x 38 x 211 mm	28	DP100-2811	15.5 / 7.02
UFS-100-28	14 <sup>7</sup> / <sub>8</sub> x 12 <sup>3</sup> / <sub>8</sub> x 8 <sup>3</sup> / <sub>8</sub> inch 377.5 x 315 x 213 mm	7 <sup>3</sup> / <sub>16</sub> x 1 <sup>1</sup> / <sub>2</sub> x 8 <sup>5</sup> / <sub>16</sub> inch 182 x 38 x 211 mm	16	DP100-1312	9.2 / 4.19
UFS-100-38	22 <sup>1</sup> / <sub>16</sub> x 12 <sup>3</sup> / <sub>8</sub> x 8 <sup>3</sup> / <sub>8</sub> inch 559.5 x 315 x 213 mm	7 <sup>3</sup> / <sub>16</sub> x 1 <sup>1</sup> / <sub>2</sub> x 8 <sup>5</sup> / <sub>16</sub> inch 182 x 38 x 211 mm	24	DP100-2112	13.5 / 6.13
UFS-100-48	29 <sup>13</sup> / <sub>16</sub> x 12 <sup>3</sup> / <sub>8</sub> x 8 <sup>3</sup> / <sub>8</sub> inch 741 x 315 x 213 mm	7 <sup>3</sup> / <sub>16</sub> x 1 <sup>1</sup> / <sub>2</sub> x 8 <sup>5</sup> / <sub>16</sub> inch 182 x 38 x 211 mm	32	DP100-2812	17.5 / 7.94

\* Security lock device is ordered separately.