



Upright Freezer Drawer Racks

For 100-Place Slide Boxes



L1 x W1 x H1: 8.25" x 7.0" x 1.34"
210 x 178 x 34mm

- ◆ Made of corrosion resistant stainless steel.
- ◆ Drawer racks allow easy access to and storage of samples without having to remove the whole rack from the freezer.
- ◆ Security lock devices are available for each rack.



Part Number	Dimension L x H x W	Compartment Size l x h x w	Capacity	Lock Device Option*	Shipping Weight(lb/Kg)
UFDS-100-24	14 ⁷ / ₈ x 6 ¹⁵ / ₁₆ x 8 ⁵ / ₈ inch 378 x 176.5 x 219 mm	7 ¹ / ₈ x 1 ¹ / ₂ x 8 ¹ / ₄ inch 180.5 x 38.4 x 210 mm	8	DD100-06	9.1 / 4.12
UFDS-100-34	22 ³ / ₁₆ x 6 ¹⁵ / ₁₆ x 8 ⁵ / ₈ inch 563.5 x 176.5 x 219 mm	7 ¹ / ₈ x 1 ¹ / ₂ x 8 ¹ / ₄ inch 180.5 x 38.4 x 210 mm	12	DD100-06	12.8 / 5.80
UFDS-100-44	29 ⁵ / ₁₆ x 6 ¹⁵ / ₁₆ x 8 ⁵ / ₈ inch 745 x 176.5 x 219 mm	7 ¹ / ₈ x 1 ¹ / ₂ x 8 ¹ / ₄ inch 180.5 x 38.4 x 210 mm	16	DD100-06	16.2 / 7.35
UFDS-100-25	15 ¹ / ₁₆ x 8 ¹¹ / ₁₆ x 8 ⁵ / ₈ inch 382 x 220 x 219 mm	7 ¹ / ₈ x 1 ¹ / ₂ x 8 ¹ / ₄ inch 180.5 x 38.4 x 210 mm	10	DD100-08	11.4 / 5.18
UFDS-100-35	22 ³ / ₁₆ x 8 ¹¹ / ₁₆ x 8 ⁵ / ₈ inch 563.5 x 220 x 219 mm	7 ¹ / ₈ x 1 ¹ / ₂ x 8 ¹ / ₄ inch 180.5 x 38.2 x 210 mm	15	DD100-08	16.1 / 7.30
UFDS-100-45	29 ⁵ / ₁₆ x 8 ¹¹ / ₁₆ x 8 ⁵ / ₈ inch 745 x 220 x 219 mm	7 ¹ / ₈ x 1 ¹ / ₂ x 8 ¹ / ₄ inch 180.5 x 38.2 x 210 mm	20	DD100-08	20.5 / 9.20
UFDS-100-26	15 ¹ / ₁₆ x 10 ³ / ₈ x 8 ⁵ / ₈ inch 382 x 263.5 x 219 mm	7 ¹ / ₈ x 1 ¹ / ₂ x 8 ¹ / ₄ inch 180.5 x 38.2 x 210 mm	12	DD100-09	13.1 / 5.92
UFDS-100-36	22 ³ / ₁₆ x 10 ³ / ₈ x 8 ⁵ / ₈ inch 563.5 x 263.5 x 219 mm	7 ¹ / ₈ x 1 ¹ / ₂ x 8 ¹ / ₄ inch 180.5 x 38.2 x 210 mm	18	DD100-09	18.4 / 8.34
UFDS-100-46	29 ⁵ / ₁₆ x 10 ³ / ₈ x 8 ⁵ / ₈ inch 745 x 263.5 x 219 mm	7 ¹ / ₈ x 1 ¹ / ₂ x 8 ¹ / ₄ inch 180.5 x 38.2 x 210 mm	24	DD100-09	23.6 / 10.55
UFDS-100-27	15 ¹ / ₁₆ x 12 ¹ / ₁₆ x 8 ⁵ / ₈ inch 382 x 307 x 219 mm	7 ¹ / ₈ x 1 ¹ / ₂ x 8 ¹ / ₄ inch 180.5 x 38.2 x 210 mm	14	DD100-11	14.7 / 6.65
UFDS-100-37	22 ³ / ₁₆ x 12 ¹ / ₁₆ x 8 ⁵ / ₈ inch 563.5 x 307 x 219 mm	7 ¹ / ₈ x 1 ¹ / ₂ x 8 ¹ / ₄ inch 180.5 x 38.2 x 210 mm	21	DD100-11	20.7 / 9.40
UFDS-100-47	29 ⁵ / ₁₆ x 12 ¹ / ₁₆ x 8 ⁵ / ₈ inch 745 x 307 x 219 mm	7 ¹ / ₈ x 1 ¹ / ₂ x 8 ¹ / ₄ inch 180.5 x 38.2 x 210 mm	28	DD100-11	26.1 / 11.85
UFDS-100-28	15 ¹ / ₁₆ x 13 ¹³ / ₁₆ x 8 ⁵ / ₈ inch 382 x 351 x 219 mm	7 ¹ / ₈ x 1 ¹ / ₂ x 8 ¹ / ₄ inch 180.5 x 38.2 x 210 mm	16	DD100-13	17.1 / 7.74
UFDS-100-38	22 ³ / ₁₆ x 13 ¹³ / ₁₆ x 8 ⁵ / ₈ inch 563.5 x 351 x 219 mm	7 ¹ / ₈ x 1 ¹ / ₂ x 8 ¹ / ₄ inch 180.5 x 38.2 x 210 mm	24	DD100-13	24.0 / 10.90
UFDS-100-48	29 ⁵ / ₁₆ x 13 ¹³ / ₁₆ x 8 ⁵ / ₈ inch 745 x 351 x 219 mm	7 ¹ / ₈ x 1 ¹ / ₂ x 8 ¹ / ₄ inch 180.5 x 38.2 x 210 mm	32	DD100-13	30.4 / 13.80

* Security lock device is ordered separately.