



## Silver Series Manual Defrost -45°C Chest Freezers

The LabRepCo FUTURA Silver Series manual defrost -45 °C chest freezers are designed for laboratory & pharmaceutical applications that require more rigorous temperature control and construction features. These durable, easy to use units feature an intelligent microprocessor controller with digital display and adjustable temperature control.

### Standard Product Features:

- Digital temperature display
- Intelligent microprocessor controller with digital display and adjustable temperature control
- Visual and audible hi or low temperature, sensor failure, door open alarms
- Keyed door locks
- Magnetically sealed door with rapid cool-down mode after door opening for quick temperature recovery
- Cabinet constructed of white powder-coated 304 steel with CFC-free, high density insulation (2.4" thick)
- Easy swivel casters
- Doors: 1 Lid
- Defrost: Manual
- Built and shipped from ISO 9001:2008 manufacturing and distribution facilities
- Each cabinet undergoes three quality control checks prior to shipping
- One wire basket (L2X-7-C45 and L2X-11-C45 only)
- Certification: UL / cUL
- Warranty: Two-year parts & compressor



**L2X-3-C45**



**L2X-7-C45**



**L2X-11-C45**

Model	CAPACITY (cu. ft.)	Exterior Dimensions (HxWxD)	Temperature Range (Setpoint)	Shipping Weight
L2X-3-C45	3.5	32.5" x 27.5" x 20.5"		102 lbs
L2X-7-C45	7	35.5" x 37" x 26.5"	-10 °C to -45°C Preset -40 °C	113.5 lbs
L2X-11-C45	11	35.5" x 50" x 26.5"		176 lbs