



TOP VIEW



FRONT VIEW



SIDE VIEW

**GENERAL DESCRIPTION**

4' Paramount Ductless Enclosure, Tall & Deep  
Model: 6963403

**Exterior Superstructure**

- \* Tempered glass
- \* Dry powder epoxy-coated aluminum

**Interior Superstructure**

- \* Dry powder epoxy-coated Aluminum and Steel Sides
- \* 1/4" Thick Tempered Safety Glass

**Filters**

- \* Accommodates 6 Filters
- \* Gas-assisted filter clamping mechanism

**Electrical Fixtures**

- \* Pre-wired Fluorescent Light with Timers

**Control Panel / LCD Display**

- \* Face Velocity
- \* Filter Status and Types
- \* Touchpad Control of Enclosure Operation
- \* Set Filter Life and Check Timers
- \* Programmable Alarm Set Points
- \* Event timer

**Wiring**

- \* Wired at Factory
- \* 8/2 Power Cord with Plug

**SASH**

- \* Gas-Assisted
- \* Pivots to full open for cleaning and material transfer
- \* 1/4" Tempered Safety Glass
- \* 12" Operating Height
- \* Epoxy Coated Aluminum Frame

**Baffles**

- \* Perforated, Epoxy Coated, Pivots Downward for Cleaning
- \* Air Foil
- \* Extruded aluminum, Clean Sweep Airfoil

**Motor**

- \* ECM Impeller, Smart-Flow System to Adjust Motor Speed

**SENSOR**

- \* Safety First Vapor Sensor
- \* Primary and Secondary Back-Up Sensors

**UTILITY PORT**

- \* (2) 1.5" diameter utility ports with plugs

**Assembly and Crating**

- \* Pre-Assembled
- \* Palletized

**Exhaust**

- \* Exhaust Shroud, recirculates filtered air back into room

FUME HOODS for:



Labconco Corporation  
8811 Prospect Avenue  
Kansas City, Mo 64132-2696  
(800) 821-5525  
www.labconco.com

No.	Date
Date:	4-7-2009
Drawing #	6963403
Drawn By:	Luke Savage
Checked By:	
Approved By:	
Sheet No.	

Log Number: \_\_\_\_\_ Order Number: \_\_\_\_\_