

Purpose-Built for Critical Healthcare and Research Applications

**For a Limited Time Only, Get 20% Off All
Corepoint Scientific Blood Bank Refrigerators & Freezers!**



NSBR492WSG/0

NSBR231WSWCR/0

Blood Bank Refrigerator	Ft³	Set Point	Defrost	Doors	Drawers	Width"	Depth"	Height"	Chart Recorder
NSBR231WSG-0	23	4°C	Cycle	1 Glass	7 Adj. Sliding	26 ¾	36 ¾	81 ½	No
NSBR231WSGCR-0	23	4°C	Cycle	1 Glass	7 Adj. Sliding	26 ¾	36 ¾	81 ½	Yes
NSBR231WSW-0	23	4°C	Cycle	1 Solid	7 Adj. Sliding	26 ¾	35 ¾	81 ½	No
NSBR231WSWCR-0	23	4°C	Cycle	1 Solid	7 Adj. Sliding	26 ¾	35 ¾	81 ½	Yes
NSBR492WSG-0	49	4°C	Cycle	2 Glass	14 Adj. Sliding	54	36 ¾	81 ½	No
NSBR492WSGCR-0	49	4°C	Cycle	2 Glass	14 Adj. Sliding	54	36 ¾	81 ½	Yes
NSBR492WSW-0	49	4°C	Cycle	2 Solid	14 Adj. Sliding	54	35 ¾	81 ½	No
NSBR492WSWCR-0	49	4°C	Cycle	2 Solid	14 Adj. Sliding	54	35 ¾	81 ½	Yes

[View Refrigerators](#)



NSBF231WSW/0

NSBF492WSWCR/0

Plasma Freezer	Ft³	Set Point	Defrost	Doors	Drawers	Width"	Depth"	Height"	Chart Recorder
NSBF231WSW-0	23	-35°C	Auto	1 Solid	7 Adj. Sliding	26 ¾	35 ¾	81 ½	No
NSBF231WSWCR-0	23	-35°C	Auto	1 Solid	7 Adj. Sliding	26 ¾	35 ¾	81 ½	Yes
NSBF492WSW-0	49	-35°C	Auto	2 Solid	14 Adj. Sliding	54	35 ¾	81 ½	No
NSBF492WSWCR-0	49	-35°C	Auto	2 Solid	14 Adj. Sliding	54	35 ¾	81 ½	Yes

[View Freezers](#)

www.labrepco.com

LABRepCo
An Employee-Owned Company

Prices shown online are already discounted, no promo code needed. Discount excludes shipping/freight charges. Promo expires December 31st, 2023.