

TW Series
TAYLOR WHARTON
by IC BIOMEDICAL LN₂ FREEZERS

TW35



Specifically designed for the every day requirements of laboratory and biorepository operations

- Universal canister compatibility
- Fits in standard shipping boxes and containers
- Rugged ribbed high-strength aluminum body, magneformed necktube design and more durable paint
- Superior vacuum performance with insulation technology that delivers maximum holding times
- Lockable lid with enlarged opening for seals or securing ties

SPECIFICATIONS	TW35
LN ₂ Capacity (L)	35
Static Hold Time (Days)	136
Evap. Rate (Liters/Day)	0.25
DIMENSIONS	
Necktube (in/mm)	4.7 / 119.4
Overall Dia. (in/mm)	18.8 / 477.5
Height (in/mm)	27 / 685.8
WEIGHT	
Empty	35 / 15.9
Full (lbs/kg)	96 / 43.5
CAPACITY	
Number of Canisters	6
Canister Diameter (in/mm)	3.7 x 11 / 94 x 279.4
1.2ml & 2.0 ml vials (5/cane)	1050
1.2ml & 2.0 ml vials (6/cane)	1260
.5 cc Straws (10/cane)	3300
.5 cc Straws- Bulk (2 Levels)	9600
.25 cc Straws (28/cane)	9240
.25 cc Straws- Bulk (2 Levels) in Goblets	21564



855.750.8191 | sales@ICBiomedical.com
www.ICBiomedical.com
Made In USA



MADE IN USA

© 2022, IC Biomedical, LLC