



phasetwo®

Cryogenic Shipping Case

For safe transportation at cryogenic temperatures.

Experience peace of mind with phasetwo's innovative Shipping Cases – meticulously crafted to prioritize your product's safety throughout its journey. Our lightweight yet durable cases shield valuable cargo during transit with tailored designs incorporating inner shock-absorbing materials for a snug fit of the Cryogenic Shipping Case.

The Shipping Case offers security through a lock and optional monitoring device for temperature and tracking. Seamlessly manage logistics with dedicated space for documents and labels that's easy to clean. Intelligent features include custom molded lid ribs for secure stacking, permanent orientation markings for precise transit positioning, and built-in D-rings for anchoring during transit and storage. Trust in phasetwo's commitment to your product's protection.

FEATURES

- Lockable Lid
- Premium materials provide industry leading LN2 performance
- Custom design made for stacking





phasetwo®

GENERAL SPECIFICATIONS

Case Dimensions	
Outer (Case Shell w/ Hardware)	
Height (in/mm)	22(559)
Width (in/mm)	17.18(436)
Depth (in/mm)	18.96(482)
Inner Usable	
Diameter (in/mm)	9.25(235)
Height (in/mm)	16.125(410)

Case Materials	Lid	Base
Material	Linear Low-Density Polyethylene (LLDPE)	Linear Low-Density Polyethylene (LLDPE)
Weight (lbs/kg)	3.95(1.72)	11.35(5.148)
Lid Hinges	Southco N6-4E-466-24 Stainless Steel	N/A
Latches / Keepers	Southco 37-10-086-10 Rubber and Zinc Plated Steel	F Style Southco 37-10-086-10 Zinc Plated Steel
Handles	N/A	Penn Elcom H4054Z - Surface Mounted Flip Handle
Tie-Down Ring	N/A	McMaster-Carr 3076T34 Zinc Plated Steel 400 lbs Capacity
Foam	2.2# Polyethylene Pad	2.2# Polyethylene Insert