

A New Level Of Security & Convenience For Your Environmental Room

Options that work together to boost convenience and security for environmental rooms.

Security Latch with Electronic Flush Strike

- Control access to your environmental room and protect inside contents
- Integrates with existing security systems
- Accepts card readers that are made for mounting in a junction box.
- Electronic door strike is UL certified and allows easy exiting
- Overlapping cover design provides increased tamper resistance



Push Bar Inside Release

- The push bar, combined with the required security latch and electronic strike option, creates secure and convenient access for your environmental room
- Made from durable stainless steel and aluminum, the push bar holds up to daily usage and works with either left or right-handed doors

Foot Treadle

- Foot treadles feature a pedal mounted to the bottom of the door that, when pressed, automatically opens the door
- It's easy to carry trays, boxes, or other objects without needing to use your hands to pull the door open
- Since entry is hands-free, foot treadles also provide a more sanitary option for door opening

