



## Neon Stool 9571, Bench-to-Counter Height, Freewheeling Casters, Black Vinyl, Foot Step, With 4D Arms

**Article #:20025409-A\***

*\*9-11 Week Factory Order*

---

### Features

- The award-winning Neon Chair is known as the best lab chair globally for its design, superior quality, ergonomic features and overall comfort.
- Innovative Seat Cushion System allows for easy removal and replacements.
- All-inclusive ergonomic functionality with premium mechanism, back height lumbar adjustment, back lock, seat height adjustment, swivel-tilt, forward tilt, tilt tension and 4-D Arms.

---

### Specifications

<b>Height Adjustment:</b>	Bench-to-Counter Height Range 25-1/2" - 36-1/4"
<b>Seat Depth Range:</b>	16-5/8" - 18-7/8"
<b>Accent Trim Color:</b>	In-Stock (Blue); Factory Orders (Orange, Green and Grey)
<b>Frame Color:</b>	Grey Back (Non-ESD Version)
<b>Seat Material:</b>	Black Vinyl and Enhanced Comfort Seat
<b>Castors:</b>	Freewheeling
<b>Base:</b>	Black Aluminum Base
<b>Cleanroom Rating:</b>	ISO Level 7
<b>Bio-Safety BSL:</b>	Level 2
<b>BIFMA Rating:</b>	Complies with Level e3 Platinum
<b>Lab Grade Rating:</b>	Complies with Scientific Equipment & Furniture Assn.- SEFA12
<b>Weight Capacity:</b>	300 lbs.
<b>Product Weight:</b>	40 lbs.
<b>Warranty:</b>	10-Year Warranty
<b>Original Design:</b>	German Engineering and Production
<b>ESD Capable:</b>	Available by Request