LASER CUT ID

Freezer Racks

Make rack identification part of your storage system.

IDENTIFY AT A GLANCE

Laser-cut identifiers make it easier to distinguish racks in your freezer.

BUILT FOR COLD STORAGE

Corrosion-resistant stainless steel designed for standard 2-inch boxes.

CHOOSE YOUR FORMAT

Upright and vertical rack options for different freezer layouts.



Upright Freezer Rack



Vertical Freezer Rack

Find your rack

Models for standard 2-inch boxes

UPRIGHT FREEZER RACKS	CONFIGURATION	BOXES	D × W × H (IN.)
CUF-442-L	4 deep × 4 high	16	22.3 × 5.5 × 9.22
CUF-452-L	4 deep × 5 high	20	22.3 × 5.5 × 11
CUF-462-L	4 deep × 6 high	24	22.3 × 5.5 × 12.59
CUF-472-L	4 deep × 7 high	28	22.3 × 5.5 × 15.38
CUF-542-L	5 deep × 4 high	20	27.7 × 5.5 × 9.22
CUF-552-L	5 deep × 5 high	25	27.5 × 5.5 × 11.03

LN2/CHEST FREEZER RACKS	CONFIGURATION	BOXES	D × W × H (IN.)
CCF-13-2-L	13 high	13	5.55 × 5.5 × 28.6
CCF-14-2-L	14 high	14	5.5 × 5.69 × 29.45
CCF-15-2-L	15 high	15	5.5 × 5.69 × 32.95

Vertical racks include locking rods. Boxes and dividers not included. Confirm freezer fit.

Ready to organize your freezer?

Explore models or contact our team for help selecting a rack.

VIEW RACK OPTIONS

800-521-0754

labrepco.com