LABELING & IDENTIFICATION

Freezer Rack Labeling Guide

Find the right identification option for freezer racks, drawers, storage boxes, and blood bag canisters.



Explore the options

02 Aluminum labels & label cards

03 Engraving, laser-cut IDs & riveted digits


04 Keychain identifiers & handheld printing

Contact our team to discuss rack compatibility and custom identification.

LABELS & CARDS

Identify racks and drawers

Choose a slide-in label, riveted label, or plastic card for your application.



Custom Aluminum Labels

LABEL-X-XX

Colored anodized aluminum labels slide into compatible label holders. Choose red, blue, green, or yellow and up to two letters or numbers.

2.5 in. W × 1.5 in. H | Ask about 3- or 4-character labels

VIEW PRODUCT >



Riveted Aluminum Labels

LABEL-HR-X-XX

Colored labels are permanently riveted to the rack. Choose up to two letters or numbers and your preferred color.

Must be ordered with the freezer rack at time of purchase.

VIEW PRODUCT >



Upright Drawer Rack Label Cards

DL-P

Plastic label cards fit the card holders on upright freezer drawer racks.

2 11/16 × 1 3/8 in. | Two packs of 50 cards

VIEW PRODUCT >

Confirm your rack model and holder compatibility before ordering.

DIRECT IDENTIFICATION

Identify storage directly

Explore identification built into racks or applied directly to metal storage.



Custom Engraving

FREEZER RACKS & METAL STORAGE BOXES

Custom identification can be engraved directly on freezer racks or aluminum storage boxes.

Contact us with your equipment and engraving requirements.

VIEW PRODUCT >



Laser Cut ID Freezer Racks

LASERCUTID-RACK

Stainless steel upright and vertical racks include a laser-cut identifier. Models hold standard 2-inch boxes.

Upright and LN2 / chest freezer configurations

VIEW PRODUCT >



Yellow Digits Riveted to Racks

BYFR

Yellow letters or numbers are riveted directly onto compatible freezer racks.

Up to 6 digits on standard racks; 1 on mini racks.

VIEW PRODUCT >



Custom Blood Bag Canisters

CUSTOM-CANISTER

Aluminum blood bag canisters offer custom engraving or a QR code for identification and tracking.

Multiple sizes available; confirm bag compatibility.

VIEW PRODUCT >

TAGS & PRINTING

Tags and on-demand labels

Identify rack handles or print labels for laboratory storage and samples.



Anodized Aluminum Tags | Letters

TAG-L-X-T

Round colored tags attach to rack handles; choose up to three letters.

2 in. diameter | Red, blue, green, or yellow

VIEW PRODUCT >



Anodized Aluminum Tags | Numbers

TAG-N-X-T

Round colored tags attach to rack handles; choose up to three digits.

2 in. diameter | Red, blue, green, or yellow

VIEW PRODUCT >



Engraved Stainless Steel Keychain

BEKC

Attach an engraved stainless steel keychain to an upright or vertical rack handle.

1.162 in. diameter | Up to three digits

VIEW PRODUCT >



Yellow Digits on Keychain

BYKC

A handle-mounted keychain with riveted yellow letters or numbers.

Up to four characters | Two-sided option available

VIEW PRODUCT >



LABeler™ Handheld Printer

L9010

Print laboratory labels and 2D barcodes using a handheld keypad and built-in cutter.

Standard tape to -100°C; ultra-low tape for LN2 storage.

VIEW PRODUCT >