ROEMER

PHARMAROAM60

The PharmaRoam line of mobile, battery-powered refrigerators and/or freezers enables healthcare providers to store vaccines and biologic products. By reducing medical waste and saving time, the PharmaRoam allows healthcare providers to transition from outdated disposable coolers to innovative, cost-effective solutions for the modern day.



THE FEATURES YOU NEED

- Outer case constructed with durable military-grade Aluminum (MIL-SPEC)
- Ergonomic spring handles automatically stow for safe shipping
- Stackable outer case allows for efficient out-of-use storage.
- Smart change technology automatically switches charge between battery and electrical power plugs without the need for user intervention
- Digital microprocessor temperature controller with built-in low battery alarm
- Watertight, sealed, an anodized payload compartment
- Complies with cold-chain protocols of VFC, FDA, CDC, AABB & AATB
- Pad lockable

SPECIFICATIONS

OVERALL DIMENSIONS	31.75" X 23.25" X 24.75"
WEIGHT	100 LBS
REFRIGERATED PAYLOAD DIMENSIONS	22" X 14" X 10"
REFRIGERATED PAYLOAD CAPACITY	9800 UNITS OF VACCINE DOSES OR 1230 SINGLE USE SYRINGES
PROGRAMMABLE TEMPERATURE RANGE	REFRIGERATOR: 2°C - 8°C FREEZER: -20°C - 0°C
TEMPERATURE UNIFORMITY	+/- 1.0° C AT SETPOINT 4.0° C
POWER INPUTS	12-24 VDC 110-240 VAC

BATTERY CAPACITY

	NEW TOTAL WEIGHT	FRIDGE BATTERY LIFE	FREEZER BATTERY LIFE
BBS50	145 LBS	36-40 HRS	8-10 HRS
BBS100	175 LBS	80-96 HRS	20-24 HRS
BBS200	245 LBS	160-196 HRS	36-48 HRS

OPTIONAL ACCESSORIES

- Cigarette lighter adapter
- Mobility options
- Temperature data logger